
# Instruction For Use [IFU]

Please select your language of choice  
from the list below:

EN

FR

DE

NL

01.132.396 V1

# Mounting Media for Coverslipping

## 1 INTENDED USE

Mounting media are designed for the purpose of permanently preserving stained cell preparations and histology tissue sections (2-5 micron thickness), derived from the human body. The preserved tissue sample is visualised microscopically to facilitate the diagnosis of the pathological state, and is intended for in-vitro examination of cellular deformity.

## 2 PRINCIPLE OF THE METHOD:

Mounting Mediums are viscous, clear liquids, consisting of a plastic dissolved in a solvent with a refractive index close to glass when dry. The mounting medium is dispensed onto a specimen in a liquid form, and the specimen is then sealed air-tight by addition of a coverslip on top of the specimen. The evaporation of the carrier solvent (Xylene, Toluene, Xylene/ Toluene substitute) results in the hardening of the mounting medium, to form a solid clear film between the slide/ specimen and the coverslip. As the mounting medium's refractive index is similar to glass, the specimen can be viewed under a microscope without any distortion of the image.

### Reagent Preparation:

Mounting Medium (DPX, ExPert and ExPert XTF) are ready to use reagents, no further preparation is required before use.

## 3 COMPONENTS:

### DPX:

**Chemical Name:** XYLENE (Ortho-xylene, meta-xylene, para-xylene & ethylbenzene)  
**CAS Number:** 1330-20-7  
**EC Number:** 215-535-7  
**Reach Reg. No.:** 01-2119488216-32  
**Conc. (%w/w):** 60 - 70%  
**Classification:** Flam. Liq. 3: H226; Acute Tox. 4: H332; Acute Tox. 4: H312; Skin Irrit. 2: H315  
 Ethylbenzene

**Chemical Name:** 100-41-4  
**CAS Number:** 202-849-4  
**Conc. (%w/w):** 25%  
**Classification:** Flam. Liq. 2: H226; Asp To. 1: H304  
 Acute Tox. 4: H332; STOT RE 2: H373  
 Aquatic Chronic 3: H412

### ExPert:

**Chemical Name:** XYLENE (Xylene o-, m-, p- or mixed isomers)  
**Index Number:** 601-022-00-9  
**CAS Number:** 1330-20-7  
**EC Number:** 215-535-7  
**Reach Reg. No.:** 01-2119488216-32  
**Conc. (%w/w):** 60 - 70%  
**Classification:** Flam. Liq. 3: H226; Acute Tox. 4: H332; Acute Tox. 4: H312; Skin Irrit. 2: H315; Eye Irrit. 2: H315; STOT SE 3: H335; STOT SE 3: H336; STOT RE 2: H373

### ExPert XTF:

**Chemical Name:** Hydrocarbons, C10-C12  
**EC Number:** 923-037-2  
**Reach Reg. No.:** Not yet assigned in the supply chain  
**Conc. (%w/w):** <45%  
**Classification:** Flam. Liq. 2: H225; Asp Tox. 1: H304; Aquatic Chronic 4: H413; EUH066

## 4 STORAGE AND SHELF LIFE AFTER FIRST OPENING:

### STORAGE AFTER FIRST OPENING:

Once opened the mounting medium can be used until the expiry date marked on the bottle, on condition that the bottle is kept tightly closed between uses to prevent evaporation of the mounting medium. Evaporation of the mounting medium may lead to a deterioration in the dispensing performance of automated coverslipping instruments, due to an increase in viscosity of the mounting medium.

### SAFETY INSTRUCTIONS:

Please refer to the Product Safety Data Sheet and/ or product label for instructions on handling, disposal, chemical warnings, precautions and product composition.

### STORAGE INSTRUCTIONS:

Opened bottles of mounting medium should be kept in a cool, dry, well ventilated area. Keep containers tightly closed. Store in correctly labelled containers. It is recommended that mounting medium is stored at a temperature between 4°C and 40°C. Keep away from direct sunlight. Keep away from incompatible materials: oxidising agents, acids and alkalis.

### HANDLING AND TREATMENT BEFORE USE:

Stained slides (cell preparations/ tissue sections) should be dehydrated routinely using graded alcohols, and then cleared using either an aromatic (Xylene or Toluene) clearing agent or XTF clearing agent and if required kept immersed in the same clearing/ agent prior to the specimen being mounted.

## 5 WARNINGS AND PRECAUTIONS:

ExPert XTF:

DPX and ExPert:

## 6 EXAMINATION & CONTROL PROCEDURES - RECOMMENDATIONS:

- Important:** When coverslipping either manually or using an automated coverslipping instrument it is important to allow a newly received bottle of mounting medium to settle for 8 hours prior to use, to allow any micro bubbles formed during transport to settle out.
- Use only in a well ventilated area. It is recommended that users manually coverslipping slides carry out the procedure within a fume hood.
- Slides should be well drained of clearing agent prior to mounting.
- Drop approximately 1 drop (0.5ml) of mounting medium onto the slide manually using a glass rod or pipette. Mounting medium may also be dispensed utilising an automated coverslipping instrument. Please consult the instrument's user manual for appropriate settings. For optimal performance, it is recommended that the prepared material is of an even thickness and that the mounting medium completely covers the material. Use only as much mounting medium as required to cover the specimen and coverslip to be added.
- If coverslipping manually gently add a clean coverslip and avoid trapping any air bubbles under the coverslip. The mounting medium should be handled with a minimum of disturbance to avoid the generation of air bubbles.
- Allow the slides to dry in a horizontal position for 1 hour prior to moving or viewing on a microscope, to avoid movement of the coverslip.
- Allow the mounted slides to dry for 48 hours (DPX and ExPert) and 72 hours (ExPert XTF) prior to filing. This ensures that slides are thoroughly dry and prevents any slides sticking together.
- The use of excessive heat to dry the slides is not recommended.
- Any excess mounting medium should be removed after positioning of the coverslip.
- Stained slides should be stored in a dark environment, away from sources of bright light. Slides mounted should be checked periodically, to confirm that the mountant is clear; with no sign of the mountant drying back or peeling off the slide.
- When using an automated coverslipping instrument, it is recommended that the device is cleaned regularly to avoid contamination of the mounting medium. (The manufacturer's user manual should be consulted for recommended cleaning procedures.)
- DPX, ExPert and ExPert XTF is a ready to use product, and any modification is done so at the user's own risk.

Specimens mounted with DPX, ExPert or ExPert XTF mounting mediums will be clearly visible, with no fogging, or distortions due to refractive index mismatch with glass microscope slides. The mounted specimen will be airtight for specimen preservation with the coverslip well bonded to the microscope slide.

## 7 PERFORMANCE CHARACTERISTICS AND METHOD LIMITATIONS:

DPX and ExPert are suitable for all staining techniques which are compatible with the use of alcohol and aromatic (Xylene/ Toluene) clearing agent. ExPert XTF is suitable for all techniques which are compatible with the use of alcohol and non-aromatic clearing agent.

It is recommended that the stained sample is cleared in the same solvent which is used in manufacture of the mounting medium to ensure optimal optical properties and transparency of the mounted specimen.

## 7 PERFORMANCE CHARACTERISTICS AND METHOD LIMITATIONS CONT.:

The coverslip and DPX, ExPert or ExPert XTF mounting medium can be removed by soaking a mounted slide in a corresponding clearing agent until the coverslip and mounting medium can be removed.

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