


Instructions

μ-Slides I Luer



The ibidi product family is comprised of a variety of μ-Slides and μ-Dishes, which have all been designed for high-end microscopic analysis of fixed or living cells. The high optical quality of the material is similar to that of glass, so you can perform all kinds of fluorescence experiments with uncompromised resolution and choice of wavelength.

The μ-Slides I Luer are designed for cell culture under perfusion and all flow applications.

Main applications are the simulation of blood vessels for arteriosclerosis research and applying defined shear stress and shear rates on cells inside the channel. The female Luers allow easy connections to tubing and pump systems. The μ-Slide I Luer comes in five versions which only differ in their channels' heights and channel volumes.

Material

ibidi μ-Slides, μ-Dishes, and μ-Plates are made of a plastic that has the highest optical quality. The bottom material exhibits extremely low birefringence and autofluorescence, similar to that of glass. Also, it is not possible to detach the bottom from the upper part. The μ-Slides, μ-Dishes, and μ-Plates are not autoclavable, since they are only temperature-stable up to 80°C/175°F. Please note that gas exchange between the medium and incubator's atmosphere occurs partially through the plastic bottom, which should not be covered.

Optical Properties ibidi Standard Bottom

Refractive index n_D (589 nm)	1.52
Abbe number	56
Thickness	No. 1.5 (180 μm)
Material	microscopy plastic

Geometry of the μ-Slides I Luer

The μ-Slides I Luer provide standard slide format according to ISO 8037/1.

General Dimensions

Number of Channels	1
Channel length	50 mm
Channel width	5.0 mm
Volume per reservoir	60 μl
Growth area	2.5 cm ² per channel
Bottom matches coverslip	No. 1.5

¹Collagen IV, BD Cat.-Nr. 35 6233, Fibronectin, BD Cat.-Nr. 354008, Poly-L-Lysin, Sigma Cat.-Nr. P4832, Poly-D-Lysin, BD Cat.-Nr. 354210

The channel volume differs, depending on the channel height (see table below).

Product name	Channel height	Channel volume
μ-Slide I ^{0.1} Luer	100 μm	25 μl
μ-Slide I ^{0.2} Luer	200 μm	50 μl
μ-Slide I ^{0.4} Luer	400 μm	100 μl
μ-Slide I ^{0.6} Luer	600 μm	150 μl
μ-Slide I ^{0.8} Luer	800 μm	200 μl

μ-Slide Surfaces

Depending on the type of cells and the special application you are using, you will need μ-Slides with different surfaces. If you do not require any special adhesion molecules for your application, the best choice will be ibiTreat, a tissue culture treated surface.

We provide precoated μ-Slides with adhesion substrates like Collagen IV, Fibronectin, Poly-L-Lysin, and Poly-D-Lysin. Such adhesion substrates have been shown to stimulate the adhesion and growth of various cell lines in μ-Slides. Only high-quality substrates are used ¹.

The uncoated μ-Slide is manufactured from hydrophobic plastic. For the cultivation of most cell lines, it is indispensable to treat the uncoated μ-Slide with biopolymers, which mediate cell adhesion and growth.

Coating your μ-Slides I Luer

The uncoated μ-Slide must be coated to promote cell adhesion. If you like to establish a certain coating for your demands we recommend to test your coating procedure on uncoated and ibiTreat μ-Slides, since we have observed that some biomolecules adhere differently to hydrophobic or hydrophilic plastic surfaces.

- Prepare your coating solution according to the manufacturer's specifications or reference.

Instructions

μ-Slides I Luer

- Apply the channel volume depending on the channel height. Leave at room temperature for at least 30 minutes.

Product name	Coating area	Growth area
μ-Slide I ^{0.1} Luer	5.1 cm ²	2.5 cm ²
μ-Slide I ^{0.2} Luer	5.2 cm ²	2.5 cm ²
μ-Slide I ^{0.4} Luer	5.4 cm ²	2.5 cm ²
μ-Slide I ^{0.6} Luer	5.6 cm ²	2.5 cm ²
μ-Slide I ^{0.8} Luer	5.8 cm ²	2.5 cm ²

- Aspirate the solution and wash with 1 ml ultra-pure water. Let dry at room temperature.

Further information about coatings are provided in [Application Note 08 "Cell culture coating"](#).

Cell Culture under Static Conditions

For many static applications with microscopic imaging, like transfection, immunofluorescence staining or cell morphology the μ-Slide I Luer is an optimal solution.

Important!

The μ-Slide I^{0.1} Luer and μ-Slide I^{0.2} Luer are not recommended for use in static cell culture!

- Trypsinize and count cells as usual. The cell density after seeding strongly depends on the channel's height. We recommend the following cell concentrations and volumes:

Product name	Volume	Cell concentration
μ-Slide I ^{0.1} Luer	25 μl	12–28 x 10 ⁵ cells/ml
μ-Slide I ^{0.2} Luer	50 μl	6–14 x 10 ⁵ cells/ml
μ-Slide I ^{0.4} Luer	100 μl	3–7 x 10 ⁵ cells/ml
μ-Slide I ^{0.6} Luer	150 μl	2–4.5 x 10 ⁵ cells/ml
μ-Slide I ^{0.8} Luer	200 μl	1.5–3.5 x 10 ⁵ cells/ml

- Apply the volume directly into the channel. The recommended cell concentration should result in a 50 % optical confluence layer after 24 hours.
- Cover reservoirs with the supplied caps. Incubate at 37°C and 5 % CO₂ as usual.
- After cell attachment fill each reservoir with 60 μl medium.

Depending on the cells we recommend exchanging the medium every day in static culture: Aspirate both reser-

voirs (not the channel). Flush fresh medium inside the channel by filling one reservoir with 120 μl medium and removing the content of the reservoir from the other well, ensuring the channel is never dry. Leave both reservoirs filled with approx. 60 μl each.

Tip:

The day before seeding the cells we recommend placing the cell medium, the μ-Slide, and the tubing into the incubator for equilibration. This will prevent the liquid inside the channel from emerging air bubbles over the incubation time.

Quick dispensing of cell suspension helps to avoid trapped air bubbles and leads to maximal homogeneity of cell distribution.

Cell Culture under Flow Conditions

Due to the Luer adapters, μ-Slide I Luer is suitable to any fluidic setup for cell cultivation under flow conditions. Cells are seeded into the channel and the flow is applied after cell attachment.

- Trypsinize and count cells as usual. The cell density after seeding strongly depends on the channel's height. We recommend the following cell concentrations and volumes:

Product name	Volume	Cell concentration
μ-Slide I ^{0.1} Luer	25 μl	5–10 x 10 ⁶ cells/ml
μ-Slide I ^{0.2} Luer	50 μl	2.5–5 x 10 ⁶ cells/ml
μ-Slide I ^{0.4} Luer	100 μl	1.2–2.5 x 10 ⁶ cells/ml
μ-Slide I ^{0.6} Luer	150 μl	0.8–1.6 x 10 ⁶ cells/ml
μ-Slide I ^{0.8} Luer	200 μl	0.6–1.2 x 10 ⁶ cells/ml

- Apply the volume directly into the channel. The recommended cell concentration should result in a 100 % optical confluence layer after some hours.
- Cover reservoirs with the supplied caps. Incubate at 37°C and 5 % CO₂ as usual.
- After cell attachment fill each reservoir with 60 μl medium.
- The μ-Slide is now ready for applying flow conditions on the adherent cells. Don't trap air bubbles when plugging in the connecting tubes.

Tip:

The day before seeding the cells we recommend placing the cell medium, the μ-Slide, and the tubing into the incubator for equilibration. This will prevent the liquid inside the channel from emerging air bubbles over the incubation time.


Quick dispensing of cell suspension helps to avoid trapped air bubbles and leads to maximal homogeneity of cell distribution.

For long term analysis of cells under flow conditions we recommend using μ-Slides with ibiTreat surface.

Application Note 13 "HUVECs under perfusion" describes a detailed protocol of a long term experiment with HUVECs and the ibidi pump system.

Detailed information about flow rates, shear stress, and shear rates is provided in **Application Note 11 "Shear stress and shear rates"** on www.ibidi.com.

Suitable flow kits (μ-Slide I Luer + tubing and adapters) are also available.



Please contact us for recommended perfusion setups. ibidi provides a variety of pump systems.

Preparation for Cell Microscopy

To analyze your cells, no special preparations are necessary. Cells can be observed live, or fixed directly in the μ-Slide on an inverted microscope. You can use any fixative of your choice. The μ-Slide material is compatible with a variety of chemicals (e.g., Acetone or Methanol). Further specifications can be found at www.ibidi.com. Due to the thin bottom of only 180 μm, high resolution microscopy is possible.

Immersion Oil

When using oil immersion objectives, use only the immersion oils specified in the table. The use of a non-recommended oil could lead to the damage of the plastic material and the objective.

Company	Product	Ordering Number
Zeiss	Immersion oil 518 F	(Zeiss) 444960
Zeiss	Immersion oil W 2010	(Zeiss) 444969
Leica	Immersion liquid	(Leica) 11513859

Instructions

μ-Slides I Luer

μ-Slide I Luer Family

The μ-Slide I Luer family is available with different channel heights and surfaces. See table below for choosing your μ-Slide I Luer.

μ-Slide I ^{0.1} Luer

Ordering Number	Treatment or Coating	Characteristics
81126	ibiTreat, tissue culture treated, sterile	hydrophilic, tissue culture treated
81122	Collagen IV, sterile	protein coating
81123	Fibronectin, sterile*	protein coating
81121	uncoated, sterile	hydrophobic

μ-Slide I ^{0.2} Luer

Ordering Number	Treatment or Coating	Characteristics
80166	ibiTreat, tissue culture treated, sterile	hydrophilic, tissue culture treated
80162	Collagen IV, sterile	protein coating
80163	Fibronectin, sterile*	protein coating
80164	Poly-L-Lysine, sterile	biopolymer coating
80165	Poly-D-Lysine, sterile*	biopolymer coating
80161	uncoated, sterile	hydrophobic

μ-Slide I ^{0.4} Luer

Ordering Number	Treatment or Coating	Characteristics
80176	ibiTreat, tissue culture treated, sterile	hydrophilic, tissue culture treated
80172	Collagen IV, sterile	protein coating
80173	Fibronectin, sterile*	protein coating
80174	Poly-L-Lysine, sterile	biopolymer coating
80175	Poly-D-Lysine, sterile*	biopolymer coating
80171	uncoated, sterile	hydrophobic

μ-Slide I ^{0.6} Luer

Ordering Number	Treatment or Coating	Characteristics
80186	ibiTreat, tissue culture treated, sterile	hydrophilic, tissue culture treated
80182	Collagen IV, sterile	protein coating
80183	Fibronectin, sterile*	protein coating
80184	Poly-L-Lysine, sterile	biopolymer coating
80185	Poly-D-Lysine, sterile*	biopolymer coating
80181	uncoated, sterile	hydrophobic

μ-Slide I ^{0.8} Luer

Ordering Number	Treatment or Coating	Characteristics
80196	ibiTreat, tissue culture treated, sterile	hydrophilic, tissue culture treated
80192	Collagen IV, sterile	protein coating
80193	Fibronectin, sterile*	protein coating
80194	Poly-L-Lysine, sterile	biopolymer coating
80195	Poly-D-Lysine, sterile*	biopolymer coating
80191	uncoated, sterile	hydrophobic

* available on request only

Selected References

- M. Anton, A. Wolf, O. Mykhaylyk, C. Koch, B. Gansbacher, and C. Plank. Optimizing Adenoviral transduction of endothelial cells under flow conditions. *Pharmaceutical Research*, 2012. doi: 10.1007/s11095-011-0631-2.
- C. Boehlke, F. Kotsis, V. Patel, S. Braeg, H. Voelker, S. Bredt, T. Beyer, H. Janusch, C. Hamann, M. Gödel, K. Müller, M. Herbst, M. Hornung, M. Doerken, M. Köttgen, R. Nitschke, P. Igarashi, G. Walz, and E. W. Kuehn. Primary Cilia Regulate mTORC1 Activity and Cell Size Through Lkb1. *Nature Cell Biology*, 2010. doi: 10.1038/ncb2117.
- H. L. Morris, C. I. Reed, J. W. Haycock, and G. C. Reilly. Mechanisms of Fluid-Flow-Induced Matrix Production in Bone Tissue Engineering. *Proceedings of the Institution of Mechanical Engineers, Part H: Journal of Engineering in Medicine*, 2010. doi: 10.1243/09544119jeim751.
- A. Sabine, Y. Agalarov, H. Maby-El Hajjami, M. Jaquet, R. Hägerling, C. Pollmann, D. Bebbber, A. Pfenniger, N. Miura, O. Dormond, J.-M. Calmes, Ralf H. Adams, T. Mäkinen, F. Kiefer, B. Kwak, and T. Petrova. Mechanotransduction, PROX1, and FOXC2 Cooperate to Control Connexin37 and Calcineurin during Lymphatic-Valve Formation. *Developmental Cell*, 2012.
- C. Schulz, E. Heiss, F. Gaertner, M. Orban, M.-L. v. Bruehl, P. Schramm, and S. Massberg. Novel Methods for Assessment of Platelet and Leukocyte Function Under Flow—Application of Epifluorescence and Two-Photon Microscopy in a Small Volume Flow Chamber Model. *The Open Biology Journal*, 2009. doi: 10.2174/1874196700902010130.
- T. G. Walsh, R. P. Murphy, P. Fitzpatrick, K. D. Rochfort, A. F. Guinan, A. Murphy, and P. M. Cummins. Stabilization of Brain Microvascular Endothelial Barrier Function by Shear Stress Involves VE cadherin Signaling Leading to Modulation of pTyr Occludin Levels. *Journal of Cellular Physiology*, 2011. doi: 10.1002/jcp.22655.

For research use only!

Further technical specifications can be found at www.ibidi.com. For questions and suggestions please contact us by e-mail info@ibidi.de or by telephone +49 (0)89/520 46 17 0. All products are developed and produced in Germany.

© ibidi GmbH, Am Klopferspitz 19, 82152 Martinsried, Germany.