MODE D'EMPLOI

1. GENERALITES

A. Electrodes de verre

Préparation:

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A. Glass electrodes

Preparation

1. Shake the electrodes gently to remove possible air bubbles ensuring that the internal buffer solution covers the whole inside surface of the membrane.
2. Clean the membrane by dipping the electrode into a 0.1 M NaOH solution for about 5 seconds immediately followed by rinsing with water.
3. a. New pH electrodes should be conditioned by soaking them in water acidified with HCl (pH 4-5) for some hours (e.g. overnight).
- b. pH and pH electrodes should be conditioned in an acidified (with HCl to pH 4-5) 0.1 M NaCl (KCl) respectively aqueous solution during 24 or preferably 48 hours.

Storage

4. a. After use the pH electrode is stored in the solution mentioned under 3.a or in distilled water.- b. pH and pH electrodes are stored in the solution mentioned under 3.b.

Reactivation

5. Drying out of the membrane or operation under extreme conditions may cause a decrease of the response time of the glass electrode.

a. Reactivation of the pH electrode is performed as indicated under point 2 followed by soaking the electrode in a 0.1 M HCl solution for several hours.

b. Reactivation of pH and pH electrodes is performed according to points 2 and 3.b.

(In case the above described treatment does not have the expected result, the pH, pH and pH electrodes may be reactivated by dipping them into a 2% HF-solution for 5 to 10 seconds, immediately followed by washing with much water. This drastic treatment should not be applied frequently since it shortens the life time of the electrodes.)

B. Reference electrodes

Preparation

1. Remove salt crystals from the exterior by washing with water.
2. Remove or puncture the stopper on the electrolyte filling opening to obtain pressure equality.
3. The reference electrode must always be filled with the appropriate electrolyte so that the galvanic element cannot dry out. If necessary refill.

The electrolyte required is:

- generally: saturated KCl-solution (aq.) (20% C)
- low temperature electrodes: saturated NaCl-solution (aq.) (20% C)
- high temperature combined electrodes: 1.0 M KNO₃ + 0.01 M KCl solution (aq.)

for pH and pH measurements: 0.1 M NH₄Cl solution (aq.)

Storage

4. a. The reference electrode can be stored in dry condition with closed electrolyte inlet to prevent loss of electrolyte through the porous plug.
- b. The reference electrode can be stored in water with open electrolyte inlet ensuring continuous outflow of electrolyte through the porous plug to prevent entry of water into the electrode.

C. Combined glass/reference electrodes

These types of electrodes should be treated according to the instructions stated under A. and B.

II. PRACTICAL INFORMATION

1. Calibration of the measuring system should be carried out daily, weekly or monthly depending on the measuring accuracy required.
2. To ensure good measuring results the electrolyte in the reference electrode should be refilled every six months.
3. A positive outflow of the electrolyte is required (important particularly in industrial processes where overpressure applied to the electrolyte may be necessary).

4. The life time of the pH glass electrode is, generally speaking, two years. High temperatures, however, decrease this time considerably.

5. The meaning of pH measurements above pH 14 and below pH 0 should be examined carefully (in such cases low pH conductivity measurements may give more reliable information).

6. pH (pH) measurement due to the disturbing influence of the pH-value of the solution under test, the presence of relatively high concentrations of other cations such as K⁺, Na⁺ and NH₄⁺ ions, and the difficulties in-

olved with the preparation of reliable standard solutions with very low Na⁺ (K⁺) concentrations, it is advisable to abstain from pH (pH) measurements above pH (pH) 5.

7. The interpretation of pH-measurements in organic solutions should be handled with care. In aqueous solutions of organic compounds it is advisable to use reference electrodes with built-in salt bridges. To reduce the liquid junction potential, the following electrolyte should consist of a salt e.g. NH₄NO₃, LiCl or KNO₃ in an organic solvent.

8. With a view to the required high insulation resistance of the glass electrode, it is necessary that the upper part of the electrode, the connection cable and the connector plug are kept dry and clean.

III. CALIBRATION OF THE MEASURING SYSTEM

1. The calibration of a pH-measuring system is described in the directions of use of the pH-meter concerned. For this procedure reliable buffer solutions should be used. Well-known and commonly used are the buffer solutions of the National Bureau of Standards (N.B.S.).

2. The electrical zero-point of pH-measuring chains is at pH 7 and the slope close to 59 mV/pH at 25°C.

3. For the calibration of pH or pH electrodes a procedure identical to the pH calibration is performed, using pH or pH standard solutions instead of pH-buffer solutions. These standard solutions have a known Na⁺ or K⁺ concentration and a buffered pH-value of approx. pH 8. Such standard solutions composed of triethanolamine (methanolamine), HCl and NaCl (resp. KCl) solutions are delivered with pH (pH) electrodes is at pH (pH) 2. The slope of pH electrodes is approx. 59 mV/pH and of pH electrodes approx. 56 mV/pH both at 25°C.

4. The response of the pH electrodes is usually linear in the range pH 6-12. However, due to the liquid junction potential of the reference electrode, deviations are possible.

5. For accurate pH measurements it is to be advised to use for the calibration procedure the following buffer solutions:

- pH 6.86 (25°C) for A.P. compensation
- pH 4.01 (25°C) for mV/pH adjustment for measurements in the solid range (slope correction)
- pH 9.18 (25°C) for mV/pH adjustment for measurements in the alkaline range (slope correction)
- Identically for pH (pH) measurements a standard solution with pH (resp. pH) 2 is used for A.P. compensation and a second standard solution corresponding to the pH (pH) range to be measured is chosen for mV/pH (resp. mV/pH) correction.

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