For USA: For laboratory use.

For EU: **IVD**

LightCycler[®] Capillaries

(100 µl) **M**^{GRADE}

 **Version 04**

August 2013

Reaction vessels for the LightCycler[®] 2.0 Instrument

REF 03 612 066 001

Store at +15 to +25°C

Table of Contents

1.	What this Product Does	3
	Contents	3
	Storage and Stability	3
	Additional Equipment and Reagents Required	3
	Intended Use	3
	Product Characteristics	4
2.	How to Use this Product	5
3.	Additional Information on this Product	8
	How this Product Works	8
	References	8
4.	Supplementary Information	9
4.1	Conventions	9
	Text Conventions	9
	Symbols	9
4.2	Changes to Previous Version	9
4.3	Ordering Information	10
4.4	Notice to Purchaser	11
4.5	Trademarks	11
4.6	Regulatory Disclaimer	11

1. What this Product Does

Contents	8 boxes, each with 96 capillaries and stoppers
Storage and Stability	If stored at +15 to +25°C, the product is stable until the expiration date printed on the label. 🕒 The product is shipped at ambient temperature.
Additional Equipment and Reagents Required	Refer to the list below for additional reagents and equipment required, but not provided, to perform PCR with the LightCycler® Capillaries (100 µl) M^{GRADE} using the LightCycler® Carousel-Based System: ▪ LightCycler® 2.0 Instrument* ▪ LightCycler® 2.0 Sample Carousel (100 µl) (provided with the LightCycler® 2.0 Instrument) ▪ Standard benchtop microcentrifuge containing a rotor for 2.0 ml reaction tubes (e.g., Biofuge 19, Heraeus Instruments or equivalent) 🕒 The LightCycler® Carousel-Based System is provided with centrifuge adapters that allow LightCycler® Capillaries (100 µl) M^{GRADE} to be centrifuged in a standard microcentrifuge rotor. or ▪ LC Carousel Centrifuge 2.0* for use with the LightCycler® 2.0 Sample Carousel (100 µl) [optional] ▪ LightCycler® Capping Tool (provided with the LightCycler® Carousel-Based Systems) ▪ LightCycler® 2.0 Capillary Releaser (provided with the LightCycler® 2.0 Instrument) <i>* available from Roche Applied Science</i>
Intended Use	The LightCycler® Capillaries (100 µl) M^{GRADE} are the appropriate reaction vessels for the LightCycler® 2.0 Instrument when working with large reaction volumes. The capillary can hold a reaction volume of up to 100 µl. High-quality borosilicate glass ensures superb PCR performance and optimal fluorescence transmission. 🕒 The LightCycler® Capillaries (100 µl) M^{GRADE} are compatible with the LightCycler® 2.0 Sample Carousel (100 µl) only. ⚠️ If the LightCycler® Capillaries (100 µl) M^{GRADE} are used in combination with the LightCycler® 2.0 Instrument and Roche Diagnostics LightCycler® kits for <i>in vitro</i> diagnostic use (not available in all countries), adhere to the instructions for use in the Instructions for Use of these kits.


1. What this Product Does, continued

Product Characteristics

A glass capillary is molded to a polypropylene reservoir, permitting easy loading of the capillary. A stopper provides a secure seal, significantly reducing the risk of contamination.

Each LightCycler® Capillary (100 µl) **M**GRADE is 51 mm long (without plastic stopper) and has an outer diameter of 3.175 mm.

Ⓚ **M**GRADE products are designed for highly sensitive bacterial and fungal-nucleic acid detection. Proprietary manufacturing processes were developed in order to eliminate contaminations that might interfere with such assays.

2. How to Use this Product

During set-up of a LightCycler® 2.0 Instrument PCR run, PCR components (reagent mix and template nucleic acid, either separately or combined) are pipetted into the capillary reservoir and then forced into the glass capillary by centrifugation. You can perform this centrifugation step in two different ways:

- by using the LightCycler® Centrifuge Adapters, in combination with a standard benchtop microcentrifuge, or
- by using the LightCycler® 2.0 Sample Carousel (100 µl), in combination with the LC Carousel Centrifuge 2.0.

⚠ Do not touch the surface of the capillaries. Always wear gloves when handling the capillaries and the stoppers.

- 1 Program the experimental protocol and define the sample numbers, names etc. Refer to the 'Software' section of the LightCycler® 2.0 Instrument Operator's Manual for detailed instructions.

⚠ When programming the experimental protocol, ensure that you select "100 µl" as the Capillary Size for the experiment.

- 2 Prepare a PCR reagent mix and omit the template nucleic acid. Details on preparing the PCR mix are described in the corresponding Instructions for Use of the LightCycler® reagent kit in use.
-

- 3 Mix by gentle vortexing.
-

- 4
 - Place the capillaries into the LightCycler® Centrifuge Adapters that have been precooled in the cooling block.
 - Pipette the PCR mix into the plastic reservoir at the top of the capillary. Add the template nucleic acid to the capillary.
-

- 5 Seal each capillary with a plastic stopper using the LightCycler® Capping Tool.

⚠ Ensure that each LightCycler® Capillary (100 µl) **M^{GRADE}** is closed tightly by checking it visually: The lower part of the plastic stopper must be completely inserted into the glass capillary.

6 Preparing a PCR run using the LightCycler® Centrifuge Adapters

- Place the centrifuge adapters each containing a LightCycler® Capillary (100 µl) **M_{GRADE}** into the rotor of a standard benchtop microcentrifuge.
 - ⚠ Use only rotors that are designed to hold 2.0 ml reaction tubes.
 - ⚠ Place the centrifuge adapters in a balanced arrangement within the centrifuge. If centrifuging an uneven number of capillaries, use an empty centrifuge adapter as counterweight.
- Centrifuge briefly, at not more than $735 \times g$ for 5 s.
- Place the capillaries into the LightCycler® 2.0 Sample Carousel, keeping the capillaries in an upright position. Make sure that all capillaries are fixed in the optimal position where the O-ring of the LightCycler® 2.0 Sample Carousel covers the lower part of the plastic reservoir. Ensure proper positioning by lightly pressing the cap until you hear a final “click” as the capillary reaches its final position.
- ⚠ When pressing the capillaries into the LightCycler® 2.0 Sample Carousel, do not press too hard. Do not use capillaries that show slight cracks or cracked slightly when inserted into the LightCycler® 2.0 Sample Carousel. Do not use capillaries that have been dropped. Capillaries with slight cracks, possibly invisible, might rupture during a run.

Preparing a PCR run using the LC Carousel Centrifuge 2.0

- Place the capillaries into the LightCycler® 2.0 Sample Carousel, keeping the capillaries in an upright position. Make sure that all capillaries are fixed in the optimal position where the O-ring of the LightCycler® 2.0 Sample Carousel covers the lower part of the plastic reservoir. Ensure proper positioning by lightly pressing the cap until you hear a final “click” as the capillary reaches its final position.
 - ⚠ When pressing the capillaries into the LightCycler® 2.0 Sample Carousel, do not press too hard. Do not use capillaries that show slight cracks or cracked slightly when inserted into the LightCycler® 2.0 Sample Carousel. Do not use capillaries that have been dropped. Capillaries with slight cracks, possibly invisible, might rupture during a run.
 - Put the loaded LightCycler® 2.0 Sample Carousel into the rotor bucket and place it in the LC Carousel Centrifuge 2.0.
 - ⚠ Only operate the LC Carousel Centrifuge 2.0 with the blue rotor and blue rotor bucket supplied together with the instrument.
 - 📖 Refer to the LC Carousel Centrifuge 2.0 Operator's Manual for detailed operating instructions.
-

- 7** Place the LightCycler® 2.0 Sample Carousel into the LightCycler® 2.0 Instrument. Ensure that the notch below sample position 1 on the LightCycler® 2.0 Sample Carousel locks into position against the pin on the thermal chamber. Visually check that the carousel is inserted correctly and fits perfectly in the thermal chamber.

 Before placing the LightCycler® 2.0 Sample Carousel into the LightCycler® 2.0 Instrument:

- Visually check that all capillaries have a uniform fill level or if capillary breakage has occurred. Any indications of an uneven fill level will point to possible leakage of a capillary. Discard the affected capillaries, as leakage could impair the experimental result especially in quantification experiments.

If capillary leakage or breakage occurs, clean the LightCycler® Centrifuge Adapters or the LightCycler® 2.0 Sample Carousel with a lab tissue soaked with decontamination solution. Use commercialized reagents such as LTK-008 (Biodelta) or DNAZap (Ambion). The LightCycler® 2.0 Sample Carousel can also be autoclaved. Clean the chamber of the benchtop centrifuge or the LC Carousel Centrifuge 2.0 with a mild commercial detergent. If necessary, use 70% ethanol for disinfecting the chamber. Refer to the LightCycler® 2.0 Instrument Operator's Manual and the LC Carousel Centrifuge 2.0 Operator's Manual for detailed instructions.

- Ensure that the thermal chamber is clean and free of any items that could interfere with the capillaries during the run.

-
- 8** Close the instrument lid. You are now ready to start the run.

-
- 9** At the end of the run, release all capillaries placed into the LightCycler® 2.0 Sample Carousel using the LightCycler® 2.0 Capillary Releaser:

Place the loaded carousel on the LightCycler® 2.0 Capillary Releaser, then push the carousel down with the flat of your hand. The capillaries will be released from the rubber O-ring and held at a slightly elevated position. You can now easily remove the capillaries from the sample carousel.

-
- 10** After use, discard the capillaries into a solid waste box.
-

3. Additional Information on this Product

How this Product Works

In the LightCycler® 2.0 Instrument, temperature cycling is achieved using air. Since air has virtually no mass, this process is significantly faster (ten times) than cycling with conventional thermal blocks. A heating coil controlled by thermocouples in a thermal chamber warms the air. PCR occurs in the specially designed LightCycler® borosilicate glass capillaries. The capillaries have a high surface-to-volume ratio, to ensure rapid equilibration between the air and the reaction components (1, 2, 3, 4).

The combination of using air for rapid thermal cycling and the high surface-to-volume ratio of the capillaries, allows a single PCR cycle to be completed in less than 95 seconds [when using the LightCycler® Capillaries (100 µl) **M^{ERAD}**]. Depending on the experimental protocol, a complete PCR run of 30 to 40 cycles can thus be performed in 45 to 65 minutes.

The optical properties of borosilicate glass make the capillaries ideally suited for use as cuvettes for fluorescence measurement (5).

References

- 1 Wittwer, CT. *et al.* (1989). Automated polymerase chain reaction in capillary tubes with hot air. *Nucleic Acids Res.* **17**, 4353-4357.
- 2 Wittwer, CT. *et al.* (1990). Minimizing the time required for DNA amplification by efficient heat transfer to small samples. *Anal. Biochem.* **186**, 328-331.
- 3 Wittwer, CT. and Garling, DJ. (1991). Rapid cycle DNA amplification: time and temperature optimization. *BioTechniques* **10**, 76-83.
- 4 Wittwer, CT. *et al.* (1994). Rapid cycle DNA amplification. In: KB. Mullis, F. Ferre and RA. Gibbs. *The Polymerase Chain Reaction*. Birkhäuser, Boston, MA, pp 174-181.
- 5 Pritham, GH. and Wittwer, CT. (1998). Continuous fluorescent monitoring of rapid cycle polymerase chain reaction. *J. Clin. Ligand Assay* **21**, 404-412.

4. Supplementary Information

4.1 Conventions

Text Conventions To make information consistent and memorable, the following text conventions are used in this document:

Text Convention	Use
Numbered stages labeled ①, ② etc.	Stages in a process that usually occur in the order listed.
Numbered instructions labeled ①, ② etc.	Steps in a procedure that must be performed in the order listed.
Asterisk *	Denotes a product available from Roche Applied Science.

Symbols

In this document and the product label, the following symbols are used to highlight important information:

Symbol	Description
	Information Note: Additional information about the current topic or procedure.
	Important Note: Information critical to the success of the procedure or use of the product.
	This product fulfills the requirements of the European Directive 98/79/EC for <i>in vitro</i> diagnostic medical devices.
	Manufacturer
	Use by
	Temperature limitation
	Batch code
	Catalog number
	<i>In vitro</i> diagnostic medical device
	Do not re-use
	Consult instructions for use.

4.2 Changes to Previous Version

- Revision of regulatory disclaimer.
- Editorial changes.
- Section Quality Control deleted.

4. Supplementary Information, continued

4.3 Ordering Information

Roche Applied Science offers a large selection of reagents and systems for life science research. For a complete overview of related products and manuals, please visit and bookmark our home page, <http://www.roche-applied-science.com> and our Special Interest Sites including:

- The LightCycler® 2.0 System: <http://www.lightcycler-online.com/>
- DNA & RNA preparation – Versatile Tools for Nucleic Acid Purification: <http://www.roche-applied-science.com/napure>
- The MagNA Pure System family for automated nucleic acid isolation: <http://www.magnapure.com>

Product	Pack Size	Cat. No.
LightCycler® 2.0 Instrument ¹⁾	1 instrument plus accessories	03 531 414 001
LC Carousel Centrifuge 2.0 ¹⁾	1 centrifuge plus rotor and bucket	03 709 507 001 (115 V) 03 709 582 001 (230 V)
LightCycler® Centrifuge Adapters	1 set	11 909 312 001
LightCycler® Capping Tool	1 tool	03 357 317 001
LightCycler® 2.0 Sample Carousel (100 µl)	1 sample carousel	03 603 954 001
LightCycler® 2.0 Capillary Releaser	1 capillary releaser	03 603 920 001

- Kits for *in vitro* diagnostic use:

For Roche Diagnostics LightCycler® kits for *in vitro* diagnostic use (not available in all countries), please contact your local Roche representative.

¹⁾ The Instruments fulfill the requirements of the European Directive 98/79/EC for *in vitro* diagnostic medical devices.

4. Supplementary Information, continued

4.4 Notice to Purchaser

For patent license limitations for individual products please refer to:
www.technical-support.roche.com.

4.5 Trademarks

LIGHTCYCLER, LC and MAGNA PURE are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

4.6 Regulatory Disclaimer

For USA: For laboratory use.

For EU: **IVD**

Contact and Support

If you have questions or experience problems with this or any Roche Applied Science (RAS) product, please contact our Technical Support staff. Our scientists commit themselves to providing rapid and effective help.

We also want you to contact us if you have suggestions for enhancing RAS product performance or using our products in new or specialized ways. Such customer information has repeatedly proven invaluable to RAS and the world-wide research community.

To ask questions, solve problems, suggest enhancements or report new applications, please visit our **Online Technical Support** Site at:

www.roche-applied-science.com/support

To call, write, fax, or email us, visit the Roche Applied Science home page, www.roche-applied-science.com, and select your home country. Country-specific contact information will be displayed.

On the Roche Applied Science home page select **Printed Materials** to find:

- in-depth Technical Manuals
- Lab FAQs: Protocols and references for life science research
- Material Safety Data Sheets
- Package Inserts and Product Instructions

or to request hard copies of printed materials.

Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim, Germany,
+49 621 7590 www.roche-applied-science.com

Distribution in USA: Roche Diagnostics Corporation
9115 Hague Road
PO Box 50457
Indianapolis, IN 46250, USA
US Customer Technical Support: 0800-262-4911

Roche Diagnostics GmbH
Sandhofer Strasse 116
D-68305 Mannheim
Germany

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LightCycler[®] Capillaries (100 µl) **M**^{GRADE}

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Referenser

- 1 Wittwer, CT. *et al.* (1989). Automated polymerase chain reaction in capillary tubes with hot air. *Nucleic Acids Res.* **17**, 4353–4357.
- 2 Wittwer, CT. *et al.* (1990). Minimizing the time required for DNA amplification by efficient heat transfer to small samples. *Anal. Biochem.* **186**, 328–331.
- 3 Wittwer, CT. and Garling, DJ. (1991). Rapid cycle DNA amplification: time and temperature optimization. *BioTechniques* **10**, 76–83.
- 4 Wittwer, CT. *et al.* (1994). Rapid cycle DNA amplification. In: KB. Mullis, F. Ferre and RA. Gibbs. *The Polymerase Chain Reaction*. Birkhäuser, Boston, MA, pp 174–181.
- 5 Pritham, GH. and Wittwer, CT. (1998). Continuous fluorescent monitoring of rapid cycle polymerase chain reaction. *J. Clin. Ligand Assay* **21**, 404–412.

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- 1 Wittwer, CT. *et al.* (1989). Automated polymerase chain reaction in capillary tubes with hot air. *Nucleic Acids Res.* **17**, 4353-4357.
- 2 Wittwer, CT. *et al.* (1990). Minimizing the time required for DNA amplification by efficient heat transfer to small samples. *Anal. Biochem.* **186**, 328-331.
- 3 Wittwer, CT. and Garling, DJ. (1991). Rapid cycle DNA amplification: time and temperature optimization. *BioTechniques* **10**, 76-83.
- 4 Wittwer, CT. *et al.* (1994). Rapid cycle DNA amplification. In: KB, Mullis, F. Ferre and RA. Gibbs. *The Polymerase Chain Reaction*. Birkhäuser, Boston, MA, pp 174-181.
- 5 Pritham, GH. and Wittwer, CT. (1998). Continuous fluorescent monitoring of rapid cycle polymerase chain reaction. *J. Clin. Ligand Assay* **21**, 404-412.

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Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim, Germany, +49 621 7590 www.roche-applied-science.com