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Gamme de volumes disponibles:

0,5 – 10 µl
10 – 100 µl
30 – 300 µl

Chemical Resistance ErgoOne Nose Cones

	Concentration [%]	Single Channel		Multichannel
		Micro Cone	Standard Cone	
Acids and Alcalis				
Ammonia Solution	25,0	R	R	R
Ammonia Solution	2,0	R	R	R
Acetic Acid	96,0	R	R	R
Acetic Acid	12,0	R	R	R
Sodium Hydroxide	20,0	R	R	NR
Sodium Hydroxide	4,0	R	R	R
Perchloric Acid	10,0	R	R	R
Nitric Acid	65,0	R	R	NR
Nitric Acid	6,3	R	R	R
Hydrochloric Acid	32,0	NR	R	R
Hydrochloric Acid	4,0	NR	R	R
Sulfuric Acid	95,0	R	R	R
Sulfuric Acid	16,0	R	R	R
Trichloroacetic Acid	40,0	R	R	R
Trichloroacetic Acid	10,0	R	R	R
Trifluoroacetic Acid	10,0	R	R	R
Trifluoroacetic Acid	100,0	R	R	R
Organic Solvents				
Acetone		NT	R	R
Acetonitrile		NT	R	R
Cleaning Solvent		NT	R	R
Chloroform		NT	R	R
Diethyl Ether		NT	R	R
DMSO (Dimethyl Sulfoxide)	10,0	NT	R	R
DMSO (Dimethyl Sulfoxide)	50,0	NT	R	R
DMSO (Dimethyl Sulfoxide)	100,0	NT	R	R
Acetic Ether		NT	R	R
Ethanol		NT	R	R
Formaldehyde	40,0	NT	R	R
Isoamyl Alcohole		NT	R	R
Isopropanol		NT	R	R
Methanol		NT	R	R
Petroleum Ether		NT	R	R
Phenol (water saturated)		NT	R	R
Carbon Tetrachloride		NT	R	R
Toluol		NT	R	R

R = resistant

NR = not resistant

NT = Not tested

Chemical Resistance ErgoOne Nose Cones

	Concentration	Single Channel		Multichannel
		Micro Cone	Standard Cone	
Cleaning and decontamination products				
Count off (Decontamination Product)		NT	R	R
Dismozone pure (Peroxide based)		NT	R	R
DNA away		NT	R	R
Helipur (Phenol based)		NT	R	R
Hexaquart S (QAV based)		NT	R	R
Hi - TOR Plus (QAV based)		NT	R	R
Korsolex basic (Aldehyde based)		NT	R	R
Meliseptol (Alcohol based)		NT	R	R
Sodium Hypochlorite	12,0	NT	R	R
RNAse away		NT	R	R
Sterillium (Alcohol based)		NT	R	R
Hydrogen Peroxide	35,0	NT	R	R
Cidex Activated Dialdehyde Solution (Aldehyde based)		NT	R	R
Other				
Cesium Chloride	saturated	NT	R	R
EDTA; pH 8	1,8 g/ml	NT	R	R
Ficoll	1,077 g/l	NT	R	R
Formamide	50,0%	NT	R	R
Glutaraldehyde	25,0%	NT	R	R
Glycerol	50,0%	NT	R	NR
Guanidine Thiocyanate	4 M	NT	R	R
Mineral Oil		NT	R	R
Sodium Acetate; pH 5,2	2 M	NT	R	R
Paraffin Oil		NT	R	R
SDS	1,0%	NT	R	R
TRIS-Buffer; pH 7,5	1 M	NT	R	R
Triton-X-100	1,0%	NT	R	R
Tween 20	1,0%	NT	R	R
Water		NT	R	R
Xylene		NT	R	R

R = resistant

NR = not resistant

NT = Not tested

Chemical Resistance ErgoOne Plunger (Piston)

	Concentration [%]	Single Channel			Multichannel	
		up to 20 µl	up to 1000 µl	5000 µl	up to 10 µl	up to 300 µl
Acids and Alcalis						
Ammonia Solution	25,0	R	R	R	R	R
Ammonia Solution	2,0	R	R	R	R	R
Acetic Acid	96,0	R	R	R	R	R
Acetic Acid	12,0	R	R	R	R	R
Sodium Hydroxide	20,0	R	NR	R	R	NR
Sodium Hydroxide	4,0	R	NR	R	R	NR
Perchloric Acid	10,0	R	R	R	R	R
Nitric Acid	65,0	R	NR	R	R	NR
Nitric Acid	6,3	R	R	R	R	R
Hydrochloric Acid	32,0	NR	R	R	NR	R
Hydrochloric Acid	4,0	NR	R	R	NR	R
Sulfuric Acid	95,0	R	R	R	R	R
Sulfuric Acid	16,0	R	R	R	R	R
Trichloroacetic Acid	40,0	R	R	R	R	R
Trichloroacetic Acid	10,0	R	R	R	R	R
Trifluoroacetic Acid	10,0	R	R	R	R	R
Trifluoroacetic Acid	100,0	R	R	R	R	R
Organic Solvents						
Acetone		NT	R	R	NT	R
Acetonitrile		NT	R	R	NT	R
Cleaning Solvent		NT	R	R	NT	R
Chloroform		NT	R	R	NT	R
Diethyl Ether		NT	R	R	NT	R
DMSO (Dimethyl Sulfoxide)	10,0	NT	R	R	NT	R
DMSO (Dimethyl Sulfoxide)	50,0	NT	R	R	NT	R
DMSO (Dimethyl Sulfoxide)	100,0	NT	R	R	NT	R
Acetic Ether		NT	R	R	NT	R
Ethanol		NT	R	R	NT	R
Formaldehyde	40,0	NT	R	R	NT	R
Isoamyl Alcohole		NT	R	R	NT	R
Isopropanol		NT	R	R	NT	R
Methanol		NT	R	R	NT	R
Petroleum Ether		NT	R	R	NT	R
Phenol (water saturated)		NT	R	R	NT	R
Carbon Tetrachloride		NT	R	R	NT	R
Toluol		NT	R	R	NT	R

R = resistant

NR = not resistant

NT = Not tested

Chemical Resistance ErgoOne Plunger (Piston)

	Concentration	Single Channel			Multichannel	
		up to 20 µl	up to 1000 µl	5000 µl	up to 10 µl	up to 300 µl
Cleaning and decontamination products						
Count off (Decontamination Product)		NT	R	R	NT	R
Dismozone pure (Peroxide based)		NT	R	R	NT	R
DNA away		NT	R	R	NT	R
Helipur (Phenol based)		NT	R	R	NT	R
Hexaquart S (QAV based)		NT	R	R	NT	R
Hi - TOR Plus (QAV based)		NT	R	R	NT	R
Korsolex basic (Aldehyde based)		NT	R	R	NT	R
Meliseptol (Alcohol based)		NT	R	R	NT	R
Sodium Hypochlorite	12,0%	NT	R	R	NT	R
RNAse away		NT	R	R	NT	R
Sterillium (Alcohol based)		NT	R	R	NT	R
Hydrogen Peroxide	35,0%	NT	R	R	NT	R
Cidex Activated Dialdehyde Solution (Aldehyde based)		NT	R	R	NT	R
Other						
Cesium Chloride	saturated	NT	R	R	NT	R
EDTA; pH 8	1,8 g/ml	NT	R	R	NT	R
Ficoll	1,077 g/l	NT	R	R	NT	R
Formamide	50,0%	NT	R	R	NT	R
Glutaraldehyde	25,0%	NT	R	R	NT	R
Glycerol	50,0%	NT	NR	R	NT	NR
Guanidine Thiocyanate	4 M	NT	R	R	NT	R
Mineral Oil		NT	R	R	NT	R
Sodium Acetate; pH 5,2	2 M	NT	R	R	NT	R
Paraffin Oil		NT	R	R	NT	R
SDS	1,0%	NT	R	R	NT	R
TRIS-Buffer; pH 7,5	1 M	NT	R	R	NT	R
Triton-X-100	1,0%	NT	R	R	NT	R
Tween 20	1,0%	NT	R	R	NT	R
Water		NT	R	R	NT	R
Xylene		NT	R	R	NT	R

R = resistant

NR = not resistant

NT = Not tested