## Microsart® @filter 100|250

### Instructions for use

Sterile, disposable filtration device for vacuum filtration consisting of a filter support with bayonet lock, a 47 mm membrane filter, a 100 or 250 ml funnel and a lid (the packaging type "in bags" is delivered without lid due to the need for stackability).

The system reduces the risk of secondary contaminations and allows fast filtration as well as high sample throughputs. No assembly, sterilization, or clean-up required. Ideal for microbiological and analytical testing of pharmaceuticals, beverages, water, and other liquids.

This product belongs to the Microsart® @vance product family – advanced colony counting from Sartorius Stedim Biotech

Sterile | Ready-for-use | Store dry at room temperature

Nitrocellulose membrane filters	
Hazard statements	Precautionary statements
H228: Flammable Solid  DANGER	P 210: Keep away from heat/sparks/open flames/hot surfaces. – No smoking. P 243: Take precautionary measures against static discharge. ADR 286; IATA A122

Producer: Sartorius Stedim Biotech GmbH  
August-Spindler-Strasse 11  
D-37079 Goettingen, Germany  
Phone: +49.551.308.0

## Instructions for use:

1. Preparation
  - a) Remove the Pouch with Tyvek® of one tray containing the Microsart® @filter units.
  - b) Remove one bag containing the stacked or individually packed Microsart® @filter units out of the box.
2. Place the tray or bag containing the Microsart® @filter units in the vicinity of the filter station.
3. Turn on the vacuum source and close the filter station plug valve.
4. Removing the Microsart® @filter unit
  - a) Open the bag at the indicated position and take out one Microsart® @filter unit.
  - b) Remove one Microsart® @filter unit from the tray by taking hold of the filter support.
5. Place the Microsart® filtration unit on the stainless steel 47 mm base support (1ZU---0002) so that the pins of the base support extend into the notches. Fix the Microsart® filtration unit into position by rotating clockwise until it engages.
6. Sampling. If you need to filter a small sample volume, it is recommended to add 10–20 ml of sterile buffer solution (e.g., 0.9% NaCl solution) beforehand.
7. Open the filter station plug valve and filter the sample.
8. Close the filter station plug valve. Rinse the inside of the funnel with 20–30 ml of sterile buffer solution.
9. Open the filter station plug valve and filter the rinsing solution.
10. Close the filter station plug valve. Remove the funnel of the Microsart® filtration unit from its filter support.
  - a) Flame the blunt-tipped tweezers (16625). Use these tweezers to remove the filter. Place the filter on a solid culture medium (e.g., agar medium) without entrapping any air bubbles.
  - b) When using Microsart @filter together with Microsart @media agar dishes, place the lid onto the Microsart @filter base and simply lift up the membrane filter to transfer it onto the agar without the usage of tweezers.
11. Incubate the filtered sample under the appropriate conditions.
12. Detach the filter support by rotating in a counter clockwise direction.

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Hersteller: Sartorius Stedim Biotech GmbH