Dräger DCD 5000

non IVD

REF 83 19 910

de	Speichelprobennehmer	Gebrauchsanweisung
en	Oral fluid sample collector	Instructions for use
fr	Collecteur d'échantillons de salive	Notice d'utilisation
es	Dispositivo de toma de muestras de saliva	Instrucciones de uso
it	Raccoglitore di saliva	Istruzioni per l'uso

 **Dräger Safety AG & Co. KGaA**
Revalstraße 1
23560 Lübeck, Germany
Tel +49 451 882 0
Fax +49 451 882 20 80
www.draeger.com

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en	Instructions for use
1	For your safety
1.1	Follow the instructions for use Any use of the Dräger DCD 5000 requires full understanding and strict observation of these instructions for use. The Dräger DCD 5000 is only to be used for the purposes specified here.
2	Description
2.1	Intended use The Dräger DCD 5000 is an oral fluid sample collector which holds approx. 0.4 mL of oral fluid to send the extracted sample to an analysis laboratory for further examination. The Dräger DCD 5000 has been designed for non-medical and forensic applications. It is suitable as a reference method for the laboratory confirmation of a preliminary drug test result which has been obtained, e.g. with the Dräger DrugTest 5000. It has been evaluated exemplary for the following drugs and their metabolites: amphetamine, benzodiazepine, delta-9-tetrahydrocannabinol, cocaine, benzoylcegonine, metamphetamine and morphine. The low variance of the collected sample volume and the proven high recovery rate of the Dräger DCD 5000 permit semi-quantitative conclusions about the analyte concentrations in the oral fluid at the time of taking the sample. If required, two Dräger DCD 5000 can be combined (Fig. 1a) to carry out simultaneous double sampling. After taking the sample, the two Dräger DCD 5000s can be separated again, packed in transport tubes and analysed separately. The Dräger DCD 5000 enables simple, hygienic sample extraction in the laboratory thanks to a patented, new ejection technology for the sample collector (Fig. 5).
2.2	Explanation of symbols
	Strictly follow the instructions for use!
	Expiration date
	Temperature limitation
	Contents sufficient for <n> tests
	Manufacturer

2.3 Scope of delivery

- 20 oral fluid sample collectors in transport tubes with plugs
- Instructions for use

2.4 Recommended additional material

- Latex gloves
- Methanol and/or methanolic solution for storage
- Extraction cylinder with screw cap for methanolic solution

For technical support, please contact your local dealer, Dräger or Technical Support.

3 Precautions and warnings

	NOTICE For hygienic reasons, always use gloves when handling the Dräger DCD 5000 and any samples. Do not touch the sample collector with bare hands before and after sampling and follow the current hygiene rules.
---	---

- Follow handling and disposal procedures correctly.
- In order to prevent contamination of the sample collector, only open the transport tube immediately before use.
- If the expiration date on the packaging has expired, do not use the Dräger DCD 5000. The expiration date is formatted as YYYY-MM, e.g. 2019-01 means the Dräger DCD 5000 must not be used after the end of January 2019.

4 Testing

4.1 Preparation

- Ensure that the test person has had no substances like, e.g. food, chewing gum or tobacco in the mouth for at least 10 minutes before sampling.

4.2 Sampling

4.2.1 Simultaneous collection of two oral fluid samples (Fig. 1a)

1. Remove the two oral fluid sample collectors from the transport tubes and join them at the ribbed handle piece to form a "twin sample collector". Further handling is the same as for taking a single oral fluid sample. Separate both oral fluid sample collectors after the sampling process.

4.2.2 Collection of a single oral fluid sample (Fig. 1b):

1. Remove the oral fluid sample collector from the transport tube and hand it to the test person.
2. Instruct the test person to place the oral fluid sample collector inside the cheek and move it carefully from one side of the mouth to the other (**Fig. 2**). The test person must not chew or suck on the oral fluid sample collector. Observe the sampling process!
3. A sufficient sample will be collected within 4 minutes. The sampling process can be terminated if the colour indicator of the sample collector turns blue within the 4 minutes (**Fig. 3**).
4. Take the oral fluid sample collector from the test person.
5. Insert the oral fluid sample collector with the mouthpiece first back into the open transport tube (**Fig. 4**).
6. Eject the mouthpiece (**Fig. 5**).
7. Seal the transport tube with the sealing plug and mark the sample unequivocally on the tube.

5 Storage

5.1 Storage of collected oral fluid samples

Basically, the oral fluid samples taken should be analysed in the laboratory as quickly as possible. Studies have shown that analyte loss from the absorption of the drug at vessel walls is possible if oral fluid samples are collected and stored directly in plastic vessels without the use of a suitable oral fluid sample collector.

For the Dräger DCD 5000 recovery rates have been determined for all common classes of drugs.

After storage at room temperature (15 - 30 °C) over a period of 7 days the recovery rate for delta-9-tetrahydrocannabinol is approx. 7%, for amphetamine, benzoylcegonine with cocaine, metamphetamine, morphine and metabolite approx. 60 - 80 %, for nordazepam, oxazepam, temazepam, lorazepam, bromazepam and flurazepam >80 %, for diazepam approx. 30 % and midazolam approx. 50%.

After immediate storage in a methanol solution over a period of 7 days at room temperature (15 - 30 °C) the recovery rate for delta-9-tetrahydrocannabinol is approx. 75 %, for benzoylcegonine with cocaine approx. 70%, and for all other drugs mentioned >90 %.

5.2 Storing an unused Dräger DCD 5000

Always store the Dräger DCD 5000 at 4 to 30 °C (39.2 – 86 °F) in its original packaging.


6 Despatch

If required, stick a security seal over the sealing plug, fill in the order form for the analysis laboratory and pack the transport tubes for despatch. Despatch and/or storage must be in accordance with local regulations (e.g. by using special envelopes).

Send samples to an accredited laboratory. The required address labels, order forms and suitable packaging are available from the representative of the laboratory of your choice.

7 Disposal

Dispose of the Dräger DCD 5000 in accordance with valid regulations.

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© Dräger Safety AG & Co. KGaA