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Cutting ring saw with blade for cutting gold or silver rings. Blade to cut the steel rings delivered separately.

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Store in a clean and dry place.

➤ Nettoyage, pré-désinfection / Cleaning, pre-disinfection

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Standard conditions for reusable surgical instruments.

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*The blades of cutting rings saw are not compatible with machine cleaning due to **CARBON STEEL** used to provide cutting quality for different types of rings.*

The applied cleaning method should be compatible with carbon steel and/or chrome steel instruments.

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The Sterilization shall be performed according to ISO 17665-1 validated steam sterilisation procedure (fractional vacuum procedure) in an autoclave in compliance with EN 285.

Saturated steam sterilization method	Temperature	Time	Drying time	Absolute pressure
Pre-vacuum	132 °C (270 °F)	4 minutes	30 minutes	2.9 bar, 2.8 atm, 41.6 psi
Pre-vacuum	134 °C (273 °F)	3 minutes	35 minutes	3.0 bar, 3.0 atm, 44.1 psi
Pre-vacuum	134 °C (273 °F)	5 minutes	35 minutes	3.0 bar, 3.0 atm, 44.1 psi
Pre-vacuum	134 °C (273 °F)	18 minutes	35 minutes	3.0 bar, 3.0 atm, 44.1 psi

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The blades of cutting rings saw are not compatible with autoclave sterilization, due to CARBON STEEL used to provide cutting quality for different types of rings.

The applied disinfection method should be compatible with carbon steel and/or chrome steel instruments.

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