Natus Snap Leads

Instructions for Use:

Natus Medical Incorporated
DBA Excel-Tech Ltd. (XLTEK)
2568 Bristol Circle
Oakville, Ontario, L6H 5S1 Canada
Tel: +1 905-829-5300
Website: natus.com

EU Representative / Importer
Natus Manufacturing Limited
IDA Business Park
Gort, Co. Galway, Ireland

Rx Only

Associated product part numbers:

1411517, 1411524, 1411520, 1411535, 1411530, 1411516, 1411541, 1411548 and 1411542

031530 Rev C 12/2020

English

Description:

Snap Leads connect snap electrodes with an amplifier for monitoring electromyogenic and physiologic data signals commonly recorded from chest, extremities and faces of patients.

Intended Use:

The reusable 1.5 mm button snap cables are intended to be used with compatible disposable snap electrodes for transferring surface biopotential signals to use during sleep studies, EEG, EMG, and recording evoked potential signals from a cutaneous location. For use on adults.

Intended User and Target Patient Group:

The product is for use only under the direction of trained medical professionals in acquiring electrophysiologic signals.

Clinical Benefits:

Facilitates obtaining the recording of an EEG/PSG study to detect any irregularities indicative of various brain or sleep disorders.

Contraindications and Side Effects:

There are no known contraindications or side effects for procedures performed with Natus reusable Snap Leads. See Warnings or Precautions below.

Operating Instructions:

- Connect the appropriate snap lead to the snap connector on the electrode.
- Reusable Snap Leads should be cleaned between uses. Cleaning instructions are below.

Cleaning Instructions:

- Clean with a commercial wipe such as CaviWipes™ or Sani-Cloth® to remove visible soil.
- Wipe the article using a lint-free cloth and air dry.

Understanding Cautions Statements:

 CAUTION
Refers to a hazardous situation that could result in minor or moderate injury or material damage if not avoided.
• Information on how the hazardous situation is avoided.

Warnings and Precautions:

 CAUTION
Device dropped or damaged in transit/use could lead to loss of function or delayed diagnosis.
• Inspect the device prior to each use and do not use if damaged.
Data cannot be recorded due to improper placement leading to delay in procedure.
• Assure proper connections and check signal quality before using.
Modifications to the device can affect function and performance.
• Do not modify this equipment without authorization of the manufacturer.

Environmental Specifications:

Operating Conditions:

- Temperature: +10°C (+50°F) to +30°C (+86°F)
- Relative Humidity: 20% to 80%
- Pressure: 70 kPa to 106 kPa

Storage Conditions:

- Temperature: -25°C (-13°F) to +55°C (+131°F)
- Relative Humidity: 20% to 80%
- Pressure: 50 kPa to 106 kPa

Disposal Instructions:

Natus is committed to meeting the requirements of the European Union WEEE (Waste Electrical and Electronic Equipment) Regulations 2014. These regulations state that electrical and electronic waste must be separately collected for the proper treatment and recovery to ensure that WEEE is reused or recycled safely. In line with that commitment Natus may pass along the obligation for take back and recycling to the end user, unless other arrangements have been made. Please contact us for details on the collection and recovery systems available to you in your region at natus.com

Electrical and Electronic Equipment (EEE) contains materials, components and substances that may be hazardous and present a risk to human health and the environment when WEEE is not handled correctly. Therefore, end users also have a role to play in ensuring that WEEE is reused and recycled safely. Users of electrical and electronic equipment must not discard WEEE together with other wastes. Users must use the municipal collection schemes or the producer/importers take-back obligation or licensed waste carriers to reduce adverse environmental impacts in connection with disposal of waste electrical and electronic equipment and to increase opportunities for reuse, recycling and recovery of waste electrical and electronic equipment.

Equipment marked with the below crossed-out wheeled bin is electrical and electronic equipment. The crossed-out wheeled bin symbol indicates that waste electrical and electronic equipment should not be discarded together with unseparated waste but must be collected separately.

Disclaimer:

Natus Medical Incorporated DBA Excel-Tech Ltd. (Xitek) is not responsible for injury, infection or other damage resulting from the use of this product. Any serious incident that has occurred in relation to the device should be reported to Natus Medical Incorporated DBA Excel-Tech Ltd. (Xitek) and the competent authority of the Member State in which the user and/or patient is established. Refer to the Natus website for an electronic copy of this document.

Glossary of Symbols

Symbol	Standards Reference	Standard Title of Symbol	Symbol Title as per Reference Standard	Explanation
Medical Device	-	-	An indication of Medical device	This product is a medical device.
Rx Only	21 CFR Part 801.109(b)(1)	Labeling-Prescription devices.	Prescription only	Indicates the product is authorized for sale by or on the order of a licensed healthcare practitioner.
	ISO 15223-1 Symbol 5.4.5 (Reference Annex B for the general prohibition symbol)	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Not made with Natural Rubber Latex	Indicates a medical device that is not made with natural rubber latex.
	2012/19/EU	Waste Electrical and Electronic Equipment (WEEE)	Disposal at end of operating life instructions	Indicates that electrical and electronic equipment waste should not be discarded together with unseparated waste but must be collected separately.
	ISO 15223-1 Symbol 5.1.1	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Manufacturer	Indicates the medical device manufacturer.
	ISO 15223-1 Symbol 5.1.2	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Authorized representative in the European Community	Indicates the Authorized representative in the European Community.

Symbol	Standards Reference	Standard Title of Symbol	Symbol Title as per Reference Standard	Explanation
	ISO 15223-1 Symbol 5.1.3	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Date of manufacture	Indicates the date when the medical device was manufactured.
	ISO 15223-1 Symbol 5.1.5	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Batch or Lot code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	ISO 15223-1 Symbol 5.1.6	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	ISO 15223-1 Symbol 5.4.3	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Consult instructions for use	Indicates the need for the user to consult the instructions for use.
	ISO 60601-1 Table D.1 #11	Medical electrical equipment — Part 1: General requirements for basic safety and essential performance.	Operating instructions	
	ISO 60601-1 Table D.2 #10	Medical electrical equipment — Part 1: General requirements for basic safety and essential performance.	Follow instructions for use	Refer to instruction manual/ Booklet. NOTE on ME EQUIPMENT "Follow instructions for use"
	ISO 15223-1 Symbol 5.4.4	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Caution: Read all warnings and precautions in instructions for use	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	ISO 60601-1 Table D.1 #10	Medical electrical equipment — Part 1: General requirements for basic safety and essential performance.		
	ISO 60601-1 Table D.2 #2	Medical electrical equipment — Part 1: General requirements for basic safety and Essential performance.	General warning sign	Indicates a hazard of potential personal injury to patient or operator.
	ISO 15223-1 Symbol 5.2.8	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
	ISO 15223-1 Symbol 5.3.7	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Temperature limit	Indicates the (storage) temperature limits to which the medical device can be safely exposed.
	ISO 15223-1 Symbol 5.3.8	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied.	Humidity limitation	Indicates the range of (storage) humidity to which the medical device can be safely exposed.

Product Part Numbers	Product Description
1411517	1.0m (39"), (1.5mm DIN) Orange, 2/pkg
1411524	2.0m (78"), (1.5mm DIN) Orange, 2/pkg
1411520	1.5m (60"), (1.5mm DIN) Black, 2/pkg
1411535	3.0m (120"), (1.5mm DIN) Black, 2/pkg
1411530	3.0m (120"), (1.5mm DIN) White, 2/pkg
1411516	1.0m (39"), (1.5mm DIN) Green, 1/pkg
1411541	0.6m (24"), (1.5mm DIN) Red, 1/pkg
1411548	0.6m (24"), (1.5mm DIN) Green, 1/pkg
1411542	0.75m (29.5"), (1.5mm DIN) Black, 1/pkg

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