
EN	Electronic Filtration Module (EFM) for Snore Sensor
FR	Electronic Filtration Module (Module filtration électronique - EFM) pour capteur de ronflements
DA	Electronic Filtration Module (elektronisk filtreringsmodul - EFM) til snorkesensor
DU	Electronic Filtration Module (elektronische filtratiemodule) (EFM) voor snurksensor
GER	Electronic Filtration Module (Elektronisches Filtrationsmodul) (EFM) für den Schnarchsens
IT	Electronic Filtration Module (Modulo di filtraggio elettronico - EFM) per sensore del russamento
NO	Electronic Filtration Module (EFM- Elektronisk filtreringsmodul) for snorkesensor
PT	Electronic Filtration Module (Módulo de Filtragem Eletrônico - EFM) para Sensores de Ronco
SP	Electronic Filtration Module (EFM, módulo de filtración electrónica) para sensor de ronquidos

86-41220-0040 Rev D

U.S. Patent US6702755

EN

CAUTION STATEMENTS

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician. This product should be used only under the direct supervision of a trained medical professional. This product is for diagnostic purposes only and is not to be used in life supporting or life sustaining situations.

FEATURES AND INTENDED USE:

The Dymedix Electronic Filtration Module (EFM) has been designed to connect the Dymedix Snore Sensor to the headbox during sleep disorder studies for recording onto a physiological recorder. Features and specifications of this device are subject to change without notice.

INSTRUMENT CONNECTIONS & SETTINGS

Refer to the polysomnograph operator's manual for direction in determining proper instrument connection and settings. Response is affected by variables such as sensitivity settings, high and low frequency filter settings, sensor placement, patient movement, etc.

PRECAUTIONS FOR CABLE

1. Store EFM at room temperature.
2. Do not sterilize EFM.
3. If product is damaged, discontinue use.
4. Connect the EFM to electrically isolated equipment.
5. Use all available recorded parameters in making diagnosis.
6. Non-accessible Lithium Battery cannot be recharged or replaced. The battery is recyclable and shall be disposed of in accordance with applicable local, state/provincial or federal/national regulations. Do not incinerate.

PATIENT INSTALLATION

1. Connect the Dymedix sensor to the two receiving connections on the EFM.
2. Connect the EFM to the headbox.
3. The EFM must be connected to an AC-coupled channel on the recording instrument.
4. Proceed with polysomnography as outlined by clinician procedure.
5. **BE SURE THE PATIENT IS PROPERLY GROUNDED.**

REMOVAL INSTRUCTIONS

1. Disconnect EFM leads from recording equipment.
2. The EFM may be cleaned with a hospital-approved disinfectant.
3. Do not immerse in fluids.

WARRANTY

Dymedix warrants its EFM to be free from defects in materials and workmanship and to perform satisfactorily for one year from the date of manufacture. Should the EFM fail, the sole liability of Dymedix is limited to replacement of the product with no charge for parts or labor. Under no circumstances shall Dymedix be liable for any consequential damages arising from the use of, or the inability to use this product. This warranty is rendered void for conditions resulting from: damage, marginal performance or malfunctions caused by: misuse, abuse, neglect, improper amplification settings, improper line voltage, power fluctuations, or any adverse environmental conditions, tampering, unauthorized modifications, adjustments or repairs to the product or its accessories. This warranty is in lieu of all other warranties, expressed or implied and is extended only to the original purchaser.

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Dymedix Diagnostics Inc.