Instructions for Use—English

Model 6500MA Single-Patient Use, Disposable Pulse Oximeter Sensor

Indications for Use

Nonin's Model 6500MA Single-Patient Use, Disposable Pulse Oximeter Sensor is indicated for non-invasive spot-checking and/or continuous monitoring of adult and pediatric patients, who are well or poorly perfused, weighing more than 60 pounds (27 kilograms). It is intended for use in environments including the operating room, surgical recovery, critical care, emergency room, long-term care, home use and medical environments.

⚠️ CAUTION: Federal law (USA) restricts this device to sale by or on the order of a licensed practitioner.

Warnings:

- Do not use the device in an MR environment, in an explosive atmosphere, or on infant or neonatal patients.
- This device is not defibrillation proof per IEC 60601-1.
- The use of sensor and oximeter combinations other than Nonin-branded components have not been tested for accuracy as a system and may affect performance of the system. Refer to Nonin pulse oximeter operator's manuals for a complete listing of Nonin-branded oximeters, sensors, and accessories.
- Inspect the sensor application site at least every 4 hours to ensure correct sensor alignment and skin integrity. Patient sensitivity to sensors may vary due to medical status and skin condition.
- Avoid excessive pressure to the sensor application site as this may cause damage to the skin beneath the sensor.

⚠️ Cautions:

- Do not use a damaged sensor. If the sensor is damaged, discontinue use immediately.
- Do not sterilize, autoclave, or immerse in liquid of any kind.
- Do not use caustic or abrasive cleaning agents on the sensor.
- Follow local governing ordinances or recycling instructions regarding disposal or recycling of the sensor and any components.
- A functional tester cannot be used to assess the accuracy of a pulse oximeter monitor or sensor.
- Refer to the pulse oximeter operator's manual for additional warnings and cautions.
- Factors that may degrade pulse oximeter performance include the following:
 - excessive ambient light
 - venous pulsations
 - excessive motion
 - excessive interference
 - anemia or low hemoglobin/concentrations
 - cardiovascular dyes
 - carboxyhemoglobin
 - methemoglobin
 - incorrect sensor type
 - poor pulse quality

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Symbols:

Symbol	Definition	Symbol	Definition
	Follow Instructions for Use		Lot Number
	CAUTION!		Catalogue number
	Do Not Reuse		Quantity
	CE Marking indicating conformance to CE Directive No. 93/42/EEC concerning medical devices		Date of manufacture
	Manufacturer		Country of manufacture
	Protected against vertically falling water when enclosure is tilted up to 15 degrees and ingress of solid foreign objects greater than equal to 2.5 mm in diameter per IEC 60529.		Storage/Shipping temperature range
	Authorized representative in the European Community		Indicates separate collection for waste electrical and electronic equipment (WEEE)
			Medical prescription required

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
Applying the Single-Patient Use Sensor

- Notes:
 - Proper sensor placement is critical for good performance. If the sensor is not positioned properly, light may bypass the tissue and result in SpO₂ inaccuracies.
 - The fingertip should not be inserted through the center fold opening.
 - The sensor cable may be secured with medical tape.

Application by a Healthcare Professional

- Place the patient's fingertip at the center fold, so the cable is on top of the finger and the black foam touches the skin.
- Bend sensor around the finger to conform to finger.
- Gently pinch just the sides of the sensor to secure comfortably. Do not squeeze the sensor as it may interfere with its performance.

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