⌘

Single Use Only

⚠

Attention! Consult accompanying documentation

R

ONLY

Directions for Use

CE

2797

EC

REP

Natus Manufacturing Limited
IDA Business Park, Gort
Co. Galway, Ireland
Tel: +353 (0)91 647400
Fax: +353 (0)91 630050

🏭

Natus Medical Incorporated
DBA Excel-Tech Ltd. (XLTEK)
2568 Bristol Circle Oakville,
Ontario L6H 5S1 Canada
Tel : 613-254-8877
Fax : 716-200-1091
1-888-662-7632
www.natus.com

DESCRIPTION

The BreathSensor airflow thermistor is indicated for single use only in sleep studies where breath signals are measured and recorded. The Model 970 series has four airflow thermistors for different patient sizes:

- Model 971 Adults
Model 974 Children and Small Adults
Model 976 Infants
Model 978 Premies

The Airflow Thermistors are for use with standard polysomnograph equipment.

The Airflow Thermistors are contraindicated for use on patients who exhibit allergic reactions to the adhesive tape.

INDICATIONS FOR USE


The Embla BreathSensor provides a qualitative air flow signal by oral and/or nasal temperature sensitive resistive components for recording onto a recording system in support of airflow analysis and sleep studies.

The Embla BreathSensor does not provide any diagnostic conclusion about the patient’s condition to the user.

The Embla BreathSensor user is a qualified medical practitioner who will exercise their professional judgment in using this information.


INSTRUCTIONS FOR USE

1. Select the appropriate size thermistor for the patient. Use patient size, not age, as a guide to selecting the correct model. To fit properly, the distance between the nasal sensing elements should match the distance between the patient’s nostrils.
2. Ensure that the patient’s skin is clean, dry, and free of oil or lotion.
3. Remove the airflow thermistor from the packaging. Remove the backing from the airflow thermistor and tape assembly. Figure ° shows the BreathSensor.

Warning: Make sure the nasal sensing elements do not block nasal airflow.

4. Position the airflow thermistor so that the two nasal sensing elements are aligned under the patient’s nostrils, as in Figure °. The oral sensing element, if present, should be positioned between the lips. (Model 978 does not have an oral sensing element.)

Warning: Discontinue use if skin irritation occurs.

5. Secure the thermistor by pressing the adhesive tape against the patient’s skin. Make sure that the sensing elements do not touch the patient’s skin or lip.
6. Insert the folded plastic tab into the airflow cable connector. See Figure °. Gently push the folded tab in as far as it will go. The tab fits snugly inside the connector.

Note: Embla BreathSensor Models 976 and 978 have a flat, unfolded tab with circuitry printed on both sides of it. This tab must remain flat when inserted. Do not fold the tab. See Figure °.

7. If necessary, wrap a piece of tape around the thermistor attachment and connector to avoid accidental detachment.

Caution: Use of this thermistor while the patient is in a heated incubator may result in reduced signal amplitude, and possibly no signal.

8. Make sure that the BreathSensor airflow thermistor is completely unobstructed. Any foreign object, such as a pacifier, resting on the thermistor will adversely affect its performance. Any obstruction, like a pH probe, in one or both nostrils will reduce the signal amplitude created by the BreathSensor airflow thermistor.

REMOVAL

Remove the BreathSensor airow thermistor by gently peeling it away from the patient’s skin.

RECOMENDED FILTER SETTINGS:

- High Frequency Filter shall be capable of 15 Hz

• Low Frequency Filter shall be capable of 0.1 Hz
- If Time Constant is used the TC shall be capable of 1.6 seconds

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Copyright © 2020 Natus. Alle Rechte vorbehalten.
Copyright © 2020 Natus. Todos los derechos reservados.