DIRECTIONS FOR USE

Disposable Subdermal Needle Electrodes
Disposable Monopolar Needle Electrodes

Safe removal of protection sheath: twist and pull.

Manufactured by:
Technomed Europe
Amerikalaan 71
6199 AE Maastricht-Airport
The Netherlands

CE
0344

Do Not Reuse

USA Only
Rx Only

STERILE EO

Doc. 27056-001A

ENGLISH

DISPOSABLE NEEDLE ELECTRODES

Indications for Use

Needle Electrodes for Neurological purposes are intended for use with recording/monitoring equipment for the recording of biopotential signals including electroencephalography (EEG), electromyography (EMG) and nerve potential signals, and are intended for stimulation/recording with stimulation/recording equipment for electromyography (EMG) and nerve potential signals.

Directions for use

When selecting a disposable needle electrode, it is important to use one that is of suitable length and diameter for the investigation being undertaken.

CAUTION!

- The disposable needle electrodes are sterile unless packaging is damaged or opened or use before date is passed.
- When introducing a disposable needle electrode into a patient it is recommended that it is not inserted up to the hub.
- If the disposable needle electrode bends before, during, or after insertion, do not straighten or re-insert. Bent disposable needle electrodes should be discarded and replaced with a new disposable needle electrode.
- Always discard used disposable needle electrodes in a properly marked medical biohazard sharps container.

WARNING!

If used as recording electrodes, the possible loss of electrocautery return path may result in high current density passing through these needles with possible tissue heating and damage resulting. If used as stimulation electrodes, thermal burns could result with prolonged stimulation outputs above 0.25W/cm².

For single patient use only. Do not reuse or re-sterilize. Cleaning and re-sterilization of the device can affect the safety, performance and effectiveness of the device and expose patients and users to unnecessary risk such as infection and transmissible disease.

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