Reusable Flow

Sensors

Item NA-N14XX – ver. 1.0

S.L.P. Ltd.
62 Anshewitz St.
Israel
Tel.: +972 3 6870640
Fax: +972 3 5371282
E-mail: info@sleepsense.com

S.L.P. Inc.
2511 Technology Dr., Suite 112
Elgin, IL 60124, USA
Tel.: +1 (630) 513 7479
Fax: +1 (630) 513 8478
Tel/Fax: +44 19328 59810

Medes Limited
5 Beaumont Gate,
Shenley Hill, Radlett,
Herts WD7 7AR, England
Tel.: +423-663-169205
Tel/Fax: +44 19328 59810

Manufactured by S.L.P. Ltd. for Natus Neurology Incorporated. Made in Israel.

Specifications are subject to change without notice
Les spécifications sont sujettes à modification sans préavis
Especificaciones sujetas a cambio sin previo aviso
Le specificazioni sono soggette a modifica senza preavviso
Die Spezifikationen sind ohne vorherige Meldung Änderungen unterworfen

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


Thank you for purchasing a SleepSense® sensor.

Description

Comprising three or four heat sensitive elements, the respiration flow sensor is placed under the patient’s nostrils and produces a signal that is directly proportional to the temperature changes of air inhaled and exhaled during respiration. The sensor generates a small analog electrical voltage that provides a clear, reliable indication of respiration airflow. The sensors should be connected to the system’s thermal-flow input or, when using the interface, to any AC (low level) channel on the recorder.

Indications for use

SleepSense sleep-lab sensors provide a qualitative measure of sounds, respiratory-effort, flow, body position or limb movements, for recording on an approved data acquisition system. They are intended for use on children and adult patients who are screened during sleep disorder studies at a sleep laboratory or the patient’s home.

Applying the Sensor

The respiration flow sensor is placed on the patient's upper lip, under the nostrils.
The sensors are delicate and must be handled with care.

Positioning the Sensor

1. Bend the two nasal prongs forward to prevent them from entering the nostrils. Make a round bend slightly away from the sensor's body.
2. Align the tip of the oral prong with the center of the mouth/upper lip.
3. Draw the wires on either side of the sensor over the patient’s ears and use the clasp under the chin to secure the wires in place.

Recommendations:

- Tape the wire on both cheeks to reduce strain and enhance reliability.
- Check that clear, strong signals are being transmitted before leaving the patient.

Cleaning the Sensor

- Wipe the sensor and cable with a non-corrosive (to plastic) cleanser to clean before use.
- Make sure the complete sensor assembly is thoroughly dry before reusing it.
- Do not soak the sensor in disinfectants.

Technical Specifications:

Sensor:	Thermocouple		Thermistor
Technology:	Small rounded TPE center with 3 prongs		
Recommended filter settings:	High pass: 0.1 Hz / Low pass: 70 Hz or as recommended in the system's instructions		
Signal output:	Approximately 0.2 mV for normal respiration	With Interface- App. 1mV	W/O Interface - NTC Ω change
Operating conditions:	5°C (40°F) - 40°C (104°F)		
Storage temperature:	-20°C (-4°F) - 60°C (140°F)		
Operating and storage humidity:	5% - 95% (Non-condensing)		

Warnings and Precautions

- **Caution:** US Federal Law restricts this device to sale, distribution and use by or on the order of a physician.
- SleepSense sensors are for professional use only.
- SleepSense sensors may be used only in conjunction with an approved recording system.
- SleepSense sensors are only intended for use on healthy skin.
- Discontinue use if sensor shows signs of wear, damage, or exposed metal.

Limited Warranty
SLP warrants the respiration-flow sensor to be free of defects in materials and workmanship for a period of six months from the date purchased. The sole liability of SLP and our distributor(s) is limited to replacement or repair of the product at the option of SLP with no charge for parts or labor if any part is proven to be defective in workmanship, performance, or materials during the warranty period. Under no circumstances shall SLP or our distributor(s) be liable for any loss of revenues or damage, direct, consequential, or incidental, including loss of profit, property damage, or personal injury arising from the use of, or the inability to use this product. This warranty is intended only for the original buyer and is in lieu of all other warranties or previous agreements, expressed or implied. This warranty is rendered void if the product is used for other than its intended purpose or is subject to abuse, misuse, tampering, neglect, or unauthorized modifications. Use of this product constitutes acceptance of this warranty in total.

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