EN

Interface Cable for Airflow Sensor

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CAUTION STATEMENTS

!

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

This product should be used only under the direct supervision of a trained medical professional.

This product is for diagnostic purposes only and is not to be used in life supporting or life sustaining situations.

FEATURES AND INTENDED USE:

The Dymedix Interface Cable connects the Airflow Sensor to the headbox during sleep disorder studies for recording onto a physiological recorder. Features and specifications of this device are subject to change without notice.

This Interface Cable must be used with a Dymedix Airflow sensor.

INSTRUMENT CONNECTIONS & SETTINGS

Set the filter settings on polysomnograph as close to Interface Cable specifications as allowed by manufacturer on PSG equipment. Refer to the polysomnogram operator’s manual for direction in determining proper instrument connection and settings. Response is affected by variables such as sensitivity settings, high and low frequency filter settings, sensor placement, patient movement, battery life etc. Incorrect connection could result in erratic or no signal.

PRECAUTIONS

1. Store Interface Cable at room temperature.

2. Do not sterilize Interface Cable.

3. If product is damaged, discontinue use.

4. Connect the Interface Cable to electrically isolated equipment.

5. Use all available recorded parameters in making diagnosis.

RED RECEIVING CONNECTORS

TO HEADBOX AIRFLOW INPUT

AIRFLOW SENSOR

PATIENT INSTALLATION

1. Connect the Dymedix airflow sensor to the two receiving connections on the Interface Cable matching the color coded ends.

2. Connect the Interface Cable to the headbox of the polysomnograph.

3. The Interface Cable must be connected to an AC-coupled channel on the recording instrument. Proceed with polysomnography as outlined by clinician procedure.

4. BE SURE PATIENT IS PROPERLY GROUNDED.

REMOVAL INSTRUCTIONS

1. Disconnect Interface Cable leads from recording equipment.

2. The Interface Cable may be cleaned with a non-corrosive hospital-approved disinfectant.

3. Do not immerse in liquids.

WARRANTY

Dymedix warrants its Interface Cable to be free from defects in materials and workmanship and to perform satisfactorily for one year from the date of manufacture. Should the Interface Cable fail, the sole liability of Dymedix is limited to replacement of the product with no charge for parts or labor. Under no circumstances shall Dymedix be liable for any consequential damages arising from the use of, or the inability to use this product. This warranty is rendered void for conditions resulting from: damage, marginal performance or malfunctions caused by: reuse, misuse, abuse, neglect, improper amplification settings, improper line voltage, power fluctuations, or any adverse environmental conditions, tampering, unauthorized modifications, adjustments or repairs to the product or its accessories. This warranty is in lieu of all other warranties, expressed or implied and is extended only to the original purchaser.

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Dymedix Diagnostics Inc.