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Fax: +972 3 6870640

Tel: +972 3 5371282

Manufactured by S.L.P., Ltd. for Natus Neurology Incorporated. Made in Israel.

To Reorder: Call: 1-800-356-0007 (in the USA) or your Natus Neurology Incorporated distributor

or email: supplies.nc@Natus.com

Item NA-N148XX - ver. 2.0 [Side A]

User Manual

Respiration Pressure Flow Sensors

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Thank you for purchasing a SleepSense® sensor.

### Description

This sensor is designed to detect respiratory airflow and snoring for recording onto a polygraph, EEG or similar type of physiological recorder. The internal Piezo Crystal senses minute pressure changes and filters them into two bands: one for slow waves to the flow channel, and one for fast waves to the snore channel. The sensor plugs directly into the electrode jack box, (no external power source is required), and a nasal cannula is attached to the patient.

### Indications for use

SleepSense sleep-lab sensors provide a qualitative measure of sounds, respiratory-effort, flow, body position or limb movements, for recording on an approved data acquisition system. They are intended for use on children and adult patients who are screened during sleep disorder studies at a sleep laboratory or the patient's home.

### Applying the Sensor

The detection of airflow is done by attaching either an ETCO2 or a standard O2 cannula to the sensor and the patient. The airflow sensed through the nose will then provide both a nasal airflow signal and a snore signal. Use a bacterial filter between the cannula and sensor. There are two sets of cables for connection to the recorder. One is for the airflow channel and one is for the snore channel. The two black connectors are connected internally in the sensor. The sensor can be connected to a nasal CPAP mask using standard respiratory tubing instead of using a cannula. When connecting an O2 cannula to the sensor always use the Luer® connector to connect to the sensor box. The ETCO2 cannula should be connected directly to the sensor box. The recorder's High Pass filter should be set at the lowest possible setting in order to maximize the time constant. **The 14831 sensor has two sensitivity-setting switches on the side panels. You can adjust these switches to increase or decrease the sensitivity of the snore or pressure channels if the signals are too strong or too weak for proper recording.**

### Cleaning the Sensor

- Wipe the sensor and cable with a non-corrosive (to plastic) cleanser to clean before use.
- Make sure the complete sensor assembly is thoroughly dry before reusing it.
- Do not soak the sensor in disinfectants.

### Technical Specifications:

|                                        |                                                                                                                         |
|----------------------------------------|-------------------------------------------------------------------------------------------------------------------------|
| <b>Recommended filter settings:</b>    | Flow: High pass: 0.01 Hz / Low pass: 70 Hz<br>Snore: High pass: 30 Hz<br>or as recommended in the system's instructions |
| <b>Signal output:</b>                  | Approximately 0.5 mV for normal respiration<br>or loud snoring                                                          |
| <b>Operating conditions:</b>           | 5°C (40°F) - 40°C (104°F)                                                                                               |
| <b>Storage temperature:</b>            | -20°C (-4°F) - 60°C (140°F)                                                                                             |
| <b>Operating and storage humidity:</b> | 5% - 95% (Non-condensing)                                                                                               |

Warnings and Precautions

- **Caution:** US Federal Law restricts this device to sale, distribution and use by or on the order of a physician.
- SleepSense sensors are for professional use only.
- SleepSense sensors may be used only in conjunction with an approved recording system.
- SleepSense sensors are only intended for use on healthy skin.
- Discontinue use if sensor shows signs of wear, damage, or exposed metal.

**Attention:** certain sensors contain a lithium battery- dispose of properly.

Limited Warranty  
SLP warrants the pressure flow sensor to be free of defects in materials and workmanship for a period of twelve months from the date purchased. The sole liability of SLP and our Distributor(s) is limited to replacement or repair of the product at the option of SLP with no charge for parts or labor if any part is proven to be defective in workmanship, performance, or materials during the warranty period. Under no circumstances shall SLP or our Distributor(s) be liable for any loss of revenues or damage, direct, consequential, or incidental, including loss of profit, property damage, or personal injury arising from the use of, or the inability to use this product. This Warranty is intended only for the original buyer and is in lieu of all other warranties or previous agreements, expressed or implied. This warranty is rendered void if the product is used for other than its intended purpose or is subject to abuse, misuse, tampering, neglect, or unauthorized modifications. Use of this product constitutes acceptance of this warranty in total.

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