
GE Healthcare

Airway Adapter, Disposable

en	Instructions For Use
cs	Návod k použití
da	Instruktioner
de	Gebrauchsanweisung
et	Kasutusjuhend
el	Οδηγίες χρήσης
es	Instrucciones de uso
fr	Mode d'emploi
ko	사용 지침
hr	Upute za upotrebu
it	Istruzioni per l'uso
lv	Lietošanas instrukcija
lt	Naudojimo instrukcija
hu	Használati útmutató
nl	Gebruiksaanwijzing
ja	取扱説明書
no	Bruksanvisning
pl	Instrukcja użycia
pt	Instruções de utilização
pt-br	Instruções de uso
ro	Instrucțiuni de utilizare
ru	Инструкция по применению
sk	Návod na použitie
sr	Uputstvo za upotrebu
fi	Käyttöohjeet
sv	Bruksanvisning
tr	Kullanım Talimatları
zh	使用说明

Type: Adult/Pediatric

REF 3090121002

Indications for Use

The Airway Adapter is intended for continuous monitoring of expired and inspired respiratory gases. The Airway Adapter is a single patient use device with a cumulative duration of use of up to 72 hours.

The Airway Adapter is for use with compatible GE Healthcare sidestream respiratory gas modules and patient monitors utilizing these modules.

The Airway Adapter is provided in a large size for use with adult/pediatric patients, for connection to the patient's breathing circuit, and may be used with separately supplied gas sample lines. The Airway Adapter is intended for professional use environments and applications specified for the connected respiratory gas modules and patient monitors.

The Airway Adapter is intended for use by qualified medical personnel only.

For prescription use only.

Warnings

Warning: Due to risk of cross contamination, do not re-use.

Warning: Loose or damaged adapter or connections may compromise ventilation or cause an inaccurate measurement of respiratory gases. Securely connect all components and check connections for leaks according to standard clinical procedures.

Warning: Do not use the Airway Adapter with infants as the adapter adds 6.6 ml dead space to the patient circuit.

Warning: Do not place the Airway Adapter (tee) between the suction catheter and endotracheal tube. This is to ensure that the Airway Adapter does not interfere with the functioning of the suction catheter.

Cautions

Caution: Do not attempt to clean, disinfect, sterilize or rinse the Airway Adapter with liquids.

Caution: Dispose of Airway Adapter in accordance with standard operating procedures or local regulations for the disposal of contaminated medical waste.

Notes

Note: For accuracy specifications and further instructions, warnings and cautions, please consult the gas monitor instructions for use.

Instructions for Use


Figure 1


Figure 2

1. Use the Adult/Pediatric size with endotracheal tubes of inner diameter > 4 mm
2. Connect the Airway Adapter to the proximal end of the breathing circuit between the endotracheal tube and the ventilator Y-section Airway Adapter facing upwards as shown in Figure 1.
3. Connect a separate gas sampling line (male Luer connector) to the Airway Adapter (female Luer connector) as shown in Figure 1. Rotate the Luer connector clockwise for a secure fit to ensure that there is no leakage of gases at the connection point.
4. Connect a separate gas sampling line (male Luer connector) to the respiratory gas monitor as shown in Figure 2.
5. After connection of the sampling line, check that gas values appear on the monitor display.
6. For reasons of hygiene, replace the Airway Adapter within a maximum period of 72 hours.

Environmental Specification

Operating temperature: +10 °C to +40 °C (50 °F to 104 °F)

Storage temperature: -25 °C to +60 °C (-13 °F to 140 °F)

Atmospheric pressure: 660 to 1060 hPa (495 to 795 mmHg)

Ambient Humidity: 10 to 95% RH (non-condensing)

Biocompatibility

The materials that come into contact with the patient have undergone extensive biocompatibility testing. Further information is available on request.

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GE and GE Monogram are trademarks of General Electric Company. All other brands are trademarks of bluepoint medical GmbH & Co. KG.
GE Healthcare Finland Oy
bluepoint medical GmbH & Co. KG