REV.	DESCRIPTION	ER/MR/ECO	SSO	APPROVAL	DATE (MM/DD/YY)
A	Initial release	MR82169	100024626	J. Ikonen	02/04/10
B	Correcting the layout so t all text is visible	099582	100014732	N.Heljo	06/11/10
C	Manufacturer symbol add to order information	100850	100014732	N. Heljo	10/05/10

Languages: Bulgarian, Croatian, Chinese, Czech, Danish, Dutch, English, Estonian, Finnish, French, German, Greek, Hungarian, Italian, Japanese, Norwegian, Latvian, Lithuanian, Portuguese, Romanian, Russian, Serbian, Slovakian, Spanish, Swedish, Turkish (26 languages)

Unfolded paper size (width x height): 23.23" x 21.74" (590mm x 552mm)

Folded paper size (width x height): 3.32" x 7.25" (84mm x 184mm)

Folding (see picture):

1. Start from unfolded sheet.
2. Fold vertically 6 times.
3. Fold horizontally 2 times.

English text shall start from the top left corner after folding.

Material: Paper, white bright, Offset 80g/m2 (Grade 50)


Printing color: black

Double sided printing


Printing quality: IEC 60601-1:2005 subclauses 7.1.2, 7.1.3.

Essential Classification: Minor

*English text here

1.


2.

3.

APPROVALS		DATE	 GE Healthcare	
DRAWN	M. Gellin	09/18/07		
CHECKED	T. Turkki	09/18/07		
			TITLE	
			INSTR TRUSIGNAL SPO2 INTERCONNECT CABLES	
CONFIDENTIAL: THIS PRINT AND THE DESIGN THEREIN REPRESENTED ARE THE PROPERTY OF GE HEALTHCARE AND MAY NOT BE USED IN ANY MANNER DETRIMENTAL TO IT' S INTEREST. THIS PRINT HAS BEEN ISSUED WITH THE UNDERSTANDING THAT IT WILL NOT BE USED FOR MANUFACTURE FOR OTHERS, AND WILL BE RETURNED UPON REQUEST.			DWG. NO.	REV.
			2042137-001	C
			SHEET	1 OF 3

TruSignal® SpO₂ Interconnect Cables

Rx Only (USA)
 US. Federal law restricts this device to sale by or on the order of a physician.

Intended use

The Interconnect Cable is a reusable cable intended for use for all patients for continuous non-invasive arterial oxygen saturation (SpO₂) and pulse rate monitoring when used with a compatible SpO₂ sensor.

Instructions for use

Always refer to the Instructions for use of the sensor and monitor for additional instructions.

1. Monitor compatibility

Use this cable only with patient monitors and sensors compatible with GE Healthcare oximetry. A complete list of compatible monitors available from GE Healthcare upon request.

2. Connecting the cable to a sensor (see Figure 1)

Always connect the cable to the sensor with the aid of alignment feature and push them together. Make sure the connectors are fully engaged.

To disconnect, grasp the connectors (into the cable) and pull straight out.

Figure 1

3. Connecting to the monitor

Align the connector with the matching receptacle on the monitor. Push the connector in fully.

To disconnect, grasp the connector (into the cable) and pull straight out.

TS-M3

To disconnect, press the buttons on the connector.

Maintenance

Cleaning

1. Spray or wipe with pad moisturized with 60% isopropyl or 90% ethyl alcohol or 2% glutaraldehyde. Although disinfection may occur, <2% phenol, <5% formaldehyde, <10% sodium hypochlorite, <0.3% formaldehyde, <0.75% tartaric acid or <5% quarternary ammonium compounds may be used.

2. Let dry completely before use.

Warning

- Do not sterilize or immerse the cable or any part of it in liquid.

Disposal

Dispose of this product and its packaging in accordance with the local environmental and waste disposal regulations.

Warnings

- Use this cable only with patient monitors and sensors compatible with GE Healthcare oximetry. The operator is responsible for checking the compatibility of the sensor and interconnect cable before use. Incompatible components can result in degraded performance.

- Failure to properly connect the cable to the sensor or the monitor may result in intermittent, inaccurate, or no readings.

- This device is not intended for use in a magnetic resonance imaging (MRI) environment.

- Inspect the cable for physical damage. Discard a damaged cable immediately. Never repair a damaged cable; never use a cable repaired by others.

- To avoid damaging the cable when disconnecting, pull from the connector, not the cable. Do not twist or bend the connectors.

- It is possible for any device to malfunction; therefore, always verify unusual data by performing a formal patient assessment.

- The protection against the effect of the discharge of a cardiac defibrillator is partly in the sensor and cable. Using other than GE Healthcare supplies and accessories might result in an impaired patient protection.

- Strangulation may result from baby or child entanglement in pulse oximeter cables.

- A healthcare professional prescribing the use of this device to home use shall provide sufficient guidance to the home user to operate the product safely and correctly.

Technical information

Storage temperature: –40°C to +60°C (-40°F to +140°F)
 Storage humidity: 5% to 95%.

Trademarks

Datex®, Ohmeda®, OxyTr®, TruSignal® and TruTrak® are the property of GE Healthcare Finland Oy. All other product names are the property of their respective owners.

Ordering information

For further information please visit www.gehealthcare.com or contact your local sales company.

Copyright © 2009 General Electric Company. All rights reserved.

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Copyright © 2009 General Electric Company. Всички права запазени.
Copyright © 2009 General Electric Company. Všechna práva vyhrazena.
Copyright © 2009 General Electric Company. Tous droits réservés.
Copyright © 2009 General Electric Company. Alle rettigheder forbeholdes.
Copyright © 2009 General Electric Company. Minden jog fenntartva.
Copyright © 2009 General Electric Company. Alle rechten voorbehouden.