
GE Healthcare

Airway Gas Sampling Set, Disposable

en	Instructions For Use
cs	Návod k použití
da	Instruktioner
de	Gebrauchsanweisung
et	kasutusjuhend
el	Οδηγίες χρήσης
es	Instrucciones de uso
fr	Mode d'emploi
ko	사용 지침
hr	Upute za uporabu
it	Istruzioni per l'uso
lv	Lietošanas instrukcija
lt	Naudojimo instrukcija
hu	Használati útmutató
nl	Gebruiksaanwijzing
ja	取扱説明書
no	Bruksanvisning
pl	Instrukcja użycia
pt	Instruções de utilização
pt-br	Instruções de uso
ro	Instrucțiuni de utilizare
ru	Инструкция по применению
sk	Návod na použitie
sr	Uputstvo za upotrebu
fi	Käyttöohjeet
sv	Bruksanvisning
tr	Kullanım Talimatları
zh	使用说明

Type: Adult/Pediatric, 2m/7ft

REF 8090121521

Type: Infant, 2m/7ft

REF 8090121522

Type: Adult/Pediatric, 3m/10ft

REF 8090121531

Indication for Use

The Airway Gas Sampling Sets are intended for continuous monitoring of expired and inspired respiratory gases: Carbon Dioxide (CO₂) and Oxygen (O₂), and respiratory rate.

The Airway Gas Sampling Sets consist of an airway adapter with gas sampling line for use only with compatible GE Healthcare sidestream respiratory gas modules and patient monitors utilizing these modules.

The airway adapter is provided in two sizes: large for use with Adult/Pediatric patients, and small for use with Infant patients, both via connection to the patient's breathing circuit.

The Airway Gas Sampling Sets are single patient use devices with a cumulative duration of use of up to 72 hours.

The Airway Gas Sampling Sets are intended for hospital care environments where anesthetic agents are not used, such as Intensive Care Units (ICU), Emergency Department (ED), Coronary Care Unit (CCU), Post Anesthesia Care Unit (PACU) and In-hospital Patient Transport.

The Airway Gas Sampling Sets are intended for use by qualified medical personnel only.

The Airway Gas Sampling Sets are for prescription use only.

Warnings

Warning: Do not use to measure anesthetic agent concentration.

Warning: Carefully route the sampling line to reduce the risk of patient entanglement or strangulation.

Warning: Do not cut or remove any part of the sampling line. Cutting the sampling line could lead to erroneous readings.

Warning: Due to risk of cross contamination, do not re-use.

Warning: Loose or damaged connections may compromise ventilation or cause an inaccurate measurement of respiratory gases. Securely connect all components and check connections for leaks according to standard clinical procedures.

Warning: Check CO₂ tubing regularly during use to ensure that no kinks are present. Kinked tubing may cause inaccurate CO₂ sampling.

Warning: Do not use adult or pediatric type sampling line configurations with infants, as this may add dead space to the patient circuit.

Warning: Do not use infant sampling line configurations with adults or pediatric, as this may result in excessive flow resistance.

Warning: Do not place the airway adapter (tee) between the suction catheter and endotracheal tube. This is to ensure that the airway adapter does not interfere with the functioning of the suction catheter.

Cautions

Caution: Do not attempt to clean, disinfect, sterilize or flush any part of the sampling line.

Caution: Dispose of sampling lines in accordance with standard operating procedures or local regulations for the disposal of contaminated medical waste.

Notes

Note: During nebulization or suction, in order to avoid moisture build up and sampling line occlusion, remove the sampling line connector from the monitor.

Note: For accuracy specifications and further instructions, warnings and cautions, please consult the gas monitor instructions for use.

Instruction for Use


Figure 1


Figure 2

1. Connect the Airway Adapter to the proximal end of the breathing circuit between the endotracheal tube and the ventilator Y-section as shown in Figure 1.
2. Connect the gas sample line (male luer connector) to the respiratory gas monitor as shown in Figure 2. Rotate the luer connector clockwise for a secure fit to ensure that there is no leakage of gases at the connection point.
3. After connection of the sampling line, check that gas values appear on the monitor display.
4. For reasons of hygiene, replace Airway Gas Sampling Set within a maximum period of 72 hours.

Environmental Specification

Operating temperature: +10 °C to +40 °C (50 °F to 104 °F)

Storage temperature: -25 °C to +60 °C (-13 °F to 140 °F)

Atmospheric pressure: 660 to 1060 hPa (495 to 795 mmHg)

Ambient Humidity: 10 to 95% RH (non-condensing)

Biocompatibility

The materials that come into contact with the patient have undergone extensive biocompatibility testing. Further information is available on request.

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GE and GE Monogram are trademarks of General Electric Company. All other brands are trademarks of bluepoint medical GmbH & Co. KG.
GE Healthcare Finland Oy
bluepoint medical GmbH & Co. KG