

2089610-001

2087640-001

English

Only use these instructions for flow sensors 2087640-001 and 2089610-001. Keep these instructions with the anesthesia system User's Reference Manual. Refer to the anesthesia system User's Reference Manual for additional information regarding flow sensor use and installation. Replace flow sensor as necessary (dispose of safely if damaged, checkout fails, or visual check continually fails). Under typical use the flow sensor meets specifications for a minimum of 1 year.

CAUTION

Do not use different instructions for flow sensors 2087640-001 and 2089610-001. Damage to the flow sensors or ineffective cleaning, disinfection, and sterilization could occur.

CAUTION

Only parts marked 134°C are autoclavable.

CAUTION

Do not insert any objects into the flow sensor. Damage to the flow sensor can occur. Do not touch the internal flow sensor flap as it may affect flow sensing performance.

CAUTION

Do not allow sensors to dry between soaking, agitation, and rinsing steps to prevent residue as it may affect flow sensing performance.

WARNING

Use personal protective equipment recommended by the chemical and autoclave manufacturers.

Part number	Symbol	Cleaning	Disinfection	Sterilization
2087640-001	134 °C	X		X
2089610-001		X	X	

MRI compatibility information

Flow sensors marked with the MR Conditional symbol may be used with appropriate General Electric MR (Magnetic Resonance) compatible anesthesia devices. Refer to anesthesia device user documentation for conditions of use.

Symbols on the equipment

BCG - Breathing Circuit
Gas

MR conditional

Not autoclavable


134 ° C

Autoclavable

Flow sensor cleaning

Note

This cleaning procedure is applicable to 2087640-001 and 2089610-001 flow sensors. For all soaking, agitation, and rinsing steps, use a minimum of 500 ml of liquid per sensor. Make sure that all surfaces of the flow sensor (except the connector) are submerged in the liquid and there are no visible air bubbles. Do not allow sensors to dry between steps. A Flow Sensor Cleaning Kit M1086468 is available for soaking. Soak only 6 sensors maximum using the cleaning kit as shown.


CAUTION

Do not submerge or get the flow sensor connector wet. Damage to the flow sensor will occur.

1. Mix a mild detergent solution in a clean soaking tray:
 - Use Prolystica 2x Neutral detergent concentration of 1 ml per liter of 20°C to 40°C distilled water. Thoroughly mix the solution.
2. Soak the sensor in the detergent solution for a minimum of 20 minutes.
3. Agitate the flow sensor in the detergent solution sliding it through the solution lengthwise 30 times in a 10 cm path for 30 seconds.

4. Remove the sensor from the solution. Using a syringe, flush each tube with approximately 10 ml of clean detergent solution through the flow sensor connector.


5. Visually inspect internal and external surfaces of the sensor for soil contaminant.
- If loose contaminant remains, repeat steps 3 and 4 in clean detergent solution until removed.
 - If contaminant remains adhered to surfaces, repeat steps 1 through 4 in clean detergent solution until removed.
6. Using a syringe, flush each tube with approximately 10 ml of distilled water through the flow sensor connector.
7. Soak the sensor in fresh distilled water for a minimum of 30 minutes at room temperature.

8. Using a syringe, flush each tube with approximately 10 ml of fresh distilled water through the flow sensor connector.
9. Rinse the sensor by agitating it in room temperature distilled water. Slide the flow sensor lengthwise in fresh distilled water 60 times in a 10 cm path for 60 seconds.
10. Soak the sensor in fresh distilled water for a minimum of 30 minutes at room temperature.
11. Using a syringe, flush each tube with approximately 10 ml of fresh distilled water through the flow sensor connector.
12. Using a dry syringe, push air through the flow sensor connector to dry each tube. Dry each tube until all droplets are cleared.
13. Visually inspect internal and external surfaces of the sensor for soil contaminant or residue. Repeat the cleaning procedure if contaminant or residue remains.
14. Perform the "*Flow sensor disinfection*" or "*Flow sensor sterilization (autoclave)*" procedure.

Flow sensor disinfection

Note

This disinfection procedure is applicable to 2089610-001 flow sensors. For all soaking, agitation, and rinsing steps, use a minimum of 500 ml of liquid per sensor. Make sure that all surfaces of the flow sensor (except the connector) are submerged in the liquid and there are no visible air bubbles. Do not allow sensors to dry between steps. A Flow Sensor Cleaning Kit M1086468 is available for soaking. Soak only 6 sensors maximum using the cleaning kit as shown.


CAUTION

Do not submerge or get the flow sensor connector wet. Damage to the flow sensor will occur.

1. Perform the "*Flow sensor cleaning*" procedure before disinfection.

2. Pour Sporox II disinfectant solution into a clean tray.
3. Agitate the flow sensor in the disinfectant solution sliding it through the solution lengthwise 60 times in a 10 cm path for 60 seconds.
4. Using a syringe, flush each tube with approximately 10 ml of clean disinfectant solution through the flow sensor connector.

5. **Soak the sensor in the disinfectant solution for a minimum of 45 minutes at approximately 20°C.**
6. Using a syringe, flush each tube with approximately 10 ml of distilled water through the flow sensor connector.
7. Rinse the sensor by agitating it in room temperature fresh distilled water. Slide the flow sensor lengthwise in the distilled water 60 times in a 10 cm path for 60 seconds.

8. Soak the sensor in fresh distilled water for a minimum of 30 minutes at room temperature.
9. Using a syringe, flush each tube with approximately 10 ml of fresh distilled water through the flow sensor connector.
10. Rinse the sensor by agitating it in room temperature fresh distilled water. Slide the flow sensor lengthwise in the distilled water 60 times in a 10 cm path for 60 seconds.
11. Soak the sensor in fresh distilled water for a minimum of 30 minutes at room temperature.
12. Rinse the sensor by agitating it in room temperature fresh distilled water. Slide the flow sensor lengthwise in the distilled water 60 times in a 10 cm path for 60 seconds.
13. Using a syringe, flush each tube with approximately 10 ml of fresh distilled water through the flow sensor connector.
14. Using a dry syringe, push air through the flow sensor connector to dry each tube. Dry each tube until all droplets are cleared.
15. Completely dry the sensor in an oven at 105°C to 110°C for a minimum of 30 minutes.

16. Visually inspect internal and external surfaces of the sensor for soil contaminant or residue. Repeat "*Flow sensor cleaning*" and "*Flow sensor disinfection*" if contaminant or residue remains.

Flow sensor sterilization (autoclave)

This sterilization procedure is applicable to 2087640-001 flow sensors.

CAUTION

Only parts marked 134°C are autoclavable.

1. Perform the "*Flow sensor cleaning*" procedure before sterilization.
2. Wrap the flow sensor in light-gauge sterilization wrap for autoclaving (one-ply polypropylene wrap Kimberly-Clark Kimguard KC600 or equivalent).
3. Sterilize in an autoclave at 132°C for 4 minutes, or 134°C for 3 minutes, using the pre-vacuum dynamic air removal cycle with a drying time of 30 minutes for either cycle. Do not exceed 134°C.
4. Allow the flow sensor to cool to room temperature before use.

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© Datex-Ohmeda, Inc. – All rights reserved.
GE and GE Monogram are trademarks of
Datex-Ohmeda, Inc., a General Electric Company,