
尺寸: W137mmHxH95mm

unit or other RF source. If probes must be used together with electro-surgical apparatus, the instrument cables must be connected to the patient. The probes must be connected to the patient and the patient must be grounded. The probes are connected should be checked for adequate isolation from electrical grounds at radio frequencies.

Caution:

Overheating of the probes could result in damage to internal wires and loss of electrical isolation or improper measurement readings.

Operation procedures:

To ensure the compatibility of the probe, see the attached the compatible label to ensure the device is compatible with the temperature probe.

Always use the probe plug matches the connector on the device. Then fully plug the probe into the device. The probe plug is located on the side of the probe. The probe plug is used to connect the probe to the device. The probe plug is located on the side of the probe. The probe plug is used to connect the probe to the device. The probe plug is located on the side of the probe. The probe plug is used to connect the probe to the device.

Locate the probe plug on the site where is intended to be measured. Fix the probe tip with medical tape or other applicable.

Start to measure and get the measurement readings after the response time (see above).

Probes should be cleaned before disinfection or sterilization to improve the efficiency.

Probes may be disinfected by washing with 70% isopropanol, ethanol. After washing, probes should be rinsed thoroughly with water. Brief immersion of the probe in detergent solutions is not harmful.

Never boil the T-Series reusable temperature probes.

Storage and handling

Probes should be stored in a clean, dry container and should be loosely coated and stored in room temperature. Don't wrap probes around equipment cases to avoid damaging internal wires.

Environment operating range:

Altitude: 0 to 10,000 ft

Temperature: -40°C to 150°C

Humidity: 0 to 100% RH

Warranty:

Altitude: Regular lifetime of Unimed temperature probes is two years. All T-Series reusable temperature probes carry one year warranty for workmanship, components and accurate tolerances. Damage through misuse or modification is not covered. Probe life with normal use should exceed one year.

CAUTION:

Federal U.S.L.A restricts this device to sale by or on the order of a physician.

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Description générale :

Catient	Lot Number	Manufacturer	Serialisation	Lotus Tree	Self-Instructions	Non-Destructive and Destructive Examinations
	LOT					

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