Preattached Leadwire ECG Electrodes

Instructions for Use

Metallic Preattached Leadwire Electrode, Medium

REF 989803100521 (13952D)

Metallic Preattached Leadwire Electrode, Square

REF 989803100561 (13953D)

Radiolucent Preattached Leadwire Electrode, Square

REF 989803100571 (13953E)

PHILIPS

EN - Instructions for Use.....	1	LV - Lietošanas instrukcija	10
BG - Указания за употреба	1	MK - Упатства за користење	10
CS - Návod k obsluze	2	NL - Gebruiksaanwijzing	11
DA - Brugervejledning	2	NO - Brukerhåndbok	11
DE - Gebrauchsanweisung	3	PL - Instrukcja obsługi	12
EL - Οδηγίες χρήσης	3	PT-BR - Instruções de uso	12
ES - Instrucciones de uso	4	PT-PT - Instruções de utilização	13
ET - Kasutusjuhend	4	RO - Instrucțiuni de utilizare	13
FI - Käyttöohjeet	5	RU - Инструкция по эксплуатации	14
FR - Manuel d'utilisation	5	SK - Návod na obsluhu	15
HR - Upute za uporabu	6	SL - Navodila za uporabo	15
HU - Használati útmutató	6	SR - Uputstvo za upotrebu	16
ID - Petunjuk Penggunaan	7	SV - Bruksanvisning	16
IT - Istruzioni d'uso	7	TR - Kullanım Talimatları	17
JA - ユーザーズガイド	8	UK - Інструкції з використання	17
KK - Пайдалану нұсқаулығы	8	VI - Hướng Dẫn Sử Dụng	18
KO - 사용 설명서	9	ZH-CN - 使用说明	18
LT - Naudojimo nurodymai	9	ZH-TW - 使用說明	18

		
Prep	Peel	Apply
Подготовка	Отлепняне	Прилагане
Připrava	Odlepení	Aplikace
Klargør	Træk af	Sæt på
Vorbereiten	Abziehen	Anbringen
Προετοιμασία	Αποκόλληση	Εφαρμογή
Preparar	Despegar	Aplicar
Valmistuge	Tõmmake lahti	Paigaldage
Valmistele	Irrota	Käytä
Préparer	Décoller	Appliquer
Připemite	Odlíjepite	Primijenite
Készítse elő	Húzza le	Helyezze fel
Persiapan	Kelupas	Pasang
Preparare	Staccare	Applicare
準備	電極をはがす	装着
Дайындау	Электродты ажырату	Қолдану
준비	벗기기	부착
Paruošti	Lupti	Taikyti
Sagatavot	Noņemot pārklājumu	Lietot
Подготви	Одлепи	Нанеси
Vorbereiden	Verwijderen	Aanbrengen
Klargør	Ta av	Bruk
Przygotuj skórę	Odklej	Zamocuj
Preparar	Remover	Aplicar
Preparar	Descolor	Aplicar
Pregătire	Dezlipire	Aplicare
Подготовка	Отделение	Наложение
Připrava	Odlepenie	Aplikácia
Připravite	Odlepite	Namestite
Připemite	Odlepite	Primenite
Förbered	Lossa	Applicera
Hazırlayın	Soyun	Uygulayın
Підготовка	Зняття	Застосування
Chuẩn bị	Bóc	Đán
准备	剥离	佩戴
準備	撕下	貼附

	EN: Manufacturer; BG: Производител; CS: Výrobce; DA: Producent; DE: Hersteller; EL: Κατασκευαστής; ES: Fabricante; ET: Tootja; FI: Valmistaja; FR: Fabricant; HR: Proizvođač; HU: Gyártó; ID: Produsen; IT: Fabriccante; JA: 製造業者; KK: Өндүрүшү; KO: 제조자; LT: Gamintojas; LV: Ražotājs; MK: Производител; NL: Fabrikant; NO: Produsent; PT: Productent; PT-BR: Fabricante; PT-PT: Fabricante; RO: Producător; RU: Производитель; SK: Výrobca; SL: Proizvajalec; SR: Proizvođač; SV: Tillverkare; TR: Üretici; UK: Виробник; VI: Nhà sản xuất; ZH-CN: 制造商; ZH-TW: 製造商
	EN: Electronic Instructions for Use; BG: Електронни указания за употреба; CS: Elektronický návod k obsluze; DA: Elektronisk brugsvejledning; DE: Elektronische Gebrauchsanweisung; EL: Οδηγίες χρήσης σε ηλεκτρονική μορφή; ES: Instrucciones de uso electrónicas; ET: Elektroniline kasutusjuhend; FI: Sähköiset käyttöohjeet; FR: Manuel d'utilisation électronique; HR: Elektroničke upute za uporabu; HU: Elektronikus használati útmutató; ID: Petunjuk Penggunaan Elektronik; IT: istruzioni d'uso in formato elettronico; JA: 電子版ユーザーズガイド; KK: Электрондык пайдалану нускаулыгы; KO: 전자 사용 설명서; LT: Elektroninės naudojimo instrukcijos; LV: Elektroniska lietošanas instrukcija; MK: Електронски упутства за користење; NL: Elektronische gebruiksaanwijzing; NO: Elektronisk brukerhåndbok; PL: Instrukcja obsługi w wersji elektronicznej; PT-BR: nstruções eletrônicas de uso; PT-PT: Instruções de utilização em formato eletrónico; RO: Instrucțiuni de utilizare în format electronic; RU: Электронная версия инструкции по эксплуатации; SK: Elektronický návod na použitie; SL: Elektronska navodila za uporabo; SR: Elektronsko uputstvo za upotrebu; SV: Elektronisk bruksanvisning; TR: Elektronik Kullanım Talimatları; UK: Електронні інструкції з використання; VI: Hướng dẫn sử dụng bản điện tử; ZH-CN: 电子版使用说明; ZH-TW: 電子版使用說明
	EN: Batch Code; BG: Партиден код; CS: Kód série; DA: Kode for parti; DE: Chargenbezeichnung; EL: Κωδικός παρτίδας; ES: Código de lote; ET: Partiiikood; FI: Eräkoodi; FR: Code de lot; HR: Šifra serije; HU: Tételszám; ID: Kode Batch; IT: Codice lotto; JA: 製造番号又は製造記号; KK: Партия коды; KO: 배치 코드; LT: Partijos kodas; LV: Sērijas numurs; MK: Код на серија; NL: Batchcode; NO: Lotnummer; PL: Kod partii; PT-BR: Código de lote; PT-PT: Código de lote; RO: Codul lotului; RU: Код партии; SK: Kód šarže; SL: Številka serije; SR: Broj serije; SV: Batchkod; TR: Parti Kodu; UK: Код партії; VI: Mã lô; ZH-CN: 批号; ZH-TW: 批次碼
	EN: Temperature Limits; BG: Температурни граници; CS: Limity teploty; DA: Temperaturgrænser; DE: Zulässiger Temperaturbereich; EL: Ορια θερμοκρασίας; ES: Límites de temperatura; ET: Temperatuuripiirangud; FI: Lämpötilarajat; FR: Limites de température; HR: Temperature granice; HU: Hőmérsékleti határértékek; ID: Batasan Suhu; IT: Limiti di temperatura; JA: 温度制限; KK: Температура шектери; KO: 온도 제한; LT: Temperatūros ribos; LV: Temperatūras robežvērtības; MK: Температурни граници; NL: Temperatuurlimieten; NO: Temperaturgrenser; PL: Zakresy temperatury; PT-BR: Limites de temperatura; PT-PT: Limites de temperatura; RO: Limite de temperatură; RU: Температурные ограничения; SK: Hraníčné hodnoty teploty; SL: Temperature omejitve; SR: Ograničenja temperature; SV: Temperaturgränser; TR: Sıcaklık Limitleri; UK: Обмеження температури; VI: Giới hạn nhiệt độ; ZH-CN: 温度限制; ZH-TW: 溫度限制
	EN: Humidity Limitation; BG: Ограничение на влажността; CS: Limit vlhkosti; DA: Fugtighedsgrænse; DE: Zulässiger Luftfeuchtigkeitsbereich; EL: Ορια υγρασίας; ES: Límites de humedad; ID: Batasan Kelembapan; IT: Limiti di umidità; JA: 湿度制限; KK: Ылгалдылык шери; KO: 습도 제한; LT: Drėgnės ribojimas; LV: Mitruma ierobežojums; MK: Ограничување на влажност; NL: Vochtigheidsgrens; NO: Fuktighetsgrenser; PL: Zakres wilgotności; PT-BR: Limite de umidade; PT-PT: Limitação da humidade; RO: Limitare umiditate; RU: Ограничения по влажности; SK: Hraníčné hodnoty vlhkosti; SL: Omejitve glede vlažnosti; SR: Ograničenje vlažnosti; SV: Luftfuktighetsbegränsning; TR: Nem Sınırlaması; UK: Обмеження вологості; VI: Giới hạn độ ẩm; ZH-CN: 湿度限制; ZH-TW: 濕度限制
	EN: Atmospheric Pressure Limitation; BG: Ограничение на атмосферното налягане; CS: Omezení atmosférického tlaku; DA: Grænse for atmosfærisk tryk; DE: Zulässiger Luftdruckbereich; EL: Ορια ατμοσφαιρικής πίεσης; ES: Límites de presión atmosférica; ET: Õhurõhu piirväärtused; FI: Ilmanpaineajat; FR: Limite de pression atmosphérique; HR: Ograničenje atmosferskog tlaka; HU: Légnyomás-határérték; ID: Batasan Tekanan Atmosfer; IT: Limitazione di pressione atmosferica; JA: 大気圧制限; KK: Атмосфералык кысым шери; KO: 대기압 제한; LT: Atmosferos slėgio apribojimas; LV: Atmosfēras spiediena ierobežojums; MK: Ограничување за атмосферски притисок; NL: Luchtdruk grens; NO: Grense for atmosfærisk tryk; PL: Zakres ciśnienia atmosferycznego; PT-BR: Limite da pressão atmosférica; PT-PT: Limitação da pressão atmosférica; RO: Limitația de presiune atmosferică; RU: Ограничения по атмосферному давлению; SK: Hraníčné hodnoty atmosférického tlaku; SL: Omejitve glede atmosferskega tlaka; SR: Ograničenje atmosferskog pritiska; SV: Atmosfärstrycksbegränsning; TR: Atmosferik Basınc Sınırlaması; UK: Обмеження атмосферного тиску; VI: Giới hạn áp suất khí quyển; ZH-CN: 大气压限制; ZH-TW: 大氣壓力限制
	EN: Keep Dry; BG: Да се съхранява на сухо място; CS: Uchovávejte v suchu; DA: Opbevares tørt; DE: Vor Feuchtigkeit schützen; EL: Διατηρείτε στεγνό; ES: Mantener en lugar seco; ET: Hoida kuivas; FI: Suojattava kosteudelta; FR: Conserver au sec; HR: Držite na suhom; HU: Szárzandó tartandó; ID: Jaga Tetap Kering; IT: Conservare in luogo asciutto; JA: 水濡れ不可; KK: Күрпак жерде сактанды; KO: 건조한 상태 유지; LT: Laikyti sausoje vietoje; LV: Glabāt sausā vietā; MK: Да се чува на суво место; NL: Droog houden; NO: Holdes tørt; PL: Przechowywać w suchym miejscu; PT-BR: Manter seco; PT-PT: Manter seco; RO: A se feri de umezeală; RU: Хранить в сухом месте; SK: Uchovávať v suchu; SL: Hraniti na suhom; SR: Držite na suvom; SV: Produkten ska hållas torr; TR: Kuru Olarak Muhafaza Edin; UK: Зберігати в сухому місці; VI: Bảo quản nơi khô ráo; ZH-CN: 保持干燥; ZH-TW: 保持乾燥

EN: Keep Away from Sunlight; BG: Да се пази от слънчева светлина; CS: Chraňte před slunečními paprsky; DA: Holdes væk fra sollys; DE: Vor Sonneneinstrahlung schützen; EL: Διατηρείται μακριά από ηλιακό φως; ES: Proteger de la luz solar; ET: Hoida eemal päikesevalgusest; FI: Säilytä auringolla suojattuna; FR: Conserver à l'abri du soleil; HR: Držite podalje od sunca; HU: Tartsa napfénytől távol; ID: Hindarkan dari Sinar Matahari; IT: Non esporre alla luce del sole; JA: 直射日光を避けること; KK: Күн сәулесі түспейтін жерде сақтаңыз; KO: 직사광선 노출 금지; LT: Saugoti nuo saulės spindulių; LV: Neglabājet tiešos saules staros; MK: Не изложувајте го производот на сончева светлина; NL: Uit de buurt van zonlicht houden; NO: Hold unna sollys; PL: Chronić przed bezpośrednim nasłonecznieniem; PT-BR: Manterha longe da luz do sol; PT-PT: Manter ao abrigo da luz solar; RO: A se feri de lumina soarelui; RU: Предохранять от солнечных лучей; SK: Nevystavujte slnečnému žiareniu; SL: Ne izpostavljajte sončni svetlobi; SR: Držati dalje od sunčeve svetlosti; SV: Skyddas från solljus; TR: Güneş ışığından Uzak Tutun; UK: Зберігайте в захищеному від сонця місці; VI: Tránh xa ánh sáng mặt trời; ZH-CN: 避开阳光; ZH-TW: 避免日光照射

EN: Prescription Use Only; BG: Само по лекарско предписание; CS: Používejte pouze na lékařský předpis; DA: Receptpligtigt; DE: Anwendung nur auf ärztliche Verordnung; EL: Χρήση μόνο κατόπιν ιατρικής συνταγής; ES: Solo bajo prescripción médica; ET: Kasutamise ainult retsepti alusel; FI: Käyttö vain lääkärin määräyksellä; FR: Sur prescription uniquement; HR: Samo na recept; HU: Csak orvosi rendelvényre; ID: Hanya dengan Resep Dokter; IT: Solo su prescrizione medica; JA: 医家向け医療機器; KK: Тек дөріөрдүн тагайындауы бойынша пайдаланылады; KO: 처방 필요; LT: Naudoiti tik skyrus gydytojų; LV: Lietošani tikai pēc norīkojuma; MK: Само за препиша на употреба; NL: Uitsluitend op recept; NO: Reseptbelagt; PL: Wyłączenie na zlecenie lekarza; PT-BR: Uso somente sob prescrição; PT-PT: Apenas com prescrição médica; RO: Doar la recomandarea medicului; RU: Только по назначению врача; SK: Iba na lekársky predpis; SL: Izdaja se samo na recept; SR: Samo na lekarski recept; SV: Endast på ordination av läkare; TR: Sadece Reçeteye Kullandır; UK: Використовувати лише за приписом лікаря; VI: Chỉ sử dụng theo toa; ZH-CN: 仅凭处方; ZH-TW: 僅限處方使用

EN: Country of Manufacture; BG: Държава на производство; CS: Země výroby; DA: Fremstillingsland; DE: Herstellungsland; EL: Χώρα κατασκευής; ES: País de fabricación; ET: Tootmisriik; FI: Valmistusmaa; FR: Pays de fabrication; HR: Država proizvodnje; HU: Gyártás ország; ID: Negara Produsen; IT: Paese di fabbricazione; JA: 製造国; KK: Өндүрүш өлкө; KO: 제조국; LT: Pagaminimo šalis; LV: Ražotājtvalsts; MK: Земја на производство; NL: Land van productie; NO: Produktionsland; PL: Kraj produkcji; PT-BR: País de fabricação; PT-PT: País de fabrico; RO: Țara de producție; RU: Страна-производитель; SK: Krajina výroby; SL: Država proizvodnje; SR: Zemlja proizvodnje; SV: Tillverkningsland; TR: Üretim ülkesi; UK: Країна-виробник; VI: Quốc gia sản xuất; ZH-CN: 生产国家/地区; ZH-TW: 製造國家/地區

EN: Medical Device; BG: Медицинско изделие; CS: Zdravotnický přístroj; DA: Medicinsk udstyr; DE: Medizinprodukt; EL: Ιατροτεχνολογικό προϊόν; ES: Dispositivo médico; ET: Meditsiiniseade; FI: Lääkinnällinen laite; FR: Dispositif médical; HR: Medicinski proizvod; HU: Orvostechnikai eszköz; ID: Perangkat Medis; IT: Dispositivo medico; JA: 医療機器; KK: Медициналык құрылғы; KO: 의료기기; LT: Medicinos prietaisas; LV: Medicīnas ierīce; MK: Медицински уред; NL: Medisch apparaat; NO: Medicinsk utstyr; PL: Wyrób medyczny; PT-BR: Dispositivo médico; PT-PT: Dispositivo médico; RO: Dispozitiv medical; RU: Медицинское устройство; SK: Zdravotnícka pomôcka; SL: Medicinski pripomoček; SR: Medicinski uređaj; SV: Medicinteknisk produkt; TR: Tıbbi Cihaz; UK: Медичний пристрій; VI: Thiết bị y tế; ZH-CN: 医疗设备; ZH-TW: 醫療設備

EN: Authorized Representative in the European Community; BG: Упълномощен представител в Европейската общност; CS: Autorizovaný zástupce v Evropském společenství; DA: Autoriseret repræsentant i EU; DE: Autorisierte EU-Vertretung; EL: Εξουσιοδοτημένος αντιπρόσωπος για την Ευρωπαϊκή Ένωση; ES: Representante autorizado en la Unión Europea; ET: Volitatud esindaja Euroopa Ühenduses; FI: Valtuutettu edustaja Euroopan yhteisössä; FR: Représentant autorisé pour la Communauté européenne; HR: Ovlašteni zastupnik u Europskoj uniji; HU: Hivatalos képviselő az Európai Közösségekben; ID: Perwakilan Uni Eropa Resmi; IT: Rappresentante autorizzato nella Comunità europea; JA: Authorized EU Representative; KK: Еуропалык қауымдастықтағы уәкілетті өкіл; KO: 유럽 공인 대리장; LT: Įgaliojatis atstovas Europos Bendrijoje; LV: Pilnvarotais pārstāvis Eiropas Kopienā; MK: Овластен претставник во Европската заедница; NL: Officiële vertegenwoordiger in de Europese Unie; NO: Autorisert EU-representant; PL: Autoryzowany przedstawiciel na terenie Wspólnoty Europejskiej; PT-BR: Representante autorizado na União Europeia; PT-PT: Representante autorizado na Comunidade Europeia; RO: Reprezentant autorizat în Comunitatea Europeană; RU: Уполномоченный представитель в странах Европейского сообщества; SK: Autorizovaný zástupca pre Európske spoločenstvo; SL: Pooblaščen zastopnik za Evropsko skupnost; SR: Ovlašteni predstavnik Evropske zajednice; SV: Auktoriserad representant i Europeiska unionen; TR: Yetkili Avrupa Topluluğu Temsilcisi; UK: Уповноважений представник на території ЄС; VI: Đại diện được ủy quyền ở Cộng đồng châu Âu; ZH-CN: 欧盟授权代表; ZH-TW: 授權的歐洲共同體代表

EN: Unique Device Identifier; BG: Уникален идентификатор на устройството; CS: Jedinčný identifikátor zařízení; DA: Unik enhedsidentifikator; DE: Eindeutige Geräteerkennung; EL: Αποκλειστικό αναγνωριστικό τεχνολογικού προϊόντος; ES: Identificador exclusivo del dispositivo; ET: Seadme kordumatu identifitseerimistunnus; FI: Laitteen yksilöivä tunnus; FR: Identifiant unique de l'appareil; HR: Jedinstveni identifikator uređaja; HU: Egyedi eszközazonosító; ID: Pengenal Perangkat Unik; IT: Identificatore univoco del dispositivo; JA: 機器固有識別子; KK: Құрылғының бірегей есшімдік факторы; KO: 고유 장치 식별 기호; LT: Unikalusis prietaiso identifikatorius; LV: īpašs unikālais identifikators; MK: Уникатна идентификација на уредот; NL: Unieke apparaat-ID; NO: Unik enhetsidentifikator; PL: Unikalny identyfikator wyrobu; PT-BR: Identificador único do dispositivo; PT-PT: Identificador único do dispositivo; RO: Identifier unic al dispozitivului; RU: Уникальный идентификатор устройства; SK: Jedinčný identifikátor zariadenia; SL: Edinstven identifikator pripomočka; SR: Jedinstvena identifikacija uređaja; SV: Unik produktidentifier (UDI, Unique Device Identifier); TR: Benzersiz Cihaz Tanımlayıcı; UK: Унікальний ідентифікатор пристрою; VI: Mã nhận dạng thiết bị duy nhất; ZH-CN: 医疗器械唯一标识; ZH-TW: 唯一裝置識別碼

EN: Catalog Number; BG: Каталоген номер; CS: Katalogové číslo; DA: Katalognummer; DE: Katalognummer; EL: Αριθμός καταλόγου; ES: Número de catálogo; ET: Katalooginumber; FI: Luettelonumero; FR: Numéro de référence; HR: Kataloški broj; HU: Katalógusszám; ID: Nomor Katalog; IT: Numero di catalogo; JA: カタログ番号; KK: Каталог бойынша; KO: 카탈로그 번호; LT: Katalogo numeris; LV: Kataloga numurs; MK: Каталогски број; NL: Catalogusnummer; NO: Katalognummer; PL: Numer katalogowy; PT-BR: Número de catálogo; PT-PT: Número de catálogo; RO: Număr de catalog; RU: Номер по каталогу; SK: Katalógové číslo; SL: Kataloška številka; SR: Kataloški broj; SV: Katalognummer; TR: Katalog Numarası; UK: Номер за каталогом; VI: Số catalog; ZH-CN: 目录号; ZH-TW: 目錄編號

	EN: Use By Date; BG: Срок на годност; CS: Spotřebujte před datem použitelnosti; DA: Sidste anvendelsesdato; DE: Verfallsdatum; EL: Ημερομηνία λήξης; ES: Fecha de caducidad; ET: Aegumiskupaevä; FI: Viimeinen käyttöpäivä; FR: Date limite d'utilisation; HR: Upotrijebiti do; HU: Lejárati dátum; ID: Tanggal Kedaluwarsa; IT: Utilizzare entro; JA: 使用期限; KK: Жарамдылық мерзімінің біту күні; KO: 유통 기한; LT: Naudojti iki; LV: Izlietot līdz; MK: Пок на употреба; NL: Uiterste gebruiksdatum; NO: Utløpsdato; PL: Termin przydatności; PT-BR: Uso por data; PT-PT: Data de validade; RO: A se utiliza până la data de; RU: Дата истечения срока годности; SK: Použite do; SL: Rok uporabe; SR: Upotrijebiti do datuma; SV: Utgångsdatum; TR: Son Kullanım Tarihi; UK: Термін придатності; VI: Hạn sử dụng; ZH-CN: 有效使用期限; ZH-TW: 使用期限
	EN: Do Not Re-use; BG: Да не се използва повторно; CS: Nepoužívejte opakovaně; DA: Må ikke genbruges; DE: Nicht wiederverwenden; EL: Μην επαναχρησιμοποιείτε; ES: No reutilizar; ET: Mitte korduskasutada; FI: Kertäkäyttöinen; FR: Ne pas réutiliser; HR: Nemojte ponovo upotrebljavati; HU: Ne használja fel többször; ID: Jangan Dipakai Ulang; IT: Non riutilizzare; JA: 再使用不可; KK: Кайта пайдалануу болмайды; KO: 재사용 금지; LT: Nenaudoti pakartotinai; LV: Nelietot atkārtoti; MK: Да не се користи повторно; NL: Niet hergebruiken; NO: Må ikke gjenbrukes; PL: Nie używać ponownie; PT-BR: Não reutilize; PT-PT: Não reutilizar; RO: A nu se reutiliza; RU: Не используйте повторно; SK: Nepoužívaťe opakovaně; SL: Samo za enkratno uporabo; SR: Nemojte ponovo da koristite; SV: Återanvänd inte; TR: Yeniden Kullanmayın; UK: Не можна повторно використовувати; VI: Không được tái sử dụng; ZH-CN: 不要重复使用; ZH-TW: 請勿重複使用
	EN: Duration for Use (Hours); BG: Продължителност на употребата (часове); CS: Doba používání (hodiny); DA: Varighed af brug (timer); DE: Anwendungsdauer (Stunden); EL: Διάρκεια χρήσης (ώρες); ES: Duración del uso (horas); ET: Kasutamisaeg (tunnid); FI: Käyttöaika (tunteina); FR: Durée d'utilisation (heures); HR: Trajanje uporabe (u satima); HU: Használat időtartama (óra); ID: Durasi Penggunaan (Jam); IT: Durata di utilizzo (ore); JA: 使用期間 (時間); KK: Пайдалануу узактыгы (саат); KO: 사용 기간(시간); LT: Naudojimo trukmė (valandomis); LV: Lietošanas ilgums (stundas); MK: Период на користење (часови); NL: Gebruiksduur (uren); NO: Varighet for bruk (timer); PL: Czas używania (godziny); PT-BR: Duração do uso (horas); PT-PT: Duração da utilização (horas); RO: Durata de utilizare (ore); RU: Длительность использования (часы); SK: Doba použitia (v hodinách); SL: Trajanje uporabe (ure); SR: Trajanje upotrebe (sati); SV: Varaktighet för användning (timmar); TR: Kullanım Süresi (Saat); UK: Тривалість використання (години); VI: Thời lượng sử dụng (Giờ); ZH-CN: 使用持续时间 (小时); ZH-TW: 使用持續時間 (小時)
	EN: Packaging Unit; BG: Опаковъчна единица; CS: Balící jednotka; DA: Emballageenhed; DE: Verpackungseinheit; EL: Μονάδα συσκευασίας; ES: Unidad de embalaje; ET: Arv karbis; FI: Pakkausyksikkö; FR: Unité de conditionnement; HR: Jedinica pakiranja; HU: Csomagolási egység; ID: Unit Kemasan; IT: Unità di imballaggio; JA: パッケージに含まれる個数; KK: Қаптамадағы бірлік; KO: 포장 단위; LT: Pakuotė; LV: epakojuuma vienība; MK: Пакување; NL: Verpakkingseenheid; NO: Emballasjeenhet; PL: Jednostka opakowania; PT-BR: Unidade de embalagem; PT-PT: Unidade de embalagem; RO: Unitate de ambalare; RU: Упаковочная единица; SK: Jednotka balenia; SL: Količina v embalaži; SR: Broj artikala u pakovanju; SV: Förpackningsenhet; TR: Ambalaj Birimi; UK: Пакувальна одиниця; VI: Đơn vị đóng gói; ZH-CN: 包装单位; ZH-TW: 包裝單位
	EN: Caution; BG: Внимание; CS: Varování; DA: Varsigt; DE: Vorsicht; ES: Precaución; ET: Arv karbis; FI: Varoituis; FR: Attention; HR: Oprez; HU: Vigyázat; ID: Perhatian; IT: Attenzione; JA: 注意; KK: Ескерту; KO: 주의; LT: Paspėjimas; LV: Piesardzība; MK: Внимание; NL: Let op; NO: OBS!; PL: Przestroga; PT-BR: Cuidado; PT-PT: Cuidado; RO: Atenție; RU: Внимание!; SK: Varovanie; SL: Pozor; SR: Oprez; SV: Viktigt; TR: Dikkat; UK: Увага!; VI: Thận trọng; ZH-CN: 小心; ZH-TW: 小心事項
	EN: Consult Instruction for use; BG: Консултирайте се с указанията за употреба; CS: Prostudujte si návod k obsluze; DA: Se brugsanvisningen; DE: Gebrauchsanweisung beachten; EL: Συμβουλευτείτε τις οδηγίες χρήσης; ES: Consulte las instrucciones de uso; ET: Lugege kasutusjuhendit; FI: Lue käyttöohjeet; FR: Consulter le manuel d'utilisation; HR: Pročitajte upute za uporabu; HU: Lásd a használati útmutatót; ID: Baca Petunjuk Penggunaan; IT: Consultare le Istruzioni d'uso; JA: ユーザーズガイド参照; KK: Пайдалануу нускаулыгынан қараңыз; KO: 사용 설명서 참조; LT: Žr. naudojimo instrukcijas; LV: Izlasiet lietošanas instrukciju; MK: Разгледajte ги упатствата за користење; NL: Raadpleeg de gebruiksaanwijzing; NO: Se brukenhåndboken; PL: Zapoznaj się z instrukcją obsługi; PT-BR: Consulte as instruções de uso; PT-PT: Consultar as instruções de utilização; RO: Consultați instrucțiunile de utilizare; RU: См. инструкции по эксплуатации; SK: Pozrite si návod na použitie; SL: Glejte navodila za uporabo; SR: Pogledajte uputstvo za upotrebu; SV: Läs bruksanvisningen; TR: Kullanım Talimatlarını Başvurun; UK: Див. Інструкції з використання; VI: Tham khảo hướng dẫn sử dụng; ZH-CN: 参见使用说明; ZH-TW: 請見使用說明

EN -Instructions for Use

Intended Use/Indications for Use - Electrocardiogram (ECG) electrodes are used to sense the neonatal heart rate. The electrocardiograph electrode is the electrical conductor which is applied to the surface of the body to transmit the electrical signal at the body surface to a processor that produces electrocardiogram. **Intended Users** - Intended for use by trained medical professionals and lay persons (in the case of continuous ambulatory monitoring). **Contraindications** - There are no known contraindications. **Patient Population** - Neonatal/Child/Pediatric. **Compatibility** - The operator is responsible for checking compatibility of Philips monitoring equipment and these electrodes. Incompatible components may result in degraded performance. **Patient Application** - To apply the solid gel electrodes to a patient see illustrations on page i. **Environment**: For use in healthcare environments.

Warnings - Employ procedures recommended by electrosurgical or monitoring equipment manufacturers to minimize risk of burns at the electrode site during electrosurgical procedures. • Do not use solvents to clean skin prior to ECG. Solvents trapped under electrodes may lead to abnormal skin reactions. • Do not apply over open wounds, lesions infected, irritated or inflamed areas to prevent further skin damage. • Do not allow monitoring electrodes or leadwires to touch defibrillation pads. Contact can cause electrical arcing, patient burns during defibrillation and may divert therapy current away from the heart. • Do not reuse due to the risk of cross infection, degradation of adhesive or electrical performance. • Do not use on patients with a history of skin sensitization or irritation to ECG electrodes. Report any observed adverse reactions to Philips. • Use caution when removing electrodes to prevent skin damage. • Dispose of electrodes properly. • To prevent possible patient aspiration, properly dispose electrode backing plastic liner. • Electrodes are not designed for use in an MRI environment. Potential for patient burns exists. • To avoid patient injury, carefully position leadwires to avoid entanglement, choking, strangulation or inhibiting circulation in patient extremities. • Ensure leadwire connectors never come into contact with other conductive parts or with earth due to the risk of electrical shock or burn to patient or user. • These leads are unshielded against electromagnetic interference. They are not to be used with defibrillators for monitoring during transcutaneous pacing or cardioversion. • Do not use electrodes for more than 72 hours. • Attempts to clean or sterilize the electrodes may result in bio-incompatibility, infection or product failure risks to the patient.

Cautions - Federal (US) law restricts this device to sale by or on the order of a medical practitioner. • Replace electrodes if they no longer firmly stick to patient skin. • Do not use electrodes beyond the expiration date marked on the package. • Ensure application sites are properly prepared - dry, clean and free of excessive hair. Inaccurate readings could result. • Use care when applying electrodes to avoid gel displacement. Loss of adhesion could occur. • Ensure the leadwire connector is properly inserted into the trunk cable. • These leadwire electrodes are not designed for telemetry use. • The 13952D and 13953D electrodes are not designed for use with X-ray equipment. X-ray images are compromised if the electrode is placed in the area of interest, leading to delay of diagnosis. • Do not open pouch until immediately prior to use. • Storing outside the recommended storage temperature range reduces shelf life.

Reordering Information - See front cover for reordering information. **Disposal** - Follow approved medical waste disposal methods as specified by your patient care facility or local regulations. **Incident Reporting** - Any serious incident that has occurred in relation to this device should be reported to Philips and the Competent Authority of the European Economic Area (EEA) countries, including Switzerland and Turkey, in which the user and/or patient is established. **Storage Environmental Specifications** - 12°C to 32°C (53.6°F to 89.6°F) at 15% to 90% relative humidity and 570 hPa to 1000 hPa atmospheric pressure.

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Copyright © 2021
Nederland B.V.