

ECG Monitoring/Diagnostic Snap Electrodes

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EN - Instructions for Use

Intended Use/Indications for Use

M2202A: Radiolucent and MRI-compatible ECG electrodes are designed for use in general electrocardiographic procedures in healthcare environments where ECG monitoring is deemed necessary and is ordered by a physician. Such procedures include in particular patient ECG surveillance and ECG diagnosis recording. Radiolucent and MRI-compatible ECG electrodes are single-use, non-sterile, disposable and are to be used on intact (uninjured) skin.

40420A, 40493D, 40493E, M4612A, M4613A and 989803192551: ECG electrodes are designed for use in general electrocardiographic procedures in healthcare environments where ECG monitoring is deemed necessary and is ordered by a physician. Such procedures include in particular patient ECG surveillance and ECG diagnosis recording. ECG electrodes are non-sterile and are to be used on intact (uninjured) skin.

Intended Users

Intended for use by licensed healthcare professionals.

Contraindications

There are no known contraindications.

Patient Population

Adult: 40493D, 40493E, M2202A, M4612A, M4613A

Adult/Pediatric: 989803192551, 40420A

WARNINGS

- Do not allow monitoring electrodes or lead wires to touch defibrillation pads. Such contact can cause electrical arcing, patient burns during defibrillation and may divert therapy current away from the heart.
- Do not reuse due to risk of cross infection, degradation of adhesive or electrical performance.
- Do not use on patients with a history of skin sensitization or irritation in direct electrode contact. Report any incompatibility reactions observed to Philips.
- Employ recommended procedures from electrosurgical or monitoring equipment manufacturers to minimize risk of burns at electrode site during electrosurgical procedures.
- Do not use solvents to clean skin prior to ECG. Solvents trapped under electrodes may lead to abnormal skin reactions.
- Do not apply over open wounds, lesions, infected, irritated or inflamed areas to prevent further damage of skin.
- Use caution when removing electrodes to prevent skin damage.
- Dispose of electrodes properly.
- Do not use electrodes longer than the duration indicated on the packaging pouch.
- For MRI applications, only use ECG electrodes suitable for this purpose. These are marked with the MR Conditional symbol on the packaging pouch. In this case, see “MRI Safety Information for M2202A” on page 2.
- Only the M2202A electrode is translucent in X-ray and MRI environments.
- Only the M2202A electrode is suitable for use up to and including 3.0T when used in accordance with MR equipment manufacturer instructions. Electrode temperature may increase after 15 minutes of continuous MRI scanning.
- Properly dispose electrode backing plastic liner to prevent possible patient aspiration.

Cautions

- Federal (US) law restricts this device to sale by or on the order of a medical practitioner.
- Do not use if gel is dry.
- Replace electrodes if they no longer firmly stick to patient skin.
- Use electrodes within seven days of opening pouch.
- Keep unused electrodes in original pouch which is closed by folding over.
- Ensure that application sites are properly prepared and are dry, clean and free of excessive hair. Inaccurate readings could result.
- Use care when applying electrode to avoid gel displacement. Loss of adhesion could occur.
- Do not use electrodes beyond the expiration date marked on the package.
- Use only electrodes that are suitable and appropriately labeled for X-ray imaging. These are marked with “Radiolucent” in the product description on the packaging pouch.

MRI Safety Information for M2202A

Non-clinical testing has demonstrated that the M2202A is MR Conditional. A patient with this device applied can be safely scanned in an MRI system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla, with:
 - Maximum spatial field gradient of 12,800 G/cm (128 T/m)
 - Maximum force product of 231 T²/m
 - Theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode) or 4 W/kg (First Level Controlled Operating Mode).

After 15 minutes of continuous scanning, under the scanning conditions defined above, the ECG electrodes are expected to produce a maximum temperature rise of less than:

- 1.7 °C (2 W/kg, 1.5 Tesla) or 3.4 °C (4 W/kg, 1.5 Tesla) RF-related temperature increase with a background temperature rise of approximately 1.2 °C (2 W/kg, 1.5 Tesla) or 2.5 °C (4 W/kg, 1.5 Tesla)
- 1.3 °C (2 W/kg, 3 Tesla) or 2.7 °C (4 W/kg, 3 Tesla) RF-related temperature increase with a background temperature rise of approximately 0.9 °C (2 W/kg, 3 Tesla) or 1.8 °C (4 W/kg, 3 Tesla).

In non-clinical testing, the image artifact caused by the device extends approximately 4.79 mm from the ECG electrodes when imaged with a gradient echo pulse sequence and a 3 Tesla MRI system.

Visual Inspection Prior to Use

Before applying any electrode to a patient, visually inspect the electrode. Check for any damage that could compromise accurate readings or cause patient/user injury. If visual inspection reveals any damage, follow appropriate disposal procedures. See “Disposal” on page 3.

Compatibility

The operator is responsible for checking compatibility of their Philips monitoring equipment and these electrodes. Incompatible components may result in degraded performance.

Patient Application

Apply electrodes to the patient at appropriate electrode locations that meet ECG monitoring or recording requirements (refer to the ECG instrument’s Instructions for Use for recommended electrode locations).

To apply the electrodes:

1. Prepare the electrode site according to your healthcare facility's protocol.
2. Peel the liner from each electrode as it is being applied.
3. Apply all electrodes necessary for the procedure, then connect the ECG lead set.
4. When the ECG procedure is completed remove the lead set, then peel each electrode slowly and carefully from the patient while being careful to prevent damage to the patient's skin.

Reordering Information

Part REF	Description
989803101671 (40493D)	Foam Wet Gel ECG Electrode
989803101681 (40493E)	Foam Wet Gel ECG Electrode
989803129091 (M4612A)	Solid Gel Foam Electrode
989803129101 (M4613A)	Solid Gel Foam Electrode
989803192551	Clear Tape Wet Gel ECG Electrode
989803101301 (40420A)	Hi-Tack Foam Wet Gel ECG Electrode
989803105971 (M2202A)	Radiolucent Foam Electrode

Disposal

Follow approved medical waste disposal methods as specified by your patient care facility or local regulations.

Incident Reporting

Any serious incident that has occurred in relation to this device should be reported to Philips and the Competent Authority of the European Economic Area (EEA) countries, including Switzerland and Turkey, in which the user and/or patient is established.

Specifications

	Operating, Storage, Transport
Temperature	5 °C to 30 °C (41 °F to 86 °F)
Atmospheric Pressure	570hPa to 1000hPa
Relative Humidity (non-condensing)	15% RH to 90% RH

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Copyright © 2020
All rights reserved
Nederland B.V.