
Mon-a-Therm™




General Purpose Temperature Probe 400TM

en	Instructions For Use
fr	Mode d'emploi
de	Gebrauchsanweisung
nl	Gebruiksaanwijzing
it	Istruzioni per l'uso
es	Instrucciones de uso
sv	Bruksanvisning
da	Bruksanvisning
no	Bruksanvisning
fi	Käyttöohjeet
pt	Instruções de utilização
ru	Инструкции по применению
zh	使用说明
pl	Instrukcja użytkowania
cs	Návod k použití
sk	Návod na použitie
sl	Navodilo za uporabo
hu	Használati utasítás
el	Οδηγίες χρήσης
tr	Kullanma Talimatları
ar	إرشادات الاستخدام
bg	Инструкции за употреба
ro	Instrucțiuni de utilizare
et	Kasutusjuhised
lt	Naudojimo instrukcijos
lv	Lietošanas instrukcija
hr	Upute za uporabu
sr	Uputstvo za upotrebu
mk	Упатство за употреба

REF 90044

REF 90050

5020400A 5020410A	 Diagram of cable assembly TEM_10290_A. It features a cable with a multi-pin connector on the left, a small rectangular component in the middle, and a right-angle connector on the right. TEM_10290_A
5020401A 5020411A	 Diagram of cable assembly TEM_10291_A. It features a cable with a multi-pin connector on the left, a small rectangular component in the middle, and a circular connector on the right. An arrow points to the circular connector, which has two small circles in the center. TEM_10291_A
5020405A 5020415A	 Diagram of cable assembly TEM_10292_A. It features a cable with a multi-pin connector on the left, a small rectangular component in the middle, and a circular connector on the right. An arrow points to the circular connector, which has six small circles arranged in a ring. TEM_10292_A

Mon-a-Therm™

General Purpose

Temperature Probe

400TM

REF 90044 12 Fr/Ch (4.0 mm)

REF 90050 9 Fr/Ch (3.0 mm)

Identification of a substance that is contained or present within the product or packaging.

Identification of a substance that is not contained or present within the product or packaging.

Follow instructions for use

Do not reuse. Single patient use only.

Sterile in unopened, undamaged package; sterilized in the U.S.A.

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection, or product failure risks to the patient.

This product contains DEHP. When used as indicated, very limited exposure to trace amounts of DEHP may occur. There is no clear clinical evidence that this degree of exposure increases clinical risk. However, in order to minimize risk of DEHP exposure in children and nursing or pregnant women, this product should only be used as directed.

Directions for Use

Description

The Mon-a-Therm™ General Purpose Temperature Probe is packaged individually as a sterile, single-use device and is available in a variety of sizes and sensor styles. The type of sensor and device size are designated on the unit package. The device features a Satin Slip Surface to aid in insertion. A temperature sensor is located near the distal tip. Each temperature probe electrically connects to a compatible reusable cable available from Covidien or an authorized distributor.

Indications

The Mon-a-Therm™ General Purpose Temperature Probe is indicated for use in the routine monitoring of temperature in an anesthetized patient. The device is intended primarily for insertion into the esophagus or rectum, although medical judgement may dictate the selection of other anatomical sites, such as the nasopharynx, for some patients.

Contraindications

Esophageal placement of the probe may be contraindicated in patients with upper esophageal diverticulae or stenosis, as well as in neonates or small infants undergoing operative procedures of the neck, such as tracheotomy or insertion of an internal jugular catheter. Rectal placement may also be contraindicated in neonates and small infants due to the possibility of colon perforation.

WARNINGS

- **Inadvertent tracheal or bronchial insertion may cause airway obstruction, which could be fatal if uncorrected. Invasion of the trachea may permit air to escape around the tracheal tube cuff. Do not try to seal the leak; instead, withdraw the probe and reinsert in the intended location or choose another anatomical site for placement. Probe position in the esophagus may be verified by direct laryngoscopy, if necessary.**
- **Should the scheduled procedure involve surgical dissection of neck structures, inform the surgeon if the probe has been inserted in the esophagus.**

- As with any temperature probe, use proper technique during electrosurgical procedures to help reduce the possibility that the probe will serve as an alternate path for radio frequency current and burn the surrounding tissues. At minimum, proper technique should include:
 - Locating and connecting an adequate electrosurgical dispersive electrode close to the active surgical site.
 - Avoiding draping the patient temperature monitor instrument cable over grounded surfaces or intertwining it with the cables of the electrosurgical generator.
- Do not use if package is opened or damaged.
- Use only with compatible interconnect cables and monitors. Not doing so may result in an inaccurate measured temperature.
- Do not modify this device. Modifications could result in inaccurate temperature measurement.
- Do not use this device during MRI scanning. Conducted current may cause burns.

CAUTIONS

- Intubate patient with tracheal tube before attempting to place the probe esophageally. Remove the probe prior to extubation of the tracheal tube.
- Insertion and removal should be performed using proper medical techniques.
- Do not attempt to force connection of the sensor. The connectors are keyed by shape. Forcing the connection could result in damage to the sensor or cable connector and a resulting inability to function.
- Ensure that the cable connector does not lie in an area where solutions may pool during the surgical procedure.
- The operation of some patient temperature monitor instruments may be temporarily affected during electrosurgical activations. Unusual temperature readings should be checked during a pause between activations.

Adverse Reactions

The following adverse reactions have been reported to be associated with the use of temperature probes during the insertion procedure or while the devices were in use. The order of listing is alphabetical and does not indicate frequency or severity.

Reported adverse reactions include: airway obstruction, aspiration pneumonia, bronchial insertion, electrical burns due to aberrant electrosurgical pathways, colon perforation, epistaxis, esophageal abrasion, esophageal perforation, rectal bleeding, tracheal insertion, and trauma to pharynx.

Instructions for Use

1. Carefully remove the sterile temperature probe from its protective package.
2. Lubricate the probe tip with a suitable water-soluble lubricant.
3. Connect the appropriate reusable cable by inserting the sensor connector into the mating receptacle at the end of the cable until it is firmly seated. The connectors are keyed by shape.
4. Insert the probe in the desired location for monitoring temperature using proper medical techniques. (If the probe is to be placed esophageally, intubate patient with a tracheal tube prior to insertion.)
5. Listen for airway leaks when inserting esophageally. Inadvertent invasion of the trachea may permit air to escape around the tracheal tube cuff. Do not try to seal the leak; instead, remove the probe and reinsert. Correct placement may be verified by direct laryngoscopy, if necessary.
6. Follow the directions for use for the patient temperature monitor. If problems are encountered in monitoring temperature, verify that the necessary cable connections have been made and maintained.
7. A probe placed esophageally should be removed prior to extubation of the tracheal tube. Remove the probe using proper medical techniques.
8. Disconnect the temperature probe, making sure not to pull directly on the cable wire.

Specifications

- For Direct Mode Temperature Measurement
- Thermistor: 400 Series (2,252 Ohm @ 25 °C)
- Stability Response Time: < 2 minutes
- Temperature Accuracy
 - Range: 35 °C to 42 °C
 - Accuracy: ± 0.2 °C
 - Extended Range: 25 °C to 45 °C
 - Accuracy: ± 0.3 °C

Compatible Interface Cables

See inside front cover for list and images of compatible interface cables.

Compatible Patient Monitors

All patient monitors* that meet the specifications for *Thermistor*, *Temperature Accuracy*, and *Compatible Interface Cables* as described above are compatible.

**Tested with Philips IntelliVue MP50 (Model #M8004A), M3001A Multi-Measurement Server (MMS)*

Disposal

Dispose of the temperature probe after each use in accordance with applicable regulations for biologically hazardous waste.

Storage and Shipping Conditions

- Temperature: -25 °C to 55 °C
- Humidity: 10% to 95% (non-condensing)
- Atmospheric Pressure: 12 kPa to 106 kPa

Protect from moisture.

Additional Copies of Instructions

Additional copies of these instructions can be obtained through an authorized distributor. In addition, permission is hereby granted under Covidien copyrights to purchasers of products obtained from an authorized distributor to make copies of these instructions for use by such purchasers.

Warranty

To obtain information about a warranty, if any, for this product, contact Covidien's Technical Services Department or your local Covidien representative.

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COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG.
Other brands are trademarks of a Covidien company.
Covidien Ireland Limited,