



Mon-a-Therm™

Skin Temperature Probe 400TM

REF 90045

en	Instructions For Use
fr	Mode d'emploi
de	Gebrauchsanweisung
nl	Gebruiksaanwijzing
it	Istruzioni per l'uso
es	Instrucciones de uso
sv	Bruksanvisning
da	Betjeningsvejledning
no	Bruksanvisning
fi	Käyttöohjeet
pt	Instruções de utilização
ru	Инструкции по применению
zh	使用说明
pl	Instrukcja użytkowania
cs	Návod k použití
sk	Návod na použitie
sl	Navodilo za uporabo
hu	Használati utasítás
el	Οδηγίες χρήσης
tr	Kullanma Talimatları
ar	تعليمات الاستخدام
bg	Инструкции за употреба
ro	Instrucțiuni de utilizare
et	Kasutusjuhend
lt	Naudojimo instrukcijos
lv	Lietošanas instrukcijas
hr	Upute za uporabu
sr	Uputstvo za upotrebu
mk	Упатство за употреба

5020400A 5020410A	 Diagram of cable assembly TEM_10290_A. It features a cable with a multi-pin connector on the left, a small rectangular component in the middle, and a right-angle connector on the right. TEM_10290_A
5020401A 5020411A	 Diagram of cable assembly TEM_10291_A. It features a cable with a multi-pin connector on the left, a small rectangular component in the middle, and a circular connector on the right. An arrow points from the circular connector to a detailed view of the connector's face, which has two small circular ports. TEM_10291_A
5020405A 5020415A	 Diagram of cable assembly TEM_10292_A. It features a cable with a multi-pin connector on the left, a small rectangular component in the middle, and a circular connector on the right. An arrow points from the circular connector to a detailed view of the connector's face, which has six circular ports arranged in a hexagonal pattern. TEM_10292_A

Mon-a-Therm™

Skin Temperature Probe

400TM

Identification of a substance that is contained or present within the product or packaging.

Identification of a substance that is not contained or present within the product or packaging.

Follow instructions for use

Do not reuse. Single patient use only.

Sterile in unopened, undamaged package; sterilized in the U.S.A.

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection, or product failure risks to the patient.

This product contains DEHP. The intended use limits exposure to intact skin, minimizing the risk of DEHP absorption. In order to avoid undue risk of DEHP exposure in children and nursing or pregnant women, product should only be used as directed.

Directions for Use

Description

The Mon-a-Therm™ Skin Temperature Probe is packaged individually as a sterile, single-use device and is available in a variety of sizes and sensor styles. The type of sensor is designated on the unit package. A temperature sensor is positioned between a thin film and a foam layer. The surface of the film is coated with a hypoallergenic adhesive that is protected by a release liner. The foam insulates the sensor from ambient conditions. Each temperature probe electrically connects to a compatible reusable cable available from Covidien or an authorized distributor.

Indications

The Mon-a-Therm™ Skin Temperature Probe is indicated for use in the routine monitoring of skin temperature. The device is intended to be held in place on the skin by its adhesive surface.

Contraindications

The use of a temperature probe may be contraindicated over traumatized areas.

WARNINGS

- **Avoid applying the probe to traumatized areas of skin (such as burns, cuts, lacerations, or ulcers), as the presence and removal of the probe may cause further trauma.**
- **Erroneous temperature readings may be encountered in unusually warm or cold environments. External heating sources, such as infrared lamps and electric blankets, may cause localized skin temperature increases that are not indicative of changes in body core temperature.**
- **As with any temperature probe, use proper technique during electrosurgical procedures to help reduce the possibility that the probe will serve as an alternate path for radio frequency current and burn the surrounding tissues. At minimum, proper technique should include:**
 - **Locating and connecting an adequate electrosurgical dispersive electrode close to the active surgical site.**

- **Avoiding draping the patient temperature monitor instrument cable over grounded surfaces or intertwining it with the cables of the electrosurgical generator.**
- **Do not use if package is opened or damaged.**
- **Use only with compatible interconnect cables and monitors. Not doing so may result in an inaccurate measured temperature.**
- **Do not modify this device. Modifications could result in inaccurate temperature measurement.**
- **Do not use this device during MRI scanning. Conducted current may cause burns.**

CAUTIONS

- For optimal probe adhesion and patient comfort, the skin temperature monitoring site should be free of moisture and excessive hair.
- Do not attempt to force connection of the sensor. The connectors are keyed by shape. Forcing the connection could result in damage to the sensor or cable connector and a resulting inability to function.
- Ensure that the cable connector does not lie in an area where solutions may pool during the surgical procedure.
- The operation of some patient temperature monitor instruments may be temporarily affected during electrosurgical activations. Unusual temperature readings should be checked during a pause between activations.

Adverse Reactions

The following adverse reactions have been reported to be associated with the use of temperature probes while the devices were in use. The order of listing does not indicate frequency or severity.

Reported adverse reactions include: electrical burns due to aberrant electrosurgical pathways, irritation to sensitive skin, and trauma to delicate skin when the probe is removed.

Instructions for Use

1. For optimal probe adhesion and patient comfort, select a skin temperature monitoring site that is free of moisture and excessive hair. In some cases, application of an alcohol wipe will aid in drying the site.
2. Carefully remove the sterile temperature probe from its protective package and uncoil the lead wire.
3. Connect the appropriate reusable cable by inserting the sensor connector into the mating receptacle at the end of the cable until it is firmly seated. The connectors are keyed by shape.
4. Peel off the release liner from the adhesive surface of the temperature probe and apply the probe to the patient's skin.
5. Allow about four minutes for the sensor to completely equilibrate before reading temperature.
6. Follow the directions for use for the patient temperature monitor. If problems are encountered in monitoring temperature, verify that the necessary cable connections have been made and maintained.
7. When temperature monitoring is complete, carefully peel the probe off the skin.
8. Disconnect the temperature sensor, making sure not to pull directly on the cable wire.

Specifications

- For Direct Mode Temperature Measurement
- Thermistor: 400 Series (2,252 Ohm @ 25°C)
- Stability Response Time: < 1 minute
- Ambient Temperature: 15°C to 35°C
- Temperature Accuracy
 - Range: 35°C to 42°C
 - Accuracy: $\pm 0.2^{\circ}\text{C}$
 - Extended Range: 25°C to 45°C
 - Accuracy: $\pm 0.3^{\circ}\text{C}$

Compatible Interface Cables

See inside front cover for list and images of compatible interface cables.

Compatible Patient Monitors

All patient monitors* that meet the specifications for *Thermistor*, *Temperature Accuracy*, and *Compatible Interface Cables* as described above are compatible.

**Tested with Philips IntelliVue MP50 (Model #M8004A), M3001A Multi-Measurement Server (MMS)*

Disposal

Dispose of the temperature sensor after each use in accordance with applicable regulations for biologically hazardous waste.

Storage and Shipping Conditions

- Temperature: -25°C to 55°C
- Humidity: 10% to 95% (non-condensing)
- Atmospheric Pressure: 12 kPa to 106 kPa

Protect from moisture.

Additional Copies of Instructions

Additional copies of these instructions can be obtained through your authorized distributor. In addition, permission is hereby granted under Covidien copyrights to purchasers of products obtained from an authorized distributor to make copies of these instructions for use by such purchasers.

Warranty

To obtain information about a warranty, if any, for this product, contact Covidien's Technical Services Department or your local Covidien representative.

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COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG.
Other brands are trademarks of a Covidien company.
Covidien Ireland Limited,