Nellcor™

Adult SpO₂ Sensor

en	Instructions For Use
fr	Mode d'emploi
de	Gebrauchsanweisung
nl	Gebruiksaanwijzing
es	Modo de empleo
it	Istruzioni per l'uso
sv	Bruksanvisning
pt	Instruções de uso
pl	Instrukcja do Użycia
ru	Указания по применению

en

Directions for Use

- Identification of a substance that is contained or present within the product or packaging.
- Identification of a substance that is not contained or present within the product or packaging.
- Do not use if package is opened or damaged

IP22

Protected against access to hazardous parts with a finger and protected against vertically falling water drops when sensor tilted up to 15°.

This product cannot be adequately cleaned and/or reused by the user in order to facilitate reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

This device is not made with natural rubber latex or DEHP.

Indications/Contraindications

The Nellcor™ Adult SpO₂ Sensor, model 901-A, is indicated for single patient use when used in the following settings: arterial oxygen saturation and pulse rate monitoring are required for patients weighing more than 30 kg.

Use this sensor only with Nellcor instruments and instruments containing Nellcor oximetry, or with instruments licensed to use Nellcor sensors (Nellcor-compatible instruments). Consult individual manufacturers for features and compatibility of particular instruments and sensor models.

Each manufacturer of Nellcor-compatible instruments is responsible for determining whether and under what conditions its instruments are compatible for safe and effective use with each Nellcor sensor model. This may include different specifications and/or warnings, cautions, or contraindications. Refer to the instrument operator's manual or consult manufacturer for complete instructions for use of this sensor with their Nellcor-compatible instrument.

The 901-A is contraindicated for use on patients who exhibit allergic reactions to the adhesive tape.

Instructions for Use

- 1) Remove plastic bandaging from the 901-A and locate transparent windows (a) on the adhesive side. Windows cover optical components **A**.

An index finger is the preferred 901-A location. Alternatively, apply the sensor to a small thumb, smaller finger, or great toe.


Note: When selecting a sensor site, priority should be given to an extremity free of an arterial catheter, blood pressure cuff, or intravascular infusion line.

- 2) Orient the 901-A so the dashed line in the middle of the sensor is centered on the tip of the digit **B**. Wrap adhesive flap on non-cable side around the digit. Note that the cable must be positioned on the top of the hand.

- 3) Fold cable end over top of digit so that windows are directly opposite each other. Wrap adhesive securely around sides of digit **C**.

- 4) Plug the 901-A into the oximeter and verify proper operation as described in the oximeter operator's manual.

Note: If the sensor does not track the pulse reliably, it may be incorrectly positioned—or the sensor site may be too thick, thin, or overly pigmented, or otherwise deeply colored (for example, as a result of externally applied coloring such as nail polish, dye, or pigmented cream) to permit appropriate light transmission. If any of these conditions exist, reposition the sensor or choose an alternate Nellcor sensor for use on a different site.

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Zusätzliche Exemplare dieser Gebrauchsanweisung sind kostenlos von Covidien bzw. von autorisierten Vertretern von Covidien erhältlich. Unter dem Copyright von Covidien dürfen Käufer von Produkten, die von Covidien oder zugelassenen Händlern gekauft wurden, zur eigenen Nutzung Kopien dieser Anleitung erstellen.