Nellcor™

Infant SpO₂ Sensor

en	Instructions For Use
fr	Mode d'emploi
de	Gebrauchsanweisung
nl	Gebruiksaanwijzing
es	Modo de empleo
it	Istruzioni per l'uso
sv	Svårskänning
sv	Brugsanvisning
pt	Instrução de Uso
ru	Указания по Применению

en

Directions for Use

- Identification of a substance that is contained or present within the product or packaging.
- Identification of a substance that is not contained or present within the product or packaging.
- Do not use if package is opened or damaged.

IP22

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

This device is not made with natural rubber latex or DEHP.

Indications/Contraindications

The Nellcor™ Infant SpO₂ Sensor, model 902-1, is indicated for single patient use when continuous noninvasive arterial oxygen saturation and pulse rate monitoring are required for patients weighing between 3 and 20 kg.

Use this sensor only with Nellcor instruments and instruments containing Nellcor oximetry, or with instruments licensed to use Nellcor sensors (Nellcor-compatible instruments). Consult individual manufacturers for features and compatibility of particular instruments and sensor models.


Each manufacturer of Nellcor-compatible instruments is responsible for determining whether and under what conditions its instruments are compatible for safe and effective use with each Nellcor sensor model. This may include different specifications and/or warnings, cautions, or contraindications. Refer to the instrument operator's manual or consult manufacturer for complete instructions for use of this sensor with their Nellcor-compatible instrument.

The 902-1 is contraindicated for use on patients who exhibit allergic reactions to the adhesive tape.

Instructions for Use

- Remove plastic backing from the 902-1 and locate transparent windows (b) on the adhesive side. Windows cover optical component. No exposure of the alignment marks (a) on the nonadhesive side and the dashed line (c) midway between the marks.
- Wrap the 902-1 firmly, but not too tightly around the toe.
- Wrap any excess tape loosely around toe. Use additional tape provided to secure cable across bottom of foot, loosely enough to maintain good circulation.
- Plug the 902-1 into the oximeter and verify proper operation as described in the oximeter operator's manual.

Note: If the sensor does not track the pulse reliably, it may be incorrectly positioned—or the sensor site may be too thick, thin, or deeply pigmented, or otherwise deeply colored (for example, as a result of externally applied coloring such as nail polish, dye, or pigmented cream) to permit appropriate light transmission. If any of these situations occur, reposition the sensor or choose an alternate Nellcor sensor for use on a different site.

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The 902-1 can be reused on the same patient as long as the adhesive tape attaches without slippage.

Enclosed adhesive "dots" are provided for reapplication. Place a transparent dot over each window as shown, then remove the protective paper that covers each dot. The sensor is now ready to be reapplied to the same patient.

Warnings

- Failure to apply the 902-1 properly may cause incorrect measurements.
- High oxygen levels may predispose a premature infant to develop retinopathy. Therefore, the upper alarm limit for oxygen saturation must be carefully selected in accordance with accurate pulse rate monitoring, considering the accuracy range of the oximeter being used.
- While the 902-1 is designed to reduce the effects of ambient light, excessive light may cause inaccurate measurements. In such cases, cover the sensor with opaque material.

Circulation distal to the sensor site should be checked routinely. The site must be inspected every 8 hours to ensure adhesion, skin integrity, and correct optical alignment. If skin integrity changes, move the sensor to another site.

- Intravascular dyes or externally applied coloring such as nail polish, dye, or pigmented cream may lead to inaccurate measurements.
- Excessive motion may compromise performance. In such cases, try to keep the patient still, or change the sensor site to one with less motion.
- Do not immerse in water or cleaning solutions. Do not sterilize.
- If the sensor is wrapped too tightly or supplemental tape is applied, venous pulsations may lead to inaccurate saturation measurements.
- Do not use the 902-1 or other oximetry sensors during MRI scanning. Conducted current may cause burns. Also, the 902-1 may affect the MRI image, and the MRI unit may affect the accuracy of oximetry measurements.
- Do not alter or modify the 902-1. Alterations or modifications may affect performance or accuracy.
- For additional warnings, cautions or contraindications when using this sensor with Nellcor-compatible instruments, refer to the instrument operator's manual or contact the manufacturer of the instrument.
- In the event of damage to the sterile packaging, do NOT sterilize. Follow local governing ordinances and recycling instructions regarding disposal or recycling of sensors.

A great toe is the preferred 902-1 location. Alternatively, apply the sensor to another digit of similar size, for example, a thumb.

Note: When selecting a sensor site, priority should be given to an extremity free of an arterial catheter, blood pressure cuff, or intravascular infusion line.

Orient the 902-1 so the window next to the cable is aligned on the bottom of the great toe as shown. The cable should extend towards the heel.

Wrap the 902-1 around the toe. The sensor should be wrapped around the toe so that the cable is aligned on the bottom of the great toe as shown. The cable should extend towards the heel.

Wrap any excess tape loosely around toe. Use additional tape provided to secure cable across bottom of foot, loosely enough to maintain good circulation.

Plug the 902-1 into the oximeter and verify proper operation as described in the oximeter operator's manual.

Note: If the sensor does not track the pulse reliably, it may be incorrectly positioned—or the sensor site may be too thick, thin, or deeply pigmented, or otherwise deeply colored (for example, as a result of externally applied coloring such as nail polish, dye, or pigmented cream) to permit appropriate light transmission. If any of these situations occur, reposition the sensor or choose an alternate Nellcor sensor for use on a different site.

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Copie aggiuntive di queste istruzioni sono disponibili gratuitamente contattando Covidien o i suoi distributori autorizzati. Inoltre, con il presente documento si accorda il permesso previsto dai copyright di Covidien agli Acquirenti dei prodotti ottenuti da Covidien o dai suoi distributori autorizzati di fare ulteriori copie delle istruzioni perche tali acquirenti possano utilizzarle.
COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and International registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company.
Covidien, Inc., 15 Hampshire Street, Mansfield, MA 02048 USA.
Covidien Ireland Limited, IDA Business & Technology Park, Tullamore.