en	Instructions For Use
fr	Mode d'emploi
de	Gebrauchsanweisung
nl	Gebruiksaanwijzing
it	Istruzioni per l'uso
es	Modo de empleo
pt	Instrução de Uso
sv	Brugsanvisning
da	Brugsanvisning
ru	Указания по применению
zh	使用说明

en	Instructions For Use
fr	Mode d'emploi
de	Gebrauchsanweisung
nl	Gebruiksaanwijzing
it	Istruzioni per l'uso
es	Modo de empleo
pt	Instrução de Uso
sv	Brugsanvisning
da	Brugsanvisning
ru	Указания по применению
zh	使用说明

	Identification of a substance that is contained or present within the product or packaging.
	Identification of a substance that is not contained or present within the product or packaging.
	Do not use if package is opened or damaged.
	Protected against access to hazardous parts with a finger and protected against vertically falling water drops when sensor tilted up to 15°.

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

This device is not made with natural rubber latex or DEHP.

en	Warnings
fr	Précautions
de	Vorsicht
nl	Let op
it	Avvertenze
es	Advertencias
pt	Atenção
sv	Varningar
da	Varnings
ru	Предупреждения
zh	警告

en	Warnings
fr	Précautions
de	Vorsicht
nl	Let op
it	Avvertenze
es	Advertencias
pt	Atenção
sv	Varningar
da	Varnings
ru	Предупреждения
zh	警告

en	Warnings
fr	Précautions
de	Vorsicht
nl	Let op
it	Avvertenze
es	Advertencias
pt	Atenção
sv	Varningar
da	Varnings
ru	Предупреждения
zh	警告

The Nellcor® Pediatric SpO₂ Sensor, model 903-P, is indicated for single patient use when continuous noninvasive arterial oxygen saturation and pulse rate monitoring are required for patients weighing between 10 and 50 kg.

Use this sensor only with Nellcor instruments and instruments containing Nellcor oximetry, or with instruments licensed to use Nellcor sensors (Nellcor-compatible instruments). Consult individual manufacturers for features and compatibility of particular instruments and sensor models.

Each manufacturer of Nellcor-compatible instruments is responsible for determining whether and under what conditions its instruments are compatible for safe and effective use with each Nellcor sensor model. This may include different specifications and/or warnings, cautions, or contraindications. Refer to the instrument operator's manual or consult manufacturer for complete instructions for use of this sensor with their Nellcor-compatible instrument.

The 903-P is contraindicated for use on patients who exhibit allergic reactions to the adhesive tape.

- Remove plastic backing from the 903-P and locate transparent windows (a) on the adhesive side. Windows cover optical components **A**.
An index finger is the preferred 903-P location. Alternatively, apply the sensor to a small thumb, smaller finger, or great toe.

Note: When selecting a sensor site, priority should be given to an extremity free of an arterial catheter, blood pressure cuff, or intravascular infusion line.

- Orient the 903-P so the dashed line in the middle of the sensor is centered on the tip of the digit **B**. Wrap adhesive flaps on non-cable end around the digit. Note that the cable must be positioned on the top of the hand or foot.

- Fold cable end over top of digit so that windows are directly opposite each other. Wrap adhesive securely around sides of digit **C**.

- Plug the 903-P into the oximeter and verify proper operation as described in the oximeter operator's manual.

Note: If the sensor does not track the pulse reliably, it may be incorrectly positioned—or the sensor site may be too thick, thin, or deeply pigmented, or otherwise deeply colored (for example, as a result of externally applied coloring such as nail polish, dye, or pigmented cream) to permit appropriate light transmission. If any of these situations occur, reposition the sensor or choose an alternate Nellcor sensor for use on a different site.

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For the accuracy specification range when used with Nellcor monitors, refer to information provided with the monitor, or (in the U.S., contact Covidien's Technical Services Department. Outside the U.S., contact your local Covidien representative.

For the accuracy specification range of this sensor when used with Nellcor-compatible instruments, refer to the instrument operator's manual or contact the instrument manufacturer.

Additional copies of these Instructions are available at no charge by calling Covidien or its authorized distributors. Also, permission is hereby granted under Covidien copyrights to purchasers of products obtained from Covidien or its authorized distributors to make additional copies of these instructions for use by such purchasers.

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COVIDien, COVIDien with logo and COVIDien logo are U.S. and internationally registered trademarks of Covidien LLC. Other brands are trademarks of a Covidien company.