IntelliVue MX40

Reusable Patient Cables

Instructions for Use

AAMI ECG Only

REF

AAMI 3-lead grabber.....	989803171801
AAMI 5-lead grabber.....	989803171831
AAMI 6-lead grabber.....	989803171861
AAMI 5-lead snap	989803171821

AAMI ECG + SpO₂

REF

AAMI 3-lead grabber.....	989803171811
AAMI 5-lead grabber.....	989803171851
AAMI 6-lead grabber.....	989803171871
AAMI 5-lead snap	989803171841

IEC ECG Only



REF

IEC 3-lead grabber	989803171901
IEC 5-lead grabber	989803171931
IEC 6-lead grabber	989803171961

IEC ECG + SpO₂

REF

IEC 3-lead grabber	989803171911
IEC 5-lead grabber	989803171951
IEC 6-lead grabber	989803171971

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

EN Instructions for Use

ASSOCIATED DOCUMENT

Refer to the *Cables and Lead Sets Care, Cleaning and Disinfection* document that accompanies this Instructions for Use (IFU). It provides validated cleaning/disinfecting agents and cleaning procedures for these cable products.

INTENDED USE

Philips IntelliVue MX40 Reusable ECG and ECG + SPO₂ Cables are intended for monitoring and recording of and to generate alarms for, multiple physiological parameters of adults and pediatrics in a hospital environment and during patient transport inside hospitals. Not intended for home use. Intended for use by health care professionals.

Indications for Use

Philips IntelliVue MX40 Reusable ECG and ECG + SPO₂ Cables are indicated for continuous monitoring of cardiac signals and pulse oximetry for both diagnostic and monitoring purposes. These devices are limited by the indications for use of the connected monitoring and diagnostic equipment. These devices are intended to interact with patient intact skin only.

Contraindications

There are no known contraindications.

PRODUCT DESCRIPTIONS

The MX40 Reusable Cables are available in ECG only or ECG/SpO₂ configurations as identified below.

AAMI Color-Coded			IEC Color-Coded		
REF	ECG Leads*	SpO ₂ Connector	Part Number	ECG Leads*	SpO ₂ Connector
989803171801	3-lead grabbers	-	989803171901	3-lead grabbers	-
989803171831	5-lead grabbers	-	989803171931	5-lead grabbers	-
989803171861	6-lead grabbers	-	989803171961	6-lead grabbers	-
989803171821	5-lead snaps	-			
989803171811	3-lead grabbers	8-pin female	989803171911	3-lead grabbers	8-pin female
989803171851	5-lead grabbers	8-pin female	989803171951	5-lead grabbers	8-pin female
989803171871	6-lead grabbers	8-pin female	989803171971	6-lead grabbers	8-pin female
989803171841	5-lead snaps	8-pin female			

*All ECG leads are 0.85m (34 in) long.

COMPATIBILITY

These cables and leads are compatible with Philips MX40 monitors and any other monitor/defibrillator that lists them as a compatible accessory in that product's IFU. Incompatible components may result in performance issues.

WARNINGS

- Before use, read and understand all warnings listed in this Instruction for Use, as well as the *IntelliVue MX40 Instructions for Use*. Otherwise patient injury could result.
- Do not use any lead set or cable if visual inspection reveals liquid/lint contaminants in cable connectors, or other visual cable damage.
- Always clean/disinfect patient cables when they are removed from their packaging for first-time use.
- Refer to the MX40 monitor IFU for information regarding proper lead/electrode placement that complies with standard AAMI or IEC practices.
- If patient cable includes a SpO₂ sensor connector, follow the sensor application instructions supplied in SpO₂ sensor's Instructions for Use.
- To avoid patient injury, ensure that cables are carefully positioned to avoid entanglement, choking, and strangulation.
- Patient must be properly grounded during electrosurgery procedures to prevent patient/user injury (i.e. electrical shock).
- These cables are not suitable for use in an MRI environment which could create a potential for patient burns or inaccurate measurements.
- Ensure that the patient cable connector is properly and securely connected to the MX40 patient monitor. Otherwise, incorrect patient physiological data may result.
- When switching between an ECG/SpO₂ patient cable and the trunk cable, the SpO₂ sensor must be moved separately and connected to the monitor.

CAUTIONS

- Federal (US) law restricts this device to sale by or on the order of a medical practitioner.
- Do not use in excessively wet environments or under massive influence of fluids (e.g. rain).
- Do not immerse the cable connectors in any fluid.
- Do not autoclave or use ultrasonic cleaners on the cable.
- Do not clean cable connector electrical contacts or connectors with bleach.
- Expected service life of these products is 18 months of typical clinical use.
- Do not disconnect a cable by pulling on the lead wires.
- Be sure to always visually inspect the cable before connecting it to the MX40 monitor, See *Visual Inspection Prior to Use* below.

VISUAL INSPECTION PRIOR TO USE

Before attaching any cable to a patient or monitor, visually inspect the cable's physical condition, connector gasket, and connector contacts as instructed below.

Inspect Cable Physical Condition

Before attaching any cable to a patient or monitor, perform a visual inspection to determine if the cable has reached end-of-life. Check for cracks, blisters, peeling, exposed wires, damaged connectors and similar wear or damage that could compromise accurate readings or cause patient/user injury (e.g. cuts). Whenever visual inspection reveals a cable is no longer suitable for continued use, follow appropriate product disposal procedures (see *PRODUCT DISPOSAL*).

Inspect Cable Connector Gasket

Before connecting a cable to the MX40 monitor, visually ensure that the gasket surrounding the mouth of the cable connector is seated properly (reference the graphic at the front of this document that shows proper gasket seating). If any part of the gasket is dislodged or damaged **DO NOT USE!** Instead, discard the cable as there is a potential for incorrect readings and/or patient/user injury.

Inspect Cable Connector Contacts

Note: A cable's SpO₂ sensor connector should remain covered with its tethered protective cap when not in use.

The cable's monitor connector has very small contacts that must be kept clean to establish good connection between cable and monitor. Therefore, before connecting a cable to the monitor, visually confirm there is no lint, moisture, or residue inside the cable connector. If the visual check indicates a need to clean *inside* the cable connector follow the steps below.

Important! Use the following steps only to clean *inside* the cable MX40 connector. To clean connector or cable external surfaces, reference the validated cable cleaning procedures provided in the *Cables and Lead Sets Care, Cleaning and Disinfection* document that accompanies this IFU.

To clean *inside* the cable's MX40 connector:

1. Prevent cleaning fluid from seeping into the cable connector and be careful not to dislodge the gasket that surrounds the inside of the connector.
2. Keep the cable connector opening pointed down as shown in the *Clean Connector Pins* graphic at the front of this IFU; i.e. this helps prevent cleaner residue from building-up inside the connector.
3. Use only Isopropyl alcohol on a modestly dampened lint-free cloth, or pre-moistened Isopropyl alcohol wipe.
4. Gently wipe inside the cable's MX40 connector using very light swipes.
5. Allow the connector pins to air dry prior to use.

CONNECTING AND DISCONNECTING THE MX40

Connecting: Use a “straight-on” pushing motion between connectors as shown in the *Connect Cable* graphic at the front of this document. A slight clicking sound signifies a secure connection.

Disconnecting: Grip onto the sides of the connector and use a slight pivoting or rocking motion to separate the connector from the device as shown in the *Disconnect Cable* graphic.

Caution: *To avoid damaging a cable, NEVER disconnect by holding or pulling on the connector wires as shown in the Never Pull Wires graphic.*

PATIENT APPLICATION

Refer to the MX40 monitor IFU for information regarding proper ECG lead/electrode placement that complies with standard AAMI or IEC practices.

Refer to the appropriate SpO₂ sensor IFU for proper patient application of any SpO₂ sensor being connected to a ECG/SpO₂ adapter cable.

CABLE CLEANING AND DISINFECTION

Refer to the document titled *Cables and Lead Sets Care, Cleaning, and Disinfection* that accompanies this IFU for validated cleaning/disinfecting agents and procedures. To reduce risk of cross-contamination:

- Always clean/disinfect patient cables when they are removed from their packaging for first-time use.
- Clean/disinfect reusable cables before applying them to a different patient.
- Clean/disinfect reusable cables that are applied to any patient, if they become visibly soiled.

CABLE STERILIZATION

All cable products described in this IFU are not intended to be sterilized.

PRODUCT DISPOSAL

Follow approved medical waste disposal methods as specified by your patient care facility or local regulations.

INCIDENT REPORTING

Any serious incident that has occurred in relation to this device should be reported to Philips and the competent authority of the European Economic Area (EEA) countries, including Switzerland and Turkey in which the user and/or patient is established.

ENVIRONMENTAL SPECIFICATIONS

Temperature and humidity ranges identified below must be maintained for all of the cable products described in this IFU. Otherwise product damage could result.

Temperature

- Operating Temperature Range: 0 to +37°C (32 to 99 °F) at 20% R.H. to 95% R.H.
- Storage Temperature Range: - 40 to +60°C (- 40 to 140 °F) at 15% R.H. to 90% R.H.

Humidity

- Minimum Operating: Humidity at +50°C; 15% R.H. 24hrs
- Maximum Operating: Humidity at +40°C; 95% R.H. 24hrs

Atmospheric Pressure

- Operating: 670 hPa to 1080 hPa
- Storage: 100 hPa to 1080 hPa

REORDERING INFORMATION

The products listed in the table below are CE marked and available in the European Union.

REF Part Number	Description
989803171801	ECG 3 LEAD GRABBER, AAMI, TELE
989803171831	ECG 5 LEAD GRABBER, AAMI, TELE
989803171861	ECG 6 LEAD GRABBER, AAMI, TELE
989803171821	ECG 5 LEAD SNAPS, AAMI, TELE
989803171811	ECG 3 LEAD GRABBER, AAMI + SPO2, TELE
989803171851	ECG 5 LEAD GRABBER, AAMI + SPO2, TELE
989803171871	ECG 6 LEAD GRABBER, AAMI + SPO2, TELE
989803171841	ECG 5 LEAD SNAPS, AAMI + SPO2, TELE
989803171901	ECG 3 LEAD GRABBER, IEC, TELE
989803171931	ECG 5 LEAD GRABBER, IEC, TELE
989803171961	ECG 6 LEAD GRABBER, IEC, TELE
989803171911	ECG 3 LEAD GRABBER, IEC + SPO2, TELE
989803171951	ECG 5 LEAD GRABBER, IEC + SPO2, TELE
989803171971	ECG 6 LEAD GRABBER, IEC+SPO2, TELE

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[REF] 부품 번호	설명
989803171801	ECG 3 LEAD GRABBER, AAMI, TELE
989803171831	ECG 5 LEAD GRABBER, AAMI, TELE
989803171861	ECG 6 LEAD GRABBER, AAMI, TELE
989803171821	ECG 5 LEAD SNAPS, AAMI, TELE
989803171811	ECG 3 LEAD GRABBER, AAMI + SPO2, TELE
989803171851	ECG 5 LEAD GRABBER, AAMI + SPO2, TELE
989803171871	ECG 6 LEAD GRABBER, AAMI + SPO2, TELE
989803171841	ECG 5 LEAD SNAPS, AAMI + SPO2, TELE
989803171901	ECG 3 LEAD GRABBER, IEC, TELE
989803171931	ECG 5 LEAD GRABBER, IEC, TELE
989803171961	ECG 6 LEAD GRABBER, IEC, TELE
989803171911	ECG 3 LEAD GRABBER, IEC + SPO2, TELE
989803171951	ECG 5 LEAD GRABBER, IEC + SPO2, TELE
989803171971	ECG 6 LEAD GRABBER, IEC+SPO2, TELE

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[REF] Referência	Descrição
989803171801	ECG 3 LEAD GRABBER, AAMI, TELE
989803171831	ECG 5 LEAD GRABBER, AAMI, TELE
989803171861	ECG 6 LEAD GRABBER, AAMI, TELE
989803171821	ECG 5 LEAD SNAPS, AAMI, TELE
989803171811	ECG 3 LEAD GRABBER, AAMI + SPO2, TELE
989803171851	ECG 5 LEAD GRABBER, AAMI + SPO2, TELE
989803171871	ECG 6 LEAD GRABBER, AAMI + SPO2, TELE
989803171841	ECG 5 LEAD SNAPS, AAMI + SPO2, TELE
989803171901	ECG 3 LEAD GRABBER, IEC, TELE
989803171931	ECG 5 LEAD GRABBER, IEC, TELE
989803171961	ECG 6 LEAD GRABBER, IEC, TELE
989803171911	ECG 3 LEAD GRABBER, IEC + SPO2, TELE
989803171951	ECG 5 LEAD GRABBER, IEC + SPO2, TELE
989803171971	ECG 6 LEAD GRABBER, IEC+SPO2, TELE

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Copyright © 2020
All rights reserved
Nederland B.V.