IntelliVue MX40 Single-Patient-Use Cables

Instructions for Use


ECG Only

AAMI 5-lead grabber

REF 989803172031

IEC 5-lead grabber

REF 989803172131


ECG + SpO₂

AAMI 5-lead grabber plus SpO₂

REF 989803172051


IEC 5-lead grabber plus SpO₂

REF 989803172151


PHILIPS

(EN) Instructions for Use	1	(MK) Упатство за користење	119
(BG) Указания за употреба	7	(NL) Gebruiksaanwijzing	126
(CS) Návod k obsluze.....	14	(NO) Brukerhåndbok	132
(DA) Brugervejledning.....	20	(PL) Instrukcja obsługi	138
(DE) Gebrauchsanweisung	26	(PT-BR) Instruções de uso	144
(EL) Οδηγίες χρήσης	32	(PT-PT) Instruções de utilização ...	150
(ES) Instrucciones de uso	39	(RO) Instrucțiuni de utilizare....	156
(ET) Kasutusjuhend	45	(RU) Инструкция по эксплуатации.....	162
(FI) Käyttöohjeet	51	(SK) Návod na použitie	169
(FR) Manuel d'utilisation	57	(SL) Navodila za uporabo.....	175
(HR) Upute za uporabu.....	64	(SR) Uputstvo za upotrebu	181
(HU) Használati útmutató.....	70	(SV) Bruksanvisning	187
(ID) Petunjuk Penggunaan.....	76	(TR) Kullanım Talimatları.....	193
(IT) Istruzioni d'uso	82	(UK) Інструкції з використання	199
(JA) ユーザーズガイド	88	(VI) Hướng Dẫn Sử Dụng.....	206
(KK) Пайдалану нұсқаулары	94	(ZH-CN) 使用说明书	212
(KO) 사용 설명서	101	(ZH-TW) 使用說明.....	217
(LT) Naudojimo instrukcija ...	107		
(LV) Lietošanas instrukcija....	113		

EN	Do not use if gasket around connector opening is damaged	JA	コネクタ開口部に嵌め込まれたガスケットが破損している場合は使用しないでください。
BG	Не използвайте, ако уплътнението около отвора на конектора е повредено	KK	Конектор саңылауының айналасындағы тығыздағыш зақымдалған болса, оны пайдаланбаңыз
CS	Nepoužívejte, pokud je těsnění kolem otvoru konektoru poškozené	KO	커넥터 입구 주변의 가스킷이 손상된 경우 사용하지 마십시오.
DA	Må ikke anvendes, hvis pakningen omkring konnektorråbningen er beskadiget	LT	Nenaudokite, jei pažeistas aplink jungties angą esantis tarpiklis
DE	Nicht verwenden, wenn die Dichtung rund um die Anschlussöffnung beschädigt ist	LV	Neizmantojiet, ja savienotāja atvēruma aptverošā starplika ir bojāta.
EL	Μην χρησιμοποιείτε αν η φλάντζα γύρω από το άνοιγμα του συνδέσμου έχει υποστεί ζημιά	MK	Немојте да го користите ако заптивката околу отворот на приклучокот е оштетена
ES	No utilizar si la junta de la abertura del conector presenta daños	NL	Niet gebruiken als de pakking rond de connectoropening beschadigd is
ET	Mitte kasutada, kui pistiku ava ümbritsev tihend on kahjustatud	NO	Må ikke brukes hvis pakningen rundt kontaktåpningen er skadet
FI	Älä käyttää, jos liittimen aukon ympärillä oleva tiiviste on vahingoittunut	PL	Nie używać w przypadku, gdy uszczelka wokół otworu złącza jest uszkodzona
FR	Ne pas utiliser si le joint autour de l'ouverture du connecteur est endommagé	PT-BR	Não use se a junta ao redor da abertura do conector estiver danificada
HR	Nemojte upotrebljavati ako je brtva oko otvora priključka oštećena	PT-PT	Não utilizar se a junta em redor da abertura do conector estiver danificada
HU	Ne használja, ha a csatlakozónylás körüli tömítés sérült	RO	Nu utilizați dacă garnitura din jurul deschiderii conectorului este deteriorată
ID	Jangan gunakan jika gasket di sekeliling bukaan konektor rusak	RU	Не используйте, если уплотнительная прокладка вокруг отверстия разъема повреждена
IT	Non utilizzare se la guarnizione intorno all'apertura del connettore è danneggiata	SK	Nepoužívať, ak je tesnenie kalibračnej komory okolo otvoru konektora poškodené

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EN Instructions for Use

INTENDED USE

Philips IntelliVue MX40 Single-Patient-Use ECG and ECG + SpO₂ cables are intended for monitoring and recording of and to generate alarms for, multiple physiological parameters of adults and pediatrics in a hospital environment and during patient transport inside hospitals. Not intended for home use. Intended for use by health care professionals

Indications for Use

Philips IntelliVue MX40 Single Patient Use ECG and ECG + SpO₂ cables are indicated for continuous monitoring of cardiac signals, or pulse oximetry, for diagnostic and monitoring purposes. These devices are limited by the indications for use of the connected monitoring and diagnostic equipment. These devices are intended to interact with patient intact skin only.

Contraindications

There are no known contraindications.

PRODUCT DESCRIPTIONS

The MX40 Single-Patient-Use Cables are available in the ECG only or ECG/SpO₂ configurations identified below.

5-Lead ECG Only		5-Lead ECG with attached SpO ₂ Sensor*		
REF	ECG Color Coding	REF	ECG Color Coding	SpO ₂ Sensor Type
989803172031	AAMI color coded grabber connectors	989803172051	AAMI color coded grabber connectors	Finger wrap sensor
989803172131	IEC color coded grabber connectors	989803172151	IEC color coded grabber connectors	Finger wrap sensor

*The preferred SpO₂ patient application site is any finger except a thumb.

COMPATIBILITY

These cables and leads are compatible with Philips MX40 monitors and any other monitor/defibrillator that lists them as a compatible accessory in that product's IFU. Incompatible components may result in performance issues.

Warnings

- Before use, read and understand all warnings listed in this Instruction for Use, as well as the *IntelliVue MX40 Instructions for Use*. Otherwise patient injury could result.
- Do not use any lead set or cable if visual inspection reveals liquid/lint contaminants in cable connectors, or other visual cable damage.
- These patient cables are for single patient use only. They should never be reused on a different patient due to risk of cross-infection.

- To prevent risk of contamination or infection, unpack the SpO₂ sensor just prior to use.
- To avoid patient injury, ensure that cables are carefully positioned to avoid entanglement, choking, and strangulation.
- Patient must be properly grounded during electrosurgery procedures to prevent patient or user injury (i.e. electrical shock).
- These cables are not suitable for use in an MRI environment as there is potential for patient burns or inaccurate measurements.
- Ensure that the patient cable connector is properly and securely connected to the MX40 patient monitor. Otherwise, incorrect patient physiological data may result.
- Cover the SpO₂ sensor with opaque material under conditions of strong or excessive light (infrared lamps, OR lamps, photo therapy). Failure to do so can result in inaccurate measurements
- To avoid venous pulsation, obstructed circulation, pressure marks, pressure necrosis, artifacts and inaccurate measurements, make sure you are using a sensor of the correct size and that the sensor is not too tight. If the sensor is too tight, because the application site is too large or becomes too large due to edema, the excessive pressure applied may result in venous congestion distal from the application site, leading to interstitial edema and tissue ischemia.
- At elevated ambient temperatures, patient skin could be severely burned after prolonged SpO₂ sensor application at sites that are not well perfused. To prevent this condition, be sure to check patient application sites frequently. This SpO₂ sensor operates without risk of exceeding 41°C on the skin if the initial skin temperature does not exceed 35°C.
- If the SpO₂ sensor is applied too loosely it may fall off or compromise proper alignment of sensor optics and result in inaccurate readings.
- Avoid SpO₂ sensor sites subject to excessive motion. Try to keep the patient still, or move the sensor to a site with less motion.
- Dysfunctional hemoglobin or intravascular dyes can cause inaccurate SpO₂ measurements.
- When switching between an ECG/SpO₂ patient cable and the trunk cable, the SpO₂ sensor must be moved separately and connected to the monitor.

CAUTIONS

- Federal (US) law restricts this device to sale by or on the order of a medical practitioner.
- These single-patient-use products are not intended to be disinfected or sterilized.
- Do not immerse the cable connectors in any fluid.
- Do not use in excessively wet environments or under massive influence of fluids (e.g. rain).
- Ensure that leads are securely connected to the trunk cable, and the trunk cable is properly plugged into the monitor/defibrillator.
- Do not autoclave or use ultrasonic cleaners on the cable.
- Do not clean cable connector electrical contacts or connectors with bleach.
- Expected service life of these single-patient-use products is 5 days of typical clinical use.
- Be sure to always visually inspect the cable before connecting it to the MX40 monitor, See *Visual Inspection Prior to Use* below.

VISUAL INSPECTION PRIOR TO USE

Before attaching any cable to a patient or monitor, visually inspect the cable's physical condition, connector gasket, and connector contacts as instructed below.

Inspect Cable Physical Condition

Before attaching any cable to a patient or monitor, perform a visual inspection to verify that the cable is not damaged. Check for cracks, blisters, peeling, exposed wires, damaged connectors and similar wear or damage that could compromise accurate readings or cause patient/user injury (e.g. cuts). If visual inspection reveals cable damage, follow appropriate product disposal procedures (see *PRODUCT DISPOSAL*).

Inspect Cable Connector Gasket

Before connecting a cable to the MX40 monitor, visually ensure that the gasket surrounding the mouth of the cable connector is seated properly (reference the graphic at the front of this document that shows proper gasket seating). If any part of the gasket is dislodged or damaged **DO NOT USE!** Instead, discard the cable as there is a potential for incorrect readings and/or patient/user injury.

Inspect Cable Connector Contacts

The cable's monitor connector has very small contacts that must be kept clean to establish good connection between cable and monitor. Therefore, before connecting a cable to the monitor, visually confirm there is no lint, moisture, or residue inside the cable connector. If the visual check indicates a need to clean *inside* the cable connector follow the steps below.

Important! Cables described in this IFU are single-patient-use and should never be cleaned or disinfected for use on another patient. Use the following steps only to clean *inside* the cable MX40 connector if visual inspection indicates the connector contacts need cleaning.

To clean *inside* the cable's MX40 connector:

1. Prevent cleaning fluid from seeping into the cable connector and be careful not to dislodge the gasket that surrounds the inside of the connector.
2. Keep the cable connector opening pointed down as shown in the *Clean Connector Pins* graphic at the front of this IFU; i.e. this helps prevent cleaner residue from building-up inside the connector.
3. Use only Isopropyl alcohol on a modestly dampened lint-free cloth, or pre-moistened Isopropyl alcohol wipe.
4. Gently wipe inside the cable's MX40 connector using very light swipes.
5. Allow the connector pins to air dry prior to use.

CONNECTING AND DISCONNECTING THE MX40

Connecting: Use a “straight-on” pushing motion between connectors as shown in the *Connect Cable* graphic at the front of this document. A slight clicking sound signifies a secure connection.

Disconnecting: Grip onto the sides of the connector and use a slight pivoting or rocking motion to separate the connector from the device as shown in the *Disconnect Cable* graphic.

Caution: *To avoid damaging a cable, NEVER disconnect by holding or pulling on the connector wires as shown in the Never Pull Wires graphic.*

PATIENT APPLICATION, ECG LEADS

Refer to the MX40 monitor IFU for information regarding proper ECG lead/electrode placement that complies with standard AAMI or IEC practices.

PATIENT APPLICATION, SpO₂ SENSOR

When using a 5-Lead AAMI or IEC cable having an integrated SpO₂ wrap sensor (cable product numbers 989803172051 or 989803172151) apply the sensor to the patient as follows:

1. The SpO₂ Sensor is for use on a pediatric or adult fingertip. If possible, choose a finger on an extremity free from arterial catheter, blood pressure cuff, or intravascular infusion line. Select the most appropriate finger for the sensor size.
2. Gently open the sensor and guide it over the fingertip, with the sensor cable on *top* of the finger and hand. The fingertip should touch the end of the wide section of the sensor. You may need to trim the patient's fingernail to position the sensor correctly.
3. Pull the white backings from the sensor's adhesive tape, and then wrap the tape around both the sensor and cable (do not stretch the tape).
4. Move the sensor application site every four hours, or more often if circulation or skin integrity is compromised.
5. If the sensor adhesive tape becomes compromised, use any medical grade tape to reapply the sensor securely around the fingertip.
6. The sensor may be reapplied to *the same patient* by repeating the above procedure.

PRODUCT DISPOSAL

Follow approved medical waste disposal methods as specified by your patient care facility or local regulations.

INCIDENT REPORTING

Any serious incident that has occurred in relation to this device should be reported to Philips and the competent authority of the European Economic Area (EEA) countries, including Switzerland and Turkey in which the user and/or patient is established.

ENVIRONMENTAL SPECIFICATIONS

Temperature and humidity ranges identified below must be maintained for all of the cable products described in this IFU. Otherwise product damage could result.

Temperature

- Operating Temperature Range: 0 to +37°C (32 to 99 °F) at 20% R.H. to 95% R.H.
- Storage Temperature Range *for ECG only cable products* 989803172031 and 989803172031: - 40 to +60°C (- 40 to 140 °F) at 15% R.H. to 90% R.H.
- Storage Temperature Range *for cable products with ECG leads and integrated SpO₂ sensor products* 989803172051 and 989803172151: 12°C to 35°C (53.6°F to 95°F) at 18% r.H. to 78% R.H.

Humidity

- Minimum Operating: Humidity at +50°C; 15% R.H. 24hrs
- Maximum Operating: Humidity at +40°C; 95% R.H. 24hrs

Atmospheric Pressure

- Operating: 670 hPa to 1080 hPa
- Storage: 100 hPa to 1080 hPa

SpO₂ SENSOR TECHNICAL SPECIFICATIONS

The following specifications apply only to products that have a SpO₂ finger wrap sensor integrated into the lead set; i.e. product REF 989803172051 and REF 989803172151 .

SpO₂ Sensor and Pulse Rate Accuracy

For information on SpO₂ and Pulse Rate accuracy specifications, refer to your monitoring instrument's Instructions for Use; e.g. the Philips MX40 monitor's IFU.

SpO₂ Light Emitting Diode Wavelength Ranges

Wavelength ranges for the light emitting diodes used in the SpO₂ sensor are within 600nm-1000nm, with an optical output power of less than 15mW. Knowing the wavelength range can be useful to clinicians performing photodynamic therapy.

SpO₂ Measurement Validation

The SpO₂ accuracy has been validated in human studies against arterial blood sample reference measured with a CO-oximeter. In a controlled desaturation study, healthy adult volunteers with saturation levels between 70% and 100% SaO₂ were studied. The population characteristics for those studies were:

- Approximately 50% female and 50% male ranging in age from 18-45
- Skin tone: from light to black

Because pulse oximeter equipment measurements are statistically distributed, only approximately 2/3 of pulse oximeter equipment measurements can be expected to fall within the $\pm$ Arms value measured by a CO-oximeter. Functional testers, such as an SpO₂ simulator cannot be used to assess the accuracy of pulse oximeter sensors.

REORDERING INFORMATION

The products listed in the table below are CE marked and available in the European Union.

REF Part Number	Description
989803172031	ECG 5 LEAD GRABBER, AAMI, SPU, TELE
989803172131	ECG 5 LEAD GRABBER IEC, SPU, TELE
989803172051	ECG 5LEAD GRABBER, AAMI+SPO2,SPU,TELE
989803172151	ECG 5 LEAD GRABBER IEC+SPO2 SPU TELE

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Copyright © 2020
All rights reserved
Nederland B.V.
Boeblingen GmbH