
Microstream™ Advance

Neonatal-Infant Nasal CO₂ Filter Line

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




EN Instructions for Use

Intended Use

The Microstream™ Advance Neonatal-Infant Nasal CO₂ Filter Line is intended to conduct a sample of the subject's breathing to a gas measurement device (capnograph) while simultaneously administering supplemental oxygen projected near the nose and mouth for inhalation. The device is to be used with monitors using Microstream™ technology. The product is intended for single patient use only.

Intended Users

The Microstream™ Advance Neonatal-Infant Nasal CO₂ Filter Line is intended to be used by licensed healthcare providers only.

Indications for Use

The Microstream™ Advance Neonatal-Infant Nasal CO₂ Filter Line is indicated for measurement of exhaled carbon dioxide (CO₂) for monitoring purposes. These devices are limited by the indications for use of the connected Philips monitoring equipment in healthcare facilities. These devices are intended to interact with patient intact skin and nasal mucosa.

Contraindications

There are no known contraindications.

WARNINGS

- Loose or damaged connections may compromise ventilation or cause an inaccurate measurement of respiratory gases. Securely connect all components, turning the sampling line connector clockwise into the monitor CO₂ port until it can no longer be turned. Check connections for leaks according to standard clinical procedures, and ensure CO₂ values appear.
- Do not attempt to clean, disinfect, sterilize or flush any part of the sampling line. Re-use of single-use accessories could pose a cross-contamination risk to the patient or damage the functioning of the monitor.
- Inaccurate CO₂ readings or ineffective O₂ delivery may occur if the patient has clogged nares.
- For 989803204421 only: Sampling lines are recommended for use with oxygen provided at up to 3l/min. At higher levels of oxygen provision, dilution of CO₂ readings may occur, leading to lower CO₂ values.
- Ensure no kinks are present in CO₂ and O₂ tubing, as kinked tubing may cause inaccurate CO₂ sampling or affect O₂ delivery.
- The products covered in this IFU can expose you to chemicals which are known to the State of California to cause cancer, birth defects, and other reproductive harm. For more information, go to www.P65Warnings.ca.gov.
- Carefully route patient cabling (CO₂ sampling line) to reduce the possibility of patient entanglement or strangulation.
- Use of tools that can generate sparks in highly oxygenated environment can lead to patient injury.

Caution

Federal (US) law restricts this device to sale by or on the order of a medical practitioner.

Required Compatible Components

These devices are to be used with Philips monitors using Microstream™ technology.

Applying Neonatal/Infant Nasal CO₂ Filter Line

When performing the following procedure, reference the step-by-step patient application graphics located on Page i.

- ① Place cannula in patient nares.
- ② For 989803204421 only: CO₂ Nasal Sampling from one nare, O₂ provision through other nare.
- ③ Fit tubing around ears and tighten at neck.
- ④ Attach O₂ connector to O₂ source.
- ⑤ Attach the CO₂ connector to monitor.

Operational Notes

- The response time for 2m sampling lines and 3m sampling lines is ≤ 4.0 seconds; for 4m sampling lines, it is ≤ 5.7 seconds.
 - Sampling lines with Extended Duration in their names include a moisture reduction component (Nafion®* or its equivalent) for use in higher humidity environments where long duration use of CO₂ sampling is required. During use, do not cover the moisture reduction component on these lines (situated at the patient's chest area) as covering this component may lead to faster occlusion of the sampling line.
 - Verify that CO₂ readings appear on monitor screen as expected after insertion.
- TM* Trademark of its respective owner.*

Reordering Information

Part Number	Description	Quantity
989803204441	Neonatal-Infant Nasal CO2 Filter Line, Short-term use: Procedural/Emergency Note: Length: 2m	25
989803204631	Neonatal-Infant Nasal CO2 Filter Line, Extended Duration Note: Length: 2m	25
989803204431	Neonatal-Infant Nasal CO2 Filter Line, Extended Duration Note: Length: 4m	25
989803204421	Neonatal-Infant Nasal CO2 Filter Line with O2 Tubing, Extended Duration Note: Length: 3m	25

Disposal

To avoid contaminating or infecting personnel, the service environment, or other equipment, decontaminate and dispose in accordance with your country’s laws for medical waste.

Incident Reporting

Any serious incident that has occurred in relation to this device should be reported to Philips and the competent authority of the European Economic Area (EEA) countries, including Switzerland and Turkey, in which the user and/or patient is established.

Environmental Specifications

The following table provides the temperature and humidity ranges that must be maintained for these products during storage and transport:

Specification	Storage	Transport
Temperature	-20°C to 50°C (-4°F to 122°F)	-20°C to 50°C (-4°F to 122°F)
Humidity	10% RH to 90% RH	10% RH to 90% RH

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Copyright © 2020
Microstream™ is a trademark of Oridion
Nederland B.V.
Böblingen GmbH
Medical 1987 Ltd. Oridion Medical 1987
Ltd. is a Medtronic company.