
Microstream™ Advance Adult-Intermediate Bite Block CO₂ Filter Line

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








EN Instructions for Use

Intended Use

The Microstream™ Advance Adult-Intermediate Bite Block CO₂ Filter Line product is intended to sample CO₂ and administer supplemental oxygen, for patients who can wear a 60 fr. bite block, during upper endoscopy type procedures. The set is intended for single patient use only.

Intended Users

The Microstream™ Advance Adult-Intermediate Bite Block CO₂ Filter Line is intended to be used by licensed healthcare providers only.

Indications for Use

The Microstream™ Advance Adult-Intermediate Bite Block CO₂ Filter Line is indicated for measurement of exhaled carbon dioxide (CO₂), and administration of O₂, for monitoring purposes in adult patients. These devices are limited by the indications for use of the connected Philips monitoring equipment in healthcare facilities. These devices are intended to interact with patient intact skin and nasal/oral mucosa.

Contraindications

There are no known contraindications.

WARNINGS

- Loose or damaged connections may compromise ventilation or cause an inaccurate measurement of respiratory gases. Securely connect all components, turning the sampling line connector clockwise into the monitor CO₂ port until it can no longer be turned. Check connections for leaks according to standard clinical procedures, and ensure CO₂ values appear.
- Do not attempt to clean, disinfect, sterilize or flush any part of the sampling line. Re-use of single-use accessories could pose a cross-contamination risk to the patient or damage the functioning of the monitor.
- Inaccurate CO₂ readings may occur if the patient has clogged nares.
- The product is recommended for use with oxygen provided at up to 10l/min, or 5l/min with the sampling line alone. At higher levels of oxygen provision, dilution of CO₂ readings may occur, leading to attenuated waveforms and lower CO₂ values.
- Ensure no kinks are present in CO₂ and O₂ tubing, as kinked tubing may cause inaccurate CO₂ sampling or affect O₂ delivery.
- Failure to correctly insert the oral prong into the bite block channel may cause an inaccurate measurement of respiratory gases. The oral prong must remain in the channel throughout the procedure.
- Carefully route patient cabling (CO₂ sampling line) to reduce the possibility of patient entanglement or strangulation.
- Use of tools that can generate sparks in highly oxygenated environment can lead to patient injury.

Cautions

- Federal (US) law restricts this device to sale by or on the order of a medical practitioner.
- Partial blockage of the oral airway due to endoscope positioning or during suction may cause periods of low readings and rounded waveforms. The occurrence will be more pronounced with high oxygen delivery levels.

Required Compatible Components

For use with all Philips Microstream™-enabled monitors. Each product includes:

- 1 Adult-Intermediate bite block
- 1 Microstream™ Advance oral/nasal breath sampling line with O₂ delivery.

Patient Application

When performing the following procedure, reference the step-by-step patient application graphics located on Pages i to iii:

- ① Place cannula on patient.
- ② Place and adjust tubing around patient ears.
- ③ Place bite block on patient and adjust bite block strap.
- ④ Attach bite block O₂ connector to O₂ tubing.
- ⑤ Position oral prong to fit into bite block groove.
- ⑥ Attach O₂ tubing to O₂ source.
- ⑦ Plug CO₂ connector into monitor.

Operational Notes

- The response time of the 2m sampling line is ≤3.5 seconds, 4m sampling lines ≤5.9 seconds.
- If sampling line and bite block are not connected via the O₂ luer, 100% of the O₂ flow will be delivered to the patient's nares. If sampling line and bite block are connected via the O₂ luer (as recommended), O₂ will simultaneously be delivered to the patient's nares and mouth, with 80% of the O₂ flow going to the patient's mouth and 20% dispersed through the O₂ holes below the patient's nares.
- After the bite block is removed from the patient at the end of the procedure, CO₂ measurement and O₂ delivery can continue during recovery with the sampling line alone. O₂ delivery should be ≤ 5L/min when being delivered only through the nasal cannula.

Reordering Information

Part Number	Description	Quantity
989803204451	Adult-Intermediate Bite Block CO ₂ Filter Line with O ₂ Connector, Short-term use: Procedural/Emergency Note: Length: 2m	25
989803204461	Adult-Intermediate Bite Block CO ₂ Filter Line with O ₂ Tubing, Short-term use: Procedural/Emergency Note: Length: 2m	25
989803204471	Adult-Intermediate Bite Block CO ₂ Filter Line with O ₂ Tubing, Short-term use: Procedural/Emergency Note: Length: 4m	25
989803204721	Adult-Intermediate Bite Block CO ₂ Filter Line with O ₂ Tubing, Short-term use: Procedural/Emergency Note: Length: 2m	25
989803204731	Adult-Intermediate Bite Block CO ₂ Filter Line with O ₂ Tubing, Short-term use: Procedural/Emergency Note: Length: 4m	25

Disposal

Discard any product that becomes damaged or deteriorated. Follow approved medical waste disposal methods specified by your facility or local regulations.

Incident Reporting

Any serious incident that has occurred in relation to this device should be reported to Philips and the competent authority of the European Economic Area (EEA) countries, including Switzerland and Turkey, in which the user and/or patient is established.

Environmental Specifications

The following table provides the temperature and humidity ranges that must be maintained for these products during storage and transport:

Specification	Storage	Transport
Temperature	-20°C to 55°C (-4°F to 131°F)	-20°C to 55°C (-4°F to 131°F)
Humidity	10% RH to 90% RH	10% RH to 90% RH

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Microstream™ is a trademark of
Copyright © 2021
Böblingen GmbH
Oridion Medical 1987 Ltd. Oridion Medical
1987 Ltd. is a Medtronic company.