

Disposable SpO₂ Sensors

Instructions for Use


Adult/Pediatric

M1131A

Infant

M1132A

Neonatal/Infant/Adult

M1133A

Adhesive-free

Neonatal/Infant/Adult

M1134A

PHILIPS

About This Edition

Contents

(BG) Инструкции за употреба.....	1	(LV) Lietošanas instrukcija.....	86
(CS) Návod k obsluze.....	6	(MK) Упатства за користење.....	91
(DA) Brugerhåndbog.....	11	(NL) Gebruiksaanwijzing.....	96
(DE) Gebrauchsanweisung.....	16	(NO) Bruksanvisning.....	101
(EL) Οδηγίες χρήσης.....	21	(PL) Instrukcja obsługi.....	106
(EN) Instructions for Use.....	26	(PT-BR) Instruções de uso.....	111
(ES) Instrucciones de uso.....	31	(RO) Instrucțiuni de utilizare.....	116
(ET) Kasutusjuhend.....	36	(RU) Инструкция по эксплуатации.....	121
(FI) Käyttöopas.....	41	(SK) Návod na použitie.....	126
(FR) Manuel d'utilisation.....	46	(SL) Navodila za uporabo.....	131
(HR) Upute za uporabu.....	51	(SR) Uputstvo za upotrebu.....	136
(HU) Használati útmutató.....	56	(SV) Bruksanvisning.....	141
(IT) Istruzioni d'uso.....	61	(TR) Kullanım Talimatları.....	146
(JA) ユーザーズガイド.....	66	(UK) Інструкції з використання.....	151
(KK) Пайдалану нұсқаулығы... 71		(ZH-CN) 使用说明.....	156
(KO) 사용 설명서.....	76	(ZH-TW) 使用說明.....	161
(LT) Naudojimo nurodymai.....	81		

R_x only

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Intended Use

The Philips disposable SpO₂ Sensors are for single patient use when continuous non-invasive arterial oxygen saturation and pulse rate monitoring are required. Sensors are for use by professionally trained clinicians such as nurses, physicians, and EMS personnel for the treatment of patients in and out of a hospital.

The M1131A, M1132A, and M1133A sensors are contraindicated for patients with allergic reactions to adhesive. See Application Table below.

Sensor Product	Neonate ¹ <3 Kg		Infant ² 3-10 Kg		Infant ³ 10-20 Kg		Pediatric ⁴ or Adult ⁴ >20 Kg	Adult ⁴ >40 Kg
	Foot	Hand	Finger	Toe	Big Toe	Thumb		
M1131A							✓	✓
M1132A			✓	✓				
M1133A	✓	✓			✓	✓		✓
M1134A⁵	✓	✓			✓	✓		✓

¹Preferred application site = foot; alternative = hand.

²Preferred application site = any finger; alternative = any toe.

³Preferred application site = big toe; alternative = thumb.

⁴Preferred application site = any finger except thumb.

⁵M1134A = adhesive-free sensor

WARNINGS

- The Philips disposable SpO₂ sensors are not sterile.
- To ensure cleanliness, be sure to unpack the disposable sensor just prior to use.
- Before using these sensors with any instrument, verify that the combination of instrument/sensor is specified in the instrument's Instructions for Use. Also read and understand all warnings listed in the instrument's Instructions for Use, as well as all warnings listed in this document, otherwise patient injury may result.
- Do not reuse. Single patient use only. Multi-patient use may result in patient infection or cross-infection. Also, do not use any sensor showing signs of damage or alteration for further patient monitoring; instead, dispose of it using proper disposal procedures (see below).
- Heavy dirt can be removed with a dry or damp cloth (M1131A Sensor only).
- Protect the connector from contact with any liquid.
- At elevated ambient temperatures, patient skin could be severely burned after prolonged sensor application at sites that are not well perfused. To prevent this condition, check patient application sites frequently. All listed sensors operate without risk of exceeding 41° C on the skin if the initial skin temperature does not exceed 35° C.
- Make sure you apply the sensor according to the instructions above at the preferred application site. Failure to do so can cause inaccurate measurements.
- To avoid venous pulsation, obstructed circulation, pressure marks, pressure necrosis, artifacts and inaccurate measurements, be sure you use a sensor of the correct size and that sensor is not too tight. If sensor is too loose, it might compromise the optical alignment or fall off. If sensor is too tight, because application site is too large or becomes too large due to edema, excessive pressure may be applied. This can result in venous congestion distal from the application site, leading to interstitial edema, hypoxia and tissue malnutrition.

- Where possible the application site for the sensor should be an extremity free of arterial catheters, blood pressure cuffs, or intravascular infusion lines.
- Avoid sites subject to excessive motion. Try to keep the patient still, or move the sensor to a site with less motion.
- Be aware that an unapplied or loose sensor may produce inaccurate readings.
- Using a M1134A sensor may compromise measurement performance on active patients (sensor may become dislodged more easily). On patients with significant motion use M1133A sensors instead.
- Make sure the application site for the sensor is not deeply pigmented, or deeply colored (for example, nail polish, artificial nails, dye or pigmented cream may cause inaccurate measurements). In any of these cases reposition the sensor or choose an alternative sensor for use at a different site.
- Dysfunctional hemoglobin or intravascular dyes can cause inaccurate measurements.
- Excessive light may cause inaccurate measurements, in such cases, cover the sensor with opaque material.
- Use correct sensor for your application.
- Only use with instruments validated by Philips and at recommended application area.
- Use 9-pin connector only with corresponding SpO₂ connection/adaptor cable on the instrument side!
- Do not use sensors during MRI scanning. Conducted current may cause burns and wrong readings.
- Avoid sensor from being put in patients mouth.
- Pulse oximetry measurements are statistically distributed. Two-thirds of all pulse oximetry measurements can be expected to fall within the stated accuracy (refer to *Specifications* later in this IFU, or the *SpO₂ Sensor Compatibility Sheet* for stated sensor accuracy).
- Inspect the sensor application site every 2 to 3 hours to ensure skin integrity, correct optical alignment, and circulation distal to the sensor site. Skin irritations may occur if the sensor is attached to one location for too long. Move the sensor application site every 4 hours, or more often if circulation or skin integrity is compromised. If the light source is not directly opposite the light detector, reapply the sensor or choose an alternative sensor for use at a different site.

Symbol Definition

Please refer to Pages i -vii.

Compatibility

These SpO₂ sensors are compatible with the Philips M8105A (IntelliVue MP5) monitor plus any other Philips device that lists them as accessories in the device's IFU. Incompatible components can result in degraded performance. Consult Monitor's IFU for further information.

Patient Application

Prior to patient application, visually inspect the sensor. Make sure cable is not damaged, connector pins are not bent, and sensor is not torn. Refer to instructions below for proper application. **Note:** Read the **Warnings** on pages 26 and 27 before applying your Philips Sensor.




Applying the M1131A Sensor

Step	Action
1	The M1131A Sensor is for use on a pediatric or adult fingertip (see page 26). If possible, choose a finger on an extremity free from arterial catheter, blood pressure cuff, or intravascular infusion line. Select the most appropriate finger for the sensor size.
2	Gently open the sensor and guide it over the fingertip, with the sensor cable on <i>top</i> of the finger and hand: • The fingertip should touch the end of the wide section of the sensor. • You may need to trim the patient's fingernail to position the sensor correctly. 
3	Pull the white backings from the sensor's adhesive tape, then: • Wrap the tape around both the sensor and cable (do not stretch the tape). • Plug the sensor cable into the monitoring instrument's 9-pin female adapter cable connector. 
4	Move the sensor application site every four hours, or more often if circulation or skin integrity is compromised.
5	If the sensor adhesive tape becomes compromised, use any medical grade tape to reapply the sensor securely around the fingertip.

Reapplying the M1131A Sensor

The sensor may be reapplied *to the same patient* by repeating the procedure above. If the adhesive tape supplied with the sensor is compromised, use any medical grade tape to secure the sensor in place.

Applying M1132A, M1133A, M1134A Sensors

Step	Action
1	Select an appropriate sensor for the application (see <i>Intended Use</i> on page 26 for patient applications). If possible, select a patient extremity free from arterial catheter, blood pressure cuff, or intravascular infusion line.
2	M1132A/M1133A: Peel the protective strip off the sensor, revealing two optical windows. The M1132A and M1133A sensors have adhesive around these windows to provide a more secure attachment to the patient. 
3	 Position the optical window that aligns with the sensor cable over the selected patient extremity, with the cable on <i>top</i> of the hand or <i>bottom</i> of the foot toward heel or toes, as applicable. 
4	Wrap sensor around extremity firmly, but not too tight, so that both optical windows oppose each other. Use the Velcro® at the end of the sensor to secure it in place. If necessary, secure the cable to the hand or foot with medical tape.
5	Plug the sensor cable into the monitoring instrument's 9-pin female adapter cable connector.
6	Move the sensor application site every four hours, or more often if circulation or skin integrity is compromised.

Reapplying the M1132A, M1133A, M1134A Sensor

These sensors may be reapplied *to same patient* up to 10 times using above procedure. Before reapplying the M1132A or M1133A sensor, gently wipe the optics side of the sensor with alcohol to rejuvenate the adhesive, allowing the alcohol to dry before reapplying. Assure the patient has no skin reaction from the dried alcohol.

Specifications

SpO₂ and Pulse Rate Accuracy

Refer to the Instructions for Use supplied with the monitoring instrument, or the *SpO₂ Sensor Compatibility Sheet* supplied with the SpO₂ sensors.

Sensor Accuracy

- **M1131A** (adult and pediatric) sensor accuracy when used in combination with the Philips M8105A (IntelliVue MP5) monitor is $\pm 3\%$ SpO₂ in the range from 70% to 100% SaO₂.

- **M1132A** (infant) and **M1133A, M1134A** (adult, infant and neonate) sensor accuracy is $\pm 2\%$ SpO₂ in the range from 70% to 100% SaO₂. When these sensors are used on neonate subjects as recommended, the literature suggests the accuracy range might increase by an additional +1% when performing an accuracy determination in a hospital environment.

For more information on SpO₂ and Pulse Rate accuracy specifications, reference your monitoring instrument's Instructions for Use.

Light Emitting Diode Wavelength Ranges

Wavelength ranges for the light emitting diodes used in these sensors are within 600nm-1000nm, with an optical output power of less than 15mW. Knowing the wavelength range can be useful to clinicians performing photodynamic therapy.

Measurement Validation

The SpO₂ accuracy has been validated in human studies against arterial blood sample reference measured with a CO-oximeter. In a controlled desaturation study, healthy adult volunteers with saturation levels between 70% and 100% SaO₂ were studied. The population characteristics for those studies were:

- Approximately 50% female and 50% male ranging in age from 18-45
- Skin tone: from light to black

Because pulse oximeter equipment measurements are statistically distributed, only approximately 2/3 of pulse oximeter equipment measurements can be expected to fall within the $\pm$ Arms value measured by a CO-oximeter. Functional testers, such as an SpO₂ simulator cannot be used to assess the accuracy of pulse oximeter sensors.

Environmental Requirements

Temperature and humidity ranges that must be maintained for all of these sensors are as listed below:

Temperature Requirements

- Operating Temperature Range: 0° C to 55° C (32° F to 131° F)
- Storage Temperature Range: 12° C to 35° C (53.6° F to 95° F)
- Transport Temperature Range: -40° C to 70° C (-40° F to 158° F)

Humidity Requirements

- Operating Humidity Range: $\leq 95\%$ RH at 40° C (104° F)
- Storage Humidity Range: 18% RH to 78% RH
- Transport Humidity Range: $\leq 90\%$ RH at 65° C (149° F)

Disposing the Sensor

Dispose the sensor properly following local regulations for hospital waste.

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© Copyright 2016
GmbH