
en

Identification of a substance that is not contained or present within the product or packaging.

Do not use if package is opened or damaged.

Directions for Use

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

This device is not made with natural rubber latex or DEHP.

Indications/Contraindications

The Nellcor™ Pediatric-Infant Foam Sensor Wrap, model FOAM P/I, is intended for use with an oxygen sensor, model OXI-P/I, or a Nellcor multisite oxygen sensor, when continuous noninvasive monitoring of arterial oxygen saturation and pulse rate is required for patients weighing between 3 kg and 40 kg.

The FOAM P/I wrap attaches the oxygen sensor to the sensor site and maintains appropriate positioning of the sensor's optical components.

The FOAM P/I is contraindicated for use on patients who exhibit allergic reactions to the pressure-sensitive adhesive.

Instructions for Use

Instructions for applying the FOAM P/I wraps to the sensor are specific to the particular Nellcor sensor used. For this reason, please refer to the Directions for Use packaged with the sensor for complete application instructions for that particular sensor.

Removal and Cleaning

Remove the wrap from the sensor by gently peeling it away.

Caution: Using excessive force when removing the FOAM P/I wrap may damage the sensor.

The FOAM P/I is intended for single patient use only. A new wrap is required for each application.

Warnings

- 1) Failure to apply the FOAM P/I and the sensor properly may cause incorrect measurements.
- 2) Reusable sensors must be moved to a new site at least every 4 hours. Because individual skin condition affects the ability of the skin to tolerate sensor placement, it may be necessary to change the sensor site more frequently with some patients. If skin integrity changes, move the sensor to another site.
- 3) Use only Nellcor wraps designed for use with Nellcor sensors. Do not use tape. Use of additional tape or various other types of wraps can cause skin damage.
- 4) When the sensor is wrapped too tightly or supplemental tape is applied, venous pulsations may lead to inaccurate saturation measurements.

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© 2011 Covidien.
COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company.
Covidien Inc, 15 Hampshire Street, Mansfield, MA 02048 USA.
EC REP Covidien Ireland Limited, IDA Business & Technology Park, Tullamore.