
Reusable SpO₂ T-Series Sensors

Instructions for Use

M1191T

Adult Finger Sensor

M1192T

Pediatric Finger Sensor

M1193T

Neonate Hand, Foot Sensor

PHILIPS

Contents

(BG) Инструкции за употреба... 1	(LV) Lietošanas instrukcija..... 120
(CS) Návod k obsluze..... 8	(MK) Упатства за користење .. 127
(DA) Brugerhåndbog 15	(NL) Gebruiksaanwijzing 134
(DE) Gebrauchsanweisung 22	(NO) Bruksanvisning..... 141
(EL) Οδηγίες χρήσης 29	(PL) Instrukcja obsługi 148
(EN) Instructions for Use 36	(PT-BR) Instruções de uso 155
(ES) Instrucciones de uso 43	(RO) Instrucțiuni de utilizare 162
(ET) Kasutusjuhend 50	(RU) Инструкция по эксплуатации 169
(FI) Käyttöopas..... 57	(SK) Návod na použitie 176
(FR) Manuel d'utilisation 64	(SL) Navodila za uporabo..... 183
(HR) Upute za uporabu..... 71	(SR) Uputstva za upotrebu..... 190
(HU) Használati útmutató..... 78	(SV) Bruksanvisning 197
(IT) Istruzioni d'uso 85	(TR) Kullanım Talimatları..... 204
(JA) ユーザーズガイド 92	(UK) Інструкції з використання..... 211
(KK) Пайдалану нұсқаулығы 99	(ZH-CN) 使用说明 218
(KO) 사용 설명서 106	(ZH-TW) 使用說明..... 225
(LT) Naudojimo nurodymai... 113	

R_x only

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Intended Use

Philips Reusable SpO₂ T-Series sensors provide continuous, non-invasive measurement of arterial oxygen saturation (provides pulse rate monitoring and plethysmograph wave) to compatible SpO₂ patient monitoring devices.

The sensors are suitable for use only by appropriately trained clinicians such as nurses, physicians, and EMS personnel for the treatment of patients in and out of a hospital, as well as during transport. All clinicians shall follow the current standard use model when applying a SpO₂ sensor to a patient; i.e. apply the M1191T Adult or M1192T Pediatric sensor to the patient's finger; apply the M1193T Neonate sensor to the patient's hand or foot (refer to Table on Page 2).

Warnings

- The SpO₂ T-Series Sensors are not sterile.
- These sensors are for use only with compatible monitoring instruments. Before using, reference your instrument's Instructions for Use, or the *Sensor & Instrument Compatibility Sheet* shipped with each sensor, to verify instrument compatibility. Be sure to use sensors only with validated instruments, otherwise patient injury may result.
- Incorrectly applied sensors may cause inaccurate readings that could lead to misdiagnosis. Always verify the sensor is properly applied.
- Pulse oximetry measurements are statistically distributed. Two-thirds of all pulse oximetry measurements can be expected to fall within the stated accuracy (refer to *Specifications* later in this IFU, or the *Sensor & Instrument Compatibility Sheet* for stated sensor accuracy).
- Connect sensor only to the SpO₂ connection or SpO₂ adapter cable of oximeter.
- Do not reuse the sensor on a *different* patient until it has been disinfected because contamination may occur. The sensor can be reused on the same patient for the patient's entire stay.
- At elevated ambient temperatures, patient skin could be severely burned after prolonged sensor application at sites that are not well perfused. To prevent this condition, be sure to check patient application sites frequently. All listed sensors operate without risk of exceeding 41°C on the skin if the initial skin temperature does not exceed 35°C.
- Be sure to apply the sensor correctly at the preferred application site, following instructions provided later in this document (refer to *Applying the Sensor*). Failure to do so can cause inaccurate measurements.
- To avoid venous pulsation, obstructed circulation, pressure marks, pressure necrosis, artifacts and inaccurate measurements, make sure you are using a sensor of the correct size and that the sensor is not too tight. If the sensor is too tight, because the application site is too large or becomes too large due to edema, excessive pressure may be applied. This can result in venous congestion distal from the application site, leading to interstitial edema and tissue ischemia.
- If the sensor is applied too loose, it may fall off or compromise proper alignment of sensor optics and result in inaccurate readings.

- Where possible, the application site for the sensor should be an extremity free of arterial catheters, blood pressure cuffs, or intravascular infusion lines.
- Avoid application sites subject to excessive motion. Try to keep the patient still, or move the sensor to a site with less motion.
- Dysfunctional hemoglobin or intravascular dyes can cause inaccurate measurements.
- Make sure the application site for the sensor is not deeply pigmented, or deeply colored (for example: nail polish, artificial nails, dye or pigmented cream may cause inaccurate measurements). In any of these cases reposition the sensor or choose an alternative sensor for use at a different site.
- Cover the sensor with opaque material under conditions of strong or excessive light (infrared lamps, OR lamps, photo therapy). Failure to do so can result in inaccurate measurements.
- Protect the connector from contact with any liquid.
- Inspect sensor application site every 2 to 3 hours to ensure skin integrity, correct optical alignment, and circulation distal to the sensor site. Skin irritation or ulceration may occur if sensor is attached to one location for too long. Move the sensor application site every 4 hours, or more often if circulation or skin integrity is compromised. If the light source is not directly opposite the light detector, reapply the sensor, or choose an alternative sensor for use at a different site.
- Do not use a sensor during MRI scanning. This may cause burns or inaccurate measurements.

T-Series Sensor/Monitor Compatibility

T-Series sensors are compatible with the Philips M8105A (IntelliVue MP5) monitor plus many *validated* Philips/Agilent/HP and OEM (Original Equipment Manufacturer) monitors. Reference the *Sensor & Instrument Compatibility Sheet* included with your T-Series Sensor for a comprehensive list of compatible monitors.

Reusable T-Series Sensors	Model	Application	Cable	Connector
	M1191T Sensor	Adult Any finger except thumb	45 cm	 9-pin D-sub Connector
	M1192T Sensor	Pediatric Any finger except thumb	45 cm	
	M1193T Sensor	Neonatal Foot or Hand	90 cm	

Accessories Supplied with Sensors

- **Sensor & Instrument Compatibility Sheet:** Lists all instruments that are compatible with the reusable T-Series Sensors.
- **Instructions for Use:** Includes T-Series Sensor product descriptions, patient applications, Warnings against improper use, cleaning and disinfection procedures, and sensor specifications.

Accessories Available for Separate Purchase

- **M1900B Adapter Cable:** 9-pin (D-Sub female) to 12-pin (round male) *adapter* cable. Connects any T-Series Sensor 9-pin connector to any compatible Philips/Agilent/HP instrument that has a 12-pin SpO₂ connection. Extends sensor cable length by 3m.
- **M1943A Adapter Cable:** 9-pin (D-sub female) to 8-pin (round male) *adapter* cable. Connects any T-Series Sensor 9-pin connector to any compatible Philips/Agilent/HP instrument that has an 8-pin SpO₂ connection. Extends sensor cable length by 1.1m.
- **M1943AL Adapter Cable:** A longer (3m) version of the M1943A *adapter* cable listed above.

Definition of Symbols

Refer to Pages i - vii.

Before Applying the Sensor

Before applying any T-Series Sensor to a patient, be sure to read and understand all Warnings listed in your SpO₂ monitoring instrument IFU, as well as the Warnings in this sensor IFU. Always apply an appropriate size sensor at the proper patient application area (finger for adult and pediatric patients; hand or foot for neonatal patients).

Inspect Sensor for Damage

- Inspect the sensor outside and inside. To check the inside, gently open the sensor cavity and check for splits on or next to the transparent silicone that covers the optical elements. Be sure there are no blisters on the silicone, and no liquid leaking from sensor optics.
- Any sensor showing signs of damage or alteration must not be used for further patient monitoring; instead, dispose of it using proper disposal procedures (see below).

Change Application Site Periodically

Move the sensor application site every 4 hours, or more often if circulation or skin integrity is compromised.

Applying the Finger Sensor

M1191T Adult Finger Sensor

For patients weighing more than 50kg.

M1192T Pediatric Finger Sensor

For patients weighing between 15kg and 50kg.

*adult application shown
with supplied wristband*



Step	Applying the M1191T (Adult) or M1192T (Pediatric) Sensor
1	Select appropriate sensor for patient size (defined above).
2	Position sensor over patient's finger, being sure the sensor cable is positioned ON TOP of the finger/hand as shown above.
3	Patient fingertip should touch but not protrude from end of sensor. If necessary, trim the fingernail to correctly position the sensor.
4	If applying a M1191T Adult finger sensor, secure the cable to BACK OF HAND with supplied M1627A wristband (wristband is supplied only with M1191T Adult finger sensors.)
5	Plug sensor into instrument (or into adapter cable if needed).
6	Inspect and change the sensor application site periodically. You may reapply the sensor to the <i>same</i> patient by repeating this procedure. Clean and disinfect the sensor before applying it to a <i>different</i> patient.

Applying the Hand/Foot Sensor

M1193T Neonatal Hand/Foot Sensor

For patients weighing 1kg-4kg.

Step	Applying the M1193T (Neonate) Sensor
1	Position the sensor over the hand or foot, with Sensor Optics aligned opposite each other. Hold the sensor in position while stretching the Strap slightly, not more than 2.5cm (1.0"). 
2	Press the stretched Strap against the Slot, thread it into the Latch, and pull for a secure fit. If long enough, thread it through the Second Latch as well. 
3	Plug the sensor cable/connector into the monitoring instrument (or the instrument's adapter cable if used).
4	Inspect and change the sensor application site periodically. You may reapply the sensor to the <i>same</i> patient by repeating this procedure. Clean and disinfect the sensor before applying it to a <i>different</i> patient.

Specifications

Sensor Accuracy

- M1191T (Adult) and M1192T (Pediatric) sensor accuracy is $\pm 3\%$ SpO₂ in the range from 70% to 100% SpO₂.
- M1193T (neonate) sensor accuracy is $\pm 3\%$ SpO₂ in the range from 70% to 100% SpO₂. The M1193T sensor has been validated on healthy adult volunteers in a controlled laboratory environment. When used on neonate subjects as recommended, the literature suggests the accuracy range might increase by an additional +1% when performing an accuracy determination in a hospital environment.

For more information on SpO₂ and Pulse Rate accuracy specifications, reference your monitoring instrument's Instructions for Use, or the *Sensor & Instrument Compatibility Sheet* included with your sensor.

Light Emitting Diode Wavelength Ranges

Wavelength ranges for the light emitting diodes used in these sensors are within 600nm-1000nm, with an optical output power of less than 15mW. Knowing the wavelength range can be useful to clinicians performing photodynamic therapy.

Measurement Validation

The SpO₂ accuracy has been validated in human studies against arterial blood sample reference measured with a CO-oximeter. In a controlled desaturation study, healthy adult volunteers with saturation levels between 70% and 100% SaO₂ were studied. The population characteristics for those studies were:

- Approximately 50% female and 50% male ranging in age from 18-45
- Skin tone: from light to black

Because pulse oximeter equipment measurements are statistically distributed, only approximately 2/3 of pulse oximeter equipment measurements are expected to fall within the $\pm$ Arms value measured by a CO-oximeter. Functional testers, such as an SpO₂ simulator cannot be used to assess the accuracy of pulse oximeter sensors.

Sensor Cleaning and Low-Level Disinfection

The reusable sensors should be cleaned and disinfected, but never sterilized. Follow the procedure outlined below.

Warnings

- Use only validated cleaning agents and disinfectants listed below; do not use any others. Failure to do so may damage the sensor or its connecting wires and shorten product lifetime or cause safety hazards.
- Select disinfectants carefully as some have very similar names but completely different compositions.
- Do not immerse sensor connector in any cleaning solutions, disinfectant, or other liquid (only sensor and cable housing may be immersed, not connector).
- Do not soak sensors in disinfectants longer than specified by the disinfectant manufacturer.
- Do not sterilize the sensors.

Validated Cleaner/Disinfectants

CLEANER, DISINFECTANT	
Brands listed with ® or TM/MC are manufacturer's trademarks	
Mild soap	Caltech-Dispatch 5200 Wipes
Isopropyl Alcohol	Viraguard®
Hydrogen Peroxide	TechSpray®
Clorox Healthcare® Bleach Germicidal Wipes	Accel TB RTU
Sodium Hypochlorite (household chlorine bleach 5% diluted 1:10 with water)	SANI-CLOTHS® BLEACH (Light Orange Top)
Metrex CaviWipes	SANI-CLOTH® HB (Blue Top)
Oxivir® Tb Wipes	Bacillol® 25
SANI-CLOTH® PLUS (Dark Orange Top)	Mikrozyd® sensitive wipes premium
SUPER SANI-CLOTH® (Purple Top)	Anios Surfa'safe
Bacillol AF	Surfanios Citro
Meliseptol	

Step	Cleaning and Low-Level Disinfection
1	Clean sensor following instructions supplied with cleaning agent.
2	Disinfect sensor following instructions supplied with disinfectant.
3	Rinse sensor in water, wiping it dry with a clean cloth, leaving to dry completely. If you notice any deterioration or damage to sensor or cable, dispose the sensor immediately.

Sensor Disposal

Dispose Deteriorated Sensors

Any sensor showing signs of physical or electrical deterioration or failure should be disinfected, decontaminated, and disposed of according to local laws regarding the disposal of hospital waste.

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© 2016
Böblingen GmbH