Neonatal and Infant Single-Patient NBP Cuffs

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



EN - Instructions for Use

Associated Documents

For more information, refer to the *Interconnect Hose for NIBP Cuffs Instructions for Use*.

Intended Use

The Philips Neonatal and Infant Single-Patient NBP Cuffs are to be used with compatible neonatal monitors (see *Compatibility*) intended for use by, or under the supervision of, a licensed physician or other healthcare provider for the noninvasive measurement of neonatal and infant human blood pressure.

Indications

Philips Neonatal and Infant Single-Patient NBP Cuffs are indicated for the measurement of non-invasive blood pressure and pulse rate for monitoring purposes. These devices are limited by the indications for use of the connected monitor in healthcare facilities. These devices are intended to interact with patient intact skin only.

Contraindications

There are no known contraindications.

Compatibility

Neonatal NBP Cuffs can be used with any monitors that lists the cuffs as accessories in that product's instructions for use.

The Neonatal NBP Cuffs can be used with the following Neonatal Blood Pressure Air Hoses:

- M1596C (1.5m Hose)
- M1597C (3.0m Hose)

Warnings

- Single-Patient cuffs are not intended for use on more than one patient.
- Use only monitor hoses identified for use with neonatal cuffs. Use of other hoses can prevent monitor from obtaining readings.
- Be sure to select the correct patient type setting on the monitor. Do not apply the higher adult inflation, overpressure limits, and measurement duration to neonatal patients. This may cause a patient injury or inaccurate measurements.
- Be sure cuff is correctly matched to patient type and size (limb circumference ranges are printed on each cuff). Cuffs that are too small may produce false high blood pressure readings; cuffs that are too large may produce false low-pressure readings.
- Kinking of the tube can result in inaccurate readings: avoid placement of cuff and tubing which could lead to kinking. Kinking may cause patient injury.
- To prevent damage to cuff (e.g. leaks), shut off auto function or remove cuff from monitor when not on patient.

- Use cuff only under direct supervision by a trained healthcare professional when attached to automated monitors without fail-to-read alarms.
- Periodically inspect cuff, connectors, and hose for signs of wear or deterioration. Replace damaged or deteriorated components, as required.
- Replace cuff if hook and loop fastener fails to hold during inflation.
- Follow your institution's guidelines for long-term patient care.
- Do not sterilize cuffs. Do not clean or disinfect. Do not use bleach.

Caution

Federal (US) law restricts this device to sale by or on the order of a medical practitioner.

Applying the Single Patient NBP Cuff

Note 1: Whenever visual inspection reveals a cuff is no longer suitable for continued use, does not perform as intended, or severe soiling occurs, dispose of the cuff. Refer to *Disposal* section.

Note 2: Refer also to the figure on page i. Numbers in bold in the following procedure correspond to the figure panel numbers.

- 1 Select the proper cuff size for the patient's limb circumference. Refer to the *Select Appropriate Cuff Size* section that follows.
- 2 Apply cuff to limb (**1**) such that the side displaying the **Caution** symbol  is toward the patient. The arrow labeled **Arteria** needs to be over the limb artery.
- 3 Wrap the cuff snugly around the patient's limb (**2, 3** and **4**). Ensure that the index edge falls within the range indicated by the range arrow (**5**).
- 4 Secure the connector(s) on the cuff to the air hose connector to prevent leaking. Refer to figures on page ii.

Select Appropriate Cuff Size

Philips Neonatal and Infant Single-Patient NBP Cuffs are available to suit a variety of patient applications and are color-coded for easy size selection.

	Neonatal # 1	Neonatal # 2	Neonatal # 3	Neonatal # 4	Infant # 5
Single-Patient Cuffs	M1866B/S	M1868B/S	M1870B/S	M1872B/S	M1873B/S
Circumference (cm)	3.1–5.7	4.3–8.0	5.8–10.9	7.1–13.1	10.0–15.0
Cuff Text Color	Orange	Light Blue	Green	Navy	Burgundy

Reordering Information

Part	Description	Part	Description
M1866B	Neonatal Single-Patient Cuff Size #1	M1870S	Neonatal Soft Single-Patient Cuff Size #3
M1866S	Neonatal Soft Single-Patient Cuff Size #1	M1872B	Neonatal Single-Patient Cuff Size #4
M1868B	Neonatal Single-Patient Cuff Size #2	M1872S	Neonatal Soft Single-Patient Cuff Size #4
M1868S	Neonatal Soft Single-Patient Cuff Size #2	M1873B	Infant Single-Patient Cuff Size #5
M1870B	Neonatal Single-Patient Cuff Size #3	M1873S	Infant Soft Single-Patient Cuff Size #5

Cleaning and Care

Neonatal NBP Cuffs are for single-patient use and are not intended to be cleaned or disinfected. If cuff is soiled, discard. Refer to *Disposal* section that follows.

Disposal

When finished using this single use cuff, follow approved disposal methods for medical waste as specified by your facility (hospital, clinic) or local regulations

Incident Reporting

Any serious incident that has occurred in relation to this device should be reported to Philips and the competent authority of the European Economic Area (EEA) countries, including Switzerland and Turkey in which the user and/or patient is established.

Specifications

Specification	Operating	Storage
Temperature	0°C to +50°C (32° to 122°F)	-20° to +60°C (-4° to +140°F)
Humidity	0 to 95% RH	0% RH to 90% RH@ 60°C (140°F)

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Copyright © 2020
Nederland B.V.