② Single Patient Use Only

STERILE EO

Sterile

Latex Free

For use with OxiMAX™ technology only

0123

CE 2005-12

Directions for Use

Indications/Contraindications

The Nellcor® adhesive reflectance oximetry sensor, model MAX-FAST, is indicated for use when continuous noninvasive arterial oxygen saturation and pulse rate monitoring are required for patients weighing greater than 40 kg.

The MAX-FAST oximetry sensor is supplied sterile for single patient use only. This sensor does not contain latex. Affixed to the sensor are three tabbed adhesive layers and a white paper backing. The white paper backing is removed prior to initial application on a patient. When the sensor is moved to a new site on the same patient, a fresh adhesive layer may be exposed by peeling the current adhesive layer away from the sensor.

This oximetry sensor is for use only with Nellcor OxiMAX™ instruments and instruments containing Nellcor OxiMAX oximetry technology. Consult individual manufacturers for features and compatibility of particular instruments and this sensor model.

Each manufacturer of Nellcor-compatible instruments is responsible for determining whether, and under what conditions, its instruments are compatible for safe and effective use with each Nellcor sensor model. This may include different specifications and/or warnings, cautions, or contraindications. Refer to the instrument operator's manual or consult the manufacturer for complete instructions for use of this sensor with their Nellcor-compatible instrument. For additional warnings, cautions or contraindications when using this sensor with Nellcor-compatible instruments, refer to the instrument operator's manual or contact the manufacturer of the instrument.

The MAX-FAST sensor is contraindicated for patients who exhibit allergic reactions to the adhesive pad; for patients who perspire profusely; or under conditions where the patient is in the Trendelenburg's position (head lower than the heart).

Instructions for Use

The MAX-FAST sensor may be reused on the same patient. The sensor may be used on the same site for a maximum of 12 hours, provided the site is inspected routinely to ensure skin integrity, correct positioning, and proper adhesion of the sensor. Because individual skin condition affects the ability of the skin to tolerate sensor placement, it may be necessary to change the sensor site more frequently with some patients.

Note: If the sensor does not track the pulse reliably, it may be incorrectly positioned—or the sensor site may be excessively wrinkled, or too deeply pigmented or otherwise deeply colored (for example, as a result of externally applied coloring such as dye or pigmented cream) to permit appropriate light transmission. If any of these situations occurs, reposition the sensor or choose an alternate Nellcor sensor for use on a different site.

Warnings

- 1) As with all medical equipment, carefully route cables to reduce the possibility of patient entanglement or strangulation.
- 2) Do not use the MAX-FAST sensor or other oximetry sensors during MRI scanning. Conducted current may cause burns. Also, the MAX-FAST sensor may affect the MRI image, and the MRI unit may affect the accuracy of oximetry measurements.
- 3) Oxygen saturation readings may be falsely low (i.e., a mixture of arterial and venous values) in the presence of tricuspid regurgitation due to retrograde venous pulsations.
- 4) Do not use a damaged sensor or pulse oximetry cable. Do not use a sensor with exposed optical components.

Precautions

- 1) Failure to apply the MAX-FAST sensor properly may cause incorrect measurements. Always follow the *Directions for Use* to assure proper placement.
- 2) While the MAX-FAST sensor is designed to reduce the effects of ambient light, excessive light might cause inaccurate measurements. In such cases, cover the sensor with opaque material.
- 3) Excessive patient motion may compromise performance. In such cases, try to keep the patient still, or choose an alternate Nellcor sensor for use on a different site.
- 4) Do not use tape with the MAX-FAST sensor. Use of additional tape or other types of adhesives may cause skin damage.
- 5) Venous pulsations, resulting from wrapping the foam headband too tightly, can lead to inaccurate saturation measurements.
- 6) Do not alter or modify the MAX-FAST sensor. Alterations or modifications may affect performance or accuracy.
- 7) For patients in a prone position, venous pooling may cause inaccurate measurement. Do not immerse in water or cleaning.

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- remove skin oils. See illustration for the recommended site. **A**
- Locate the two optical components of the sensor. **B**
- Remove the white paper backing to expose the first of three adhesive pads. **C** The MAX-FAST sensor is now ready to apply to the patient.

Note: There are three adhesive pads affixed to the sensor, each with a pull-tab for removal. When repositioning the sensor on the same patient, first expose a new adhesive pad by grasping the tab and peeling off the old adhesive pad. The sensor is ready to be reapplied to the patient.

- Press the MAX-FAST sensor firmly in place for 10 seconds, ensuring that the entire surface area of the adhesive pad makes contact with the skin.

Note: The MAX-FAST sensor must be applied to a flat, hairless surface of the forehead.

- If desired, the sensor cable can be secured to the patient's clothing or other material by using the clip located on the cable. To open, pinch the sides of the clip; release to close. **D**
- If the optional foam headband is used, carefully wrap the headband around the patient's head, securing the headband in place with the tapeless fastener. **E** Never wrap the headband tightly. Ensure that the sensor is completely covered. Ensure that the headband applies equal pressure to the entire sensor.

Note: The headband is sold separately. Contact Nellcor's Customer Service Department or your local Nellcor representative for ordering information.

- Plug the MAX-FAST oximetry sensor into the oximeter, turn on the oximeter, and verify proper operation as described in the oximeter operator's manual.

- In the event of damage to the sterile packaging, do not resterilize. Follow local governing ordinances regarding the disposal of sensors.

If you have questions regarding any of this information, contact Nellcor's Technical Services Department, or your local Nellcor representative.

Accuracy Specifications

For the accuracy specification range when used with Nellcor monitors, refer to information provided with the monitor, or, in the U.S., contact Nellcor's Technical Services Department. Outside the U.S., contact your local Nellcor representative.

For the accuracy specification range of this sensor when used with Nellcor-compatible instruments, refer to the instrument operator's manual or contact the instrument manufacturer.

Additional Copies of Instructions

Additional copies of these Instructions are available at no charge by calling Nellcor or its authorized distributors. Also, permission is hereby granted under Nellcor Puritan Bennett copyrights to Purchasers of products obtained from Nellcor or its authorized distributors to make additional copies of these instructions for use by such purchasers.

Rx ONLY


This sensor is sold under the following U.S. patents and foreign equivalents, for single patient use only: 4,621,643; 4,700,708; and 4,830,014. Any other use of this sensor is not authorized by Nellcor Puritan Bennett under the referenced patents.



© 2001 Nellcor Puritan Bennett Inc. All rights reserved.

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© 2001 Nellcor Puritan Bennett Inc. Tous droits réservés.