Medtronic

Microstream™ Luer

Adult-Pediatric Airway Adapter
Sampling Line

en	Instructions For Use
fr	Mode d'emploi
de	Gebrauchsanweisung
nl	Gebruiksaanwijzing
it	Istruzioni per l'uso
es	Instrucciones de uso
pt	Instruções de uso
ru	Инструкция по эксплуатации
da	Brugsanvisning
sv	Bruksanvisning
el	Οδηγίες Χρήσης
no	Bruksanvisning
hu	Használati útmutató
pl	Instrukcja użycia
cs	Návod k použití
fi	Käyttöohjeet
tr	Kullanma Talimatları
sk	Pokyny na používanie
sl	Navodila za uporabo
ro	Instrucţiuni de utilizare
bg	Инструкции за употреба
et	Kasutusjuhised
lt	Naudojimo instrukcijos
lv	Lietošanas norādījumi
hr	Upute za uporabu
sr	Uputstva za upotrebu


REF 010989

REF ML


Table of Contents

English	7
French	8
German	9
Dutch	10
Italian	11
Spanish	12
Portuguese	13
Russian.....	14
Danish.....	15
Swedish	16
Greek.....	17
Norwegian	18
Hungarian	19
Polish.....	20
Czech.....	21
Finnish	22
Turkish	23
Slovak	24
Slovenian.....	25
Romanian.....	26
Bulgarian	27
Estonian.....	28
Lithuanian.....	29
Latvian.....	30
Croatian.....	31
Serbian.....	32

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
Airway Adapter and Sampling Line

Instructions for Use

Single patient use, not to be re-processed

These products are intended to conduct a sample of the patient's breathing from a ventilator or anesthesia machine to a gas measurement device for measuring the percentage of CO₂ in the patient's exhaled breath.

Microstream Luer products covered in this DFU are intended to be used with capnography monitors with the following requirements:

- Have a standard luer connector that meets the requirements of ISO 594-2
- Comply with capnography standard (80601-2-55)
- Operate at a sampling flow rate of 90-180 ml/min
- Can withstand a sampling line pressure drop of up to 70mBar @ 180ml/min.

Airway Adapters for Adult to Pediatric patients:

- Added dead space <6.6 cc.
- To be used with an endotracheal tube > 4.5mm bore.
- Suitable for use with 22mm and 15mm patient tubing.

WARNINGS

- **Do not attempt to clean, disinfect, sterilize or flush any part of the sampling line. Re-use of single-use accessories could pose a cross-contamination risk to the patient or damage the functioning of the monitor.**
- **When used with a suction catheter closed suction system, do not place the airway adapter between the suction catheter and endotracheal tube. This is to ensure that the airway adapter does not interfere with the functioning of the suction catheter, and that the suction process does not damage the airway adapter. Damage of the airway adapter may harm the patient.**
- **When connecting a sampling line to the monitor, insert the sampling line connector clockwise into the monitor CO₂ port and turn the connector until it can no longer be turned, to ensure that it is connected securely to the monitor. This will assure that there is no leak of gases during measurement at the connection point and that measurement accuracy is not compromised. Ensure CO₂ values appear on the screen.**
- **Ensure no kinks are present in CO₂ and O₂ tubing, as kinked tubing may cause inaccurate CO₂ sampling or affect O₂ delivery.**
- **Carefully route the sampling line to reduce the possibility of patient entanglement or strangulation.**

CAUTIONS

- Dispose of sampling lines according to standard operating procedures or local regulations for the disposal of contaminated medical waste.

Notes

- During nebulization, lavage, or suction, follow monitor instructions for use (IFU) to suspend pump activity, in order to avoid moisture buildup and sampling line occlusion. In addition, especially if unable to suspend the CO₂ pump, disconnect the sampling line connector from the CO₂ port on the monitor.
- During setup, ensure that airway adapter can be easily attached and detached from the tubing before proceeding.
- Response time of the 2m sampling line is ≤3.5 seconds.
- Replace the sampling line according to hospital protocol or when a blockage is indicated by the device. Excessive patient secretions or a build-up of liquids in the airway tubing may occlude the sampling line, requiring more frequent replacement.
- Any serious incident related to device use that may occur should be reported immediately to the manufacturer, the local competent authority, and any other regulators as required.

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Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.
Other brands are trademarks of a Medtronic company.
© 2019 Medtronic. All rights reserved.
Oridion Medical 1987 Ltd., 7 Hamarpe St.,