Microstream™ Luer


Adult-Pediatric Intubated CO₂ Filter Line,
Short-term use: Procedural/Emergency

en	Instructions for Use
fr	Mode d'emploi
de	Gebrauchsanweisung
nl	Gebruiksaanwijzing
it	Istruzioni per l'uso
es	Instrucciones de uso
pt	Instruções de uso
ru	Инструкция по эксплуатации
da	Brugsanvisning
sv	Bruksanvisning
el	Οδηγίες χρήσης
no	Bruksanvisning
hu	Használati útmutató
pl	Instrukcja użycia
cs	Návod k použití
fi	Käyttöohjeet
tr	Kullanma Talimatları
sk	Pokyny na používanie
sl	Navodila za uporabo
ro	Instrucțiuni de utilizare
bg	Инструкции за употреба
et	Kasutusjuhised
lt	Naudojimo instrukcijos
lv	Lietošanas norādījumi
hr	Upute za uporabu
sr	Uputstva za upotrebu


Table of Contents

English	10
French	11
German	12
Dutch	13
Italian	14
Spanish	15
Portuguese	16
Russian.....	17
Danish.....	18
Swedish	19
Greek.....	20
Norwegian	21
Hungarian	22
Polish.....	23
Czech.....	24
Finnish	25
Turkish	26
Slovak	27
Slovenian.....	28
Romanian.....	29
Bulgarian	30
Estonian.....	31
Lithuanian.....	32
Latvian.....	33
Croatian.....	34
Serbian.....	35

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
Microstream™ Luer

Adult-Pediatric Intubated CO₂ Filter Line

Indications for Use:

Intended to conduct a sample of the patient's breathing from a ventilator or anesthesia machine to a gas measurement device for measuring the percentage of CO₂ in the patient's exhaled breath. The set is intended for single patient use only. Intended population: Intubated Adult-Pediatric patients.

Microstream Luer products covered in this DFU are intended to be used with capnography monitors with the following requirements:

- Have a standard luer connector that meets the requirements of ISO 594-2
- Comply with capnography standard (80601-2-55)
- Operate at a sampling flow rate of 90-180 ml/min
- Can withstand a sampling line pressure drop of up to 70mBar @ 180ml/min.

Airway adapters for adult to pediatric patients:

- Added dead space <6.6 cc.
- To be used with an endotracheal tube > 4.5mm bore.
- Suitable for use with 22mm and 15mm patient tubing.

WARNINGS

- **Loose or damaged connections may compromise ventilation or cause an inaccurate measurement of respiratory gases. Securely connect all components, turning the sampling line connector clockwise into the monitor CO₂ port until it can no longer be turned. Check connections for leaks according to standard clinical procedures, and ensure CO₂ values appear.**
- **Do not attempt to clean, disinfect, sterilize or flush any part of the sampling line. Re-use of single-use accessories could pose a cross-contamination risk to the patient or damage the functioning of the monitor.**
- **When used with a suction catheter, do not place the airway adapter between the suction catheter and endotracheal tube. This is to ensure that the airway adapter does not interfere with the functioning of the suction catheter, and that the suction process does not damage the airway adapter. Damage of the airway adapter may harm the patient.**
- **Ensure no kinks are present in CO₂ tubing, as kinked tubing may cause inaccurate CO₂ sampling.**
- **Carefully route the sampling line to reduce the possibility of patient entanglement or strangulation.**

CAUTIONS

- Dispose of sampling lines according to standard operating procedures or local regulations for the disposal of contaminated medical waste.

Notes

- During setup, ensure that airway adapter can be easily attached and detached from the tubing before proceeding.
- During nebulization, lavage, or suction, follow monitor instructions for use (IFU) to suspend pump activity, in order to avoid moisture buildup and sampling line occlusion. In addition, especially if unable to suspend the CO₂ pump, disconnect the sampling line connector from the CO₂ port on the monitor.
- The response time of the 2m sampling line is ≤3.5 seconds, 4m sampling lines ≤5.9 seconds.
- Replace the sampling line according to hospital protocol or when a blockage is indicated by the device. Excessive patient secretions or a build-up of liquids in the airway tubing may occlude the sampling line, requiring more frequent replacement.
- Any serious incident related to device use that may occur should be reported immediately to the manufacturer, the local competent authority, and any other regulators as required.

TM* Trademark of its respective owner.

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Together are trademarks of Medtronic.
Other brands are trademarks of a
© 2019 Medtronic. All rights reserved.
Oridion Medical 1987 Ltd., 7 Hamarpe St.,