Medtronic

REF OxySoftNHC


Nellcor™ OxySoft™

Neonatal-Adult SpO₂ Sensor

en	Home Use Guide
fr	Mode d'emploi à domicile
de	Bedienungsanleitung für die Heimanwendung
nl	Handleiding voor thuisgebruik
it	Guida per l'uso domiciliare
es-L	Guía de uso doméstico
sv	Manual för hemmabruk
da	Vejledning til hjemmebrug
no	Veiledning for hjemmebruk
fi	Kotikäyttöopas
pt-BR	Guia de uso doméstico
zh-CN	家庭使用指南
pl	Podręcznik stosowania w warunkach domowych
tr	Evde Kullanım Kılavuzu


1


2


3


1

4

5

Glossary

Ambient Light – Bright light sources in the area of the patient sensor, such as fluorescent lights, infrared heating lamps, and direct sunlight.

Bland-Altman – Graphs used to visualize the differences in measurements between two different instruments or two different measurement techniques to assess accuracy.

Clinician – The trained healthcare professional that assists with monitoring the patient and using the sensor in your home like a doctor or nurse.

Detector – Receives light from emitter.

Emitter – Light source from the sensor.

Flammable – Capable of catching fire and burning quickly. Some examples of flammable materials are gasoline, propane, and natural gas.

Operator – A person who operates equipment.

Oxygen Saturation (SpO₂) – A measurement of the percentage of oxygen circulating in the patient's blood.

Patient – A person receiving or registered to receive medical treatment or monitoring.

Pulse Oximetry – A monitor that is outside the body to determine the pulse rate and level of oxygen saturation.

Pulse Rate (PR) – A measurement of the number of times the heart beats per minute by monitoring blood pulsing through the body.

Retinopathy of Prematurity (ROP) – A potentially blinding eye disorder that primarily affects premature infants.

Sensor – A device used to collect and send patient information to a monitoring system such as pulse rate and amount of oxygen in the blood.

Sensor Site – The place on the patient's body where the sensor is applied, such as the finger.

Sterile – Free from bacteria or other living microorganisms; totally clean.

Note: This Instruction for Use is intended for the sensor operator. The sensor may be prescribed by the clinician (patient's doctor) for use at home by the patient or patient's caregiver. The home patient or patient's caregiver is the intended operator. For professional use, the Instruction for Use is available at <https://manuals.medtronic.com/manuals/main/region>

General Description of the Device

The Nellcor™ OxySoft™ Neonatal-Adult SpO₂ Sensor, model OxySoftN, is intended for use in non-invasive spot-check or continuous monitoring for up to 3 days of blood oxygen saturation (SpO₂) and pulse rate for neonates and infants weighing less than 3 kg or adults and pediatrics weighing more than 40 kg.

Intended Purpose

Pulse oximeters provide measurements of arterial blood oxygen saturation (SpO₂) and pulse rate.

Indications for Use (Who Can Use This Device)

The Nellcor™ OxySoftN single patient use sensor is indicated for use with neonatal, infant, pediatric and adult patients during both no motion and motion conditions and for patients who are either well or poorly perfused, in hospitals, hospital-type facilities, intra-hospital transport, mobile emergency medical applications including both ground and air transport and home environments.

Transport environments include intra-hospital transport and both ground and air emergency transport: road ambulances and fixed-wing aircraft and helicopters. The Nellcor™ OxySoft™ SpO₂ Sensor works with monitoring systems that use OxiMax™ and Nellcor™ compatible pulse oximetry systems to facilitate spot checking or continuous monitoring of a patient's arterial hemoglobin (SpO₂) and pulse rate via topical application of the sensor over a finger, foot, or hand. This device is for prescription use only.

Contraindications (Who Should Not Use This Device)

The OxySoftNHC sensor is contraindicated for use on patients who exhibit allergic reactions to synthetic fabrics or silicone.

Safety Information

Read this entire guide before using the sensor. This guide provides important information for avoiding injury and for proper use of the sensor.

WARNINGS (What You Must Do to Avoid Serious Harm)

- Only use the sensor with the compatible systems provided within this Home Use Guide. Using with unapproved or third-party systems can result in inaccurate SpO₂ readings and/or could disrupt other electronic equipment.
- Do not alter (such as folding the sensor cable) and do not make any modification to the sensor as this may damage the sensor. Damaged sensors may affect measurement accuracy and increase the risk of electric shock.
- Avoid misapplication or misuse of the sensor. Incorrect use may lead to inaccurate SpO₂ measurements or

skin damage such as redness or sores. Examples of incorrect applications include, but are not limited to:

- a. Applying the sensor too tightly or applying third-party accessories (too much pressure restricting blood flow).
 - b. Wrapping the sensor with another material.
 - c. Applying the sensor with the assistance of third-party tape or other adhesives.
 - d. Leaving the sensor in one place for longer than recommended. If you are unsure of correct sensor use, contact your clinician or Medtronic.
- To prevent skin reactions, check sensor site periodically and reposition every 4 to 6 hours. Avoid application of the sensors to areas with poor skin condition (such as reddened or damaged).
 - Circulation near the sensor site should be checked routinely. The site should be inspected every 4 to 6 hours for adhesion, application pressure, skin condition, and correct emitter and detector position. If skin condition changes, move the sensor to another site.
 - Check and adequately prepare the sensor site before placing the sensor. Incorrect selection of sensor site can cause inaccurate readings.
 - Check the sensor prior to using to ensure no damage to the sensor (such as exposed wires, emitter, and detector). Do not use a sensor that was opened or damaged during transit.
 - The sensor is designed for single patient use only and may not be used on multiple patients due to an increased risk of cross-contamination among patients.
 - This product cannot be cleaned and/or sterilized by the operator in order to safely reuse and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a skin reaction, infection, or product failure.

- Excessive motion may affect performance. In such cases, try to keep the patient still, or change the sensor site to one with less motion.
- Avoid excessive patient motion, for example during transport/vehicular motion. Such motion may affect sensor performance and accuracy.
- The operator is responsible for checking compatibility of monitor probe and cable before use.
- Certain physical conditions may affect measurement accuracy. These conditions include, but are not limited to poor circulation, skin pigmentation, skin thickness, skin temperature, current tobacco use or sickle cell disease, dyes injected into bloodstream, and use of externally applied coloring such as nail polish, sunscreen, dye, or pigmented cream. If you are unsure if any of these conditions apply to the patient, contact your clinician or Medtronic (see back page).
- Route all cables carefully to reduce the possibility of patient entanglement or strangulation.
- Keep the sensor out of reach of children and pets, to avoid accidents such as choking. Keep in an area protected from children, pets and pests to avoid damage to the sensor.
- Do not use the sensor beyond the listed Use By Date located on the pouch label.
- Do not use the sensor on an area with poor blood circulation.
- Check the sensor prior to using to ensure no adhesive discoloration has occurred.
- Failure to apply the OxySoftNHC sensor properly may cause incorrect measurements. See Instructions for Use section for proper application procedure.
- Avoid possible disruptions from electronic equipment such as, but not limited to: Cellular phones, radio transmitters, motors, telephones, lamps, electrosurgical units, defibrillators, and other devices. Disruption may cause inaccurate SpO₂ measurements.

If you are unsure if your device is working properly, contact your clinician.

- If the sensor does not track the pulse reliably, the sensor site may be too thick or too thin, or the sensor may be incorrectly positioned. If this situation occurs, reposition the sensor or choose an alternate recommended sensor site.
- Do not use the sensor around flammable substances. This may cause explosion or fire.
- While the OxySoftNHC sensor is designed to reduce the effects of ambient light, excessive light may cause inaccurate measurements. In such cases, cover the sensor with opaque material (such as a blanket).
- High oxygen levels is associated with retinopathy of prematurity in premature infants. Therefore, the upper alarm limit for oxygen saturation must be carefully selected in accordance with accepted clinical standards and considering the accuracy range of the oximeter being used.
- Do not operate or store the monitoring system outside the environmental ranges listed below. Doing so may cause it to work incorrectly.
- No servicing or maintenance required for this device.

Cautions (What You Must Do to Avoid Other Harm)

- Do not allow patient to lay on the sensor as it can cause bruising and/or skin damage or redness.
- The sensor is designed to be used and disposed after 3 days of use.
- Do not use this device without a prescription from a licensed doctor or clinician.
- Sensor is packaged as sterile. In the event of damage to the sterile packaging, do NOT re-sterilize. Follow local governing ordinances and recycling instructions regarding disposal or recycling of sensors.

- Medical equipment applied to the patient (for example, blood pressure monitor) can change the way blood flows, affecting sensor function.

Note: The materials used in the manufacture of the sensor do not contain natural rubber latex. The materials have undergone extensive testing to meet ISO 10993 to assess potential reaction to the skin. Further information is available upon request.

Instructions for Use

For a visual demonstration of sensor placement visit <http://www.medtronic.com/assetviewer.html?path=/content/dam/covidien/library/us/en/product/pulse-oximetry/pulse-oximetry-oxysoft-application-video.mp4>

To apply the OxySoftNHC sensor:

1. Select sensor application site. Choose a location away from other medical devices that may change blood flow (such as a blood pressure monitor or tube inserted into bloodstream).

Neonates: The preferred sensor site is a foot. **2**

Consult your clinician if an alternative site is needed. Consider how the cable will be positioned.

Adults: The preferred sensor site is an index finger **5**. Alternatively, other fingers may be used.

2. Clean and remove any substances such as nail polish from the sensor application site. This can affect the sensor's accuracy.
3. Remove the plastic backing from the sensor and locate the optical components beneath the sun (☀) and dot (●) icons. Note: **1** The sun icons are over the emitter (light) and the dot icons are over the detector.
Note: When selecting a sensor site, priority should be given to a hand or foot that is away from a tube inserted into bloodstream or blood pressure monitor.
4. Place the sun (☀) icons on the sensor application site and wrap the sensor around the site so the dot (●) icons oppose the sun icons. **2 4**

Neonates: Orient the sensor based on the desired cable position, either along the side of the leg or arm, or away from the patient. **3**

Adults: Place the sun (☀) icons on the nail bed. Do not place on a joint. Note that the cable will be positioned on the left side of the selected finger. **5**

5. Wrap the length of the bandage firmly around the sensor application site, but not too tightly. The sun and dot icons must oppose each other.
6. Plug the sensor into the oximeter and verify proper operation as described in the oximeter operator's manual.

Note: Oximeters could require a sensor extension cable.

To remove the OxySoftNHC Sensor:

1. Peel the blue end to gently remove the sensor.
2. Disconnect the sensor from the monitor by unplugging the sensor from the extension cable or monitor.

Reapplication of the OxySoftNHC Sensor:

When reapplying to the same patient, refer to application procedure above. It is recommended that the site be checked every 4-6 hours for reapplication. The sensor can be reused on the same patient for up to 3 days of use life if the sensor adhesive attaches without slippage.

Compatible Boards*	Compatible extension cables
PMB05N	DEC-4
10046102	DEC-8
PMB10N	DEC10U
PMC10CV283N	DOC-4, DOC-8
PMC10UB305N	DOC-10
Nell1-SR, Nell1	MC10, MC10-P2

Compatible Boards*	Compatible extension cables
Nell1MMS	MCU10-10
Nell2	OEM10
Nell3-SR, Nell3	OEM10-U
MP404	SCP10

*See monitor manufacturer for specific board details.

Note: For additional warnings, cautions or contraindications when using this sensor with Nellcor™-compatible instruments, refer to the system's operator's manual or contact the manufacturer of the system.

Accuracy Specifications

This sensor has Nellcor™ OxiMax™ technology in its design. When connected to an OxiMax™-enabled instrument, this sensor uses OxiMax™ technology to provide additional advanced sensor performance features. If not using with a Nellcor monitor, consult monitor manufacturer for features and compatibility.

Each Nellcor™ OxiMax™-compatible instrument manufacturer is responsible for determining ideal compatibility settings to provide safe and effective use of each Nellcor™ OxiMax™ sensor. This includes specifications and/or warnings, cautions, or contraindications. Refer to each system's operator's manual or consult manufacturer for complete instructions, warnings, cautions, or contraindications on the use of this sensor with their Nellcor™ OxiMax™ compatible instrument.

Note: An SpO₂ laboratory tester cannot be used to assess the accuracy of the OxySoftNHC sensor.

Sensor Saturation Accuracy	
Poor Blood Circulation ¹	70 to 100% $\pm$ 2 digits
Adult and Neonate ^{2,3,4}	70 to 100% $\pm$ 2 digits
Adult and Neonate with Motion ⁵	70 to 100% $\pm$ 3 digits
Low Oxygen Saturation ²	60 to 80% $\pm$ 3 digits

Pulse Rate Accuracy⁶	
Normal Conditions	20 – 250 BPM: $\pm$ 3 BPM
Poor Blood Circulation	20 – 250 BPM: $\pm$ 3 BPM
Motion	20 – 250 BPM: $\pm$ 5 BPM
(refer to monitor's operating manual for more information)	

¹ Reading accuracy in the presence of poor blood circulation was confirmed using signals supplied by a patient simulator. SpO₂ and pulse rate values were varied across the monitoring range over a range of weak signal conditions and compared to the known blood oxygen level and pulse rate of the input signals. These values were confirmed on the following monitors: Nellcor™ PM1000N, N600x, OxiCable

² Confirmed by clinical studies, where the measured values of the sensors were compared with known values using a different blood oxygen measurement method in healthy adult subjects over the specified oxygen saturation range. These values were confirmed on the following monitors: Nellcor™ PM1000N, N600x, OxiCable

³ Accuracy was demonstrated in a low oxygen supply study with 10 subjects, 50% were male and 50% female, and a range of skin pigmentation (20% with darker skin and 80% with lighter skin). The observed accuracy for each subgroup was within stated device accuracy.

- 4 Clinical function of the OxySoftN sensor has been demonstrated on hospitalized neonate patients. The observed SpO₂ accuracy was 2.68% in a study of 27 patients with ages of 1 to 9 days, weight from 640 to 4710 grams, where 64 samples of known oxygen level (SaO₂) was 82.1% to 99.3%.
- 5 Motion performance was confirmed during a controlled low oxygen blood study with healthy adult subjects over the specified oxygen saturation range. Subjects performed rubbing and tapping movements 1-2 cm in amplitude with distance with random spacing and speed between movements. These values were confirmed on the following monitor: Nellcor™ N600x
- 6 Pulse rate accuracy was confirmed by simulated laboratory tests. These values were confirmed on the following monitors: Nellcor™ PM1000N, N600x, OxiCable

Optical Specifications	
LED wavelength (Approximate)	Output power
Red: 650 to 670 nm	10.0 mW
Infrared: 857 to 910 nm	10.0 mW

Environmental Specifications	
Operational temperature range	0°C to 40°C
Transport and Storage between usage temperature range	-40°C to 70°C
Transient operational temperature range	-20°C to 50°C
Transport and Storage atmospheric pressure range	500 hPa to 1060 hPa

Environmental Specifications	
Operational atmospheric pressure range	620 hPa to 1060 hPa
Operational and Storage Relative humidity range	15% to 95% non-condensing

Note: These environmental specifications are specific to the sensor. Ensure sensor application within limited environmental conditions for safe performance. For information on the environmental specification of the compatible monitoring systems, see the operator's manual for that system.

Cleaning/Disinfection/Sterilization

This product cannot be cleaned and/or sterilized by the operator in-order to safely reuse and is therefore intended for single use. Attempts to clean or sterilize this device may result in skin reaction, infection, or product failure.

Disposal

Do not dispose with domestic waste. Follow local regulations when disposing or recycling this device and its component parts. A list of component materials is available by contacting Medtronic at 800.255.6774, option 1, or a local Medtronic representative.

Troubleshooting:

- Monitor the sensor site while the sensor is connected to the patient.
- Do not use the sensor past the 3-day use life period.
- Do not attempt to clean or sterilize the sensor.
- Check to see that leads are connected, wires are not damaged, and the power is supplied to the device.
- The sensor will emit a red light when working correctly.

- Contact Medtronic to report unexpected operation or events.
- For additional help, contact Medtronic at 800.255.6774. Local support information can be found at www.medtronic.com/covidien/en-us/support.html

Adverse Events

Adverse events associated with the use of this device, in order of increasing severity are electric shock, low oxygen saturation, high oxygen saturation, and death. Any serious incident that occurs in relation to the device or its accessories should be reported to Medtronic and to the national competent authority.

Additional Copies of Instructions

Additional copies of these Instructions are available at no charge by calling Medtronic or its authorized distributors. Also, permission is hereby granted under Medtronic copyrights to purchasers of products obtained from Medtronic or its authorized distributors to make additional copies of these instructions for use by such purchasers.

California Proposition 65 Statement

For information regarding California Proposition 65, please refer to www.covidien.com/caprop65.

Clinical Studies

Refer to the graphs (Bland-Altman) at the end of this booklet.

Additional graphs are available in the IFU.

NOTE: FDA-cleared prescription pulse oximeters are required to have a minimum average (mean) accuracy that is demonstrated by reduced oxygen saturation studies done on healthy patients. This testing compares the pulse oximeter saturation readings to known blood

oxygen level readings for values between 70–100%. The typical accuracy (reported as Accuracy Root Mean Square or A_{RMS}) of recently FDA-cleared pulse oximeters is within 2 to 3% of known blood oxygen values. This generally means that during testing, about 66% of SpO_2 values were within 2 or 3% of known known blood oxygen values and about 95% of SpO_2 values were within 4 to 6% of known blood gas values, respectively.

Figure 1 Accuracy graph (Bland Altman Plot): Covidien Nellcor™ monitor (**N-600x**)/**OXYSOFTN**⁴

Number on graph:

- | | | |
|---|-----------------|--|
| 1 | Vertical axis | Compared OxySoftN blood oxygen reading to known blood oxygen level ($SpO_2 - SaO_2$) |
| 2 | Horizontal axis | Known blood oxygen level used for comparison (SaO_2) |
| 3 | | 95% of measurements appear below this line |
| 4 | | Minimum average accuracy (Mean) |
| 5 | | 95% of measurements appear above this line |
| 6 | | The typical accuracy (reported as Accuracy Root Mean Square or A_{RMS}) |

Figure 2 Accuracy graph (Bland Altman Plot): Covidien Nellcor™ monitor (**N-600x**)/**OXYSOFTN Motion**⁵

Number on graph:

- | | | |
|---|-----------------|--|
| 1 | Vertical axis | Compared OxySoftN blood oxygen reading to known blood oxygen level ($SpO_2 - SaO_2$) |
| 2 | Horizontal axis | Known blood oxygen level used for comparison (SaO_2) |
| 3 | | 95% of measurements appear below this line |
| 4 | | Minimum average accuracy (Mean) |
| 5 | | 95% of measurements appear above this line |
| 6 | | The typical accuracy (reported as Accuracy Root Mean Square or A_{RMS}) |

⁴Clinical accuracy study, conducted with 10 test subjects

⁵Clinical motion accuracy study conducted with 14 test subjects

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

STERILE EO	Sterile using ethylene oxide
	Single use
	Do not use if package is damaged
	Consult instructions for use
	Caution
	____ hPa Atmospheric pressure limitation
	____ % Humidity limitation
	____ °F ____ °C Temperature limitation
	Keep dry
REF	Catalogue number
	Manufacturer
	Use by date

Page Left Intentionally Blank

Part No. PT00158241 Rev A 2022-05

Medtronic and Medtronic logo are trademarks of Medtronic.

^{TM*} Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.

Covidien is a Medtronic company.

© 2020 Medtronic.

 Covidien llc, 15 Hampshire Street,
Mansfield, MA 02048 USA.

Medtronic B.V., Earl Bakkenstraat 10,
6422 PJ Heerlen, The Netherlands.

www.medtronic.com

Rx
ONLY

C € 0123