

Nellcor™ Adult SpO₂ Sensor Non-Adhesive

en	Instructions For Use
fr	Mode d'emploi
de	Gebrauchsanweisung
nl	Gebruiksaanwijzing
es	Instrucciones de uso
it	Istruzioni per l'uso
sv	Bruksanvisning
da	Brugsanvisning
pt	Instruções de uso
ru	Инструкция по применению
zh	使用说明

[REF] SC-A
[REF] SC-A-I

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

This device is not made with natural rubber latex or DEHP.

Directions for Use

Indications/Contraindications

The Nellcor™ Adult SpO₂ Sensor Non-Adhesive, model SC-A, is indicated for single patient use when continuous noninvasive arterial oxygen saturation and pulse rate monitoring are required for adults weighing more than 40 kg.

Use this sensor only with Nellcor instruments and instruments containing Nellcor OxiMax technology, or with instruments licensed to use Nellcor OxiMax sensors (Nellcor-compatible instruments). When connected to an OxiMax-compatible instrument, this sensor uses OxiMax technology to provide additional advanced sensor performance features. Consult individual manufacturers for features and compatibility of particular instruments and sensor models.

Each manufacturer of Nellcor-compatible instruments is responsible for determining whether and under what conditions its instruments are compatible for safe and effective use with each Nellcor sensor model. This may include different specifications and/or warnings, cautions, or contraindications. Refer to the instrument operator's manual or consult manufacturer for complete instructions for use of this sensor with their Nellcor-compatible instrument.

Note: Non-adhesive sensor consists of sensor, cable, wrap and loop fastener.

Instructions for Use



- 1) Locate transparent windows on the sensor. Windows cover optical components. This is the side that contacts the patient. Note corresponding alignment marks (a) on the outside and five horizontal lines (b) between the marks **A**.
- 2) Orient the SC-A so the horizontal lines are on the lateral edge of the site:

The preferred site is an index finger. Alternatively, other fingers may be used. The window next to the cable goes on the nail side, distal to the first joint. Do not place on a joint. Note that the cable must be positioned on the *top* of the hand **B**.

Note: When selecting a sensor site, priority should be given to an extremity free of an arterial catheter, blood pressure cuff, or intravascular infusion line.

- 3) Wrap the SC-A firmly, but not too tightly around the finger. Windows must oppose each other. Secure the hook and loop fastener.
- 4) Use the cable wrap—attached to the sensor cable—to secure the sensor cable to the patient **C**. Wrap it around the finger and secure with the hook and loop fastener.
- 5) Plug the SC-A into the oximeter and verify proper operation as described in the oximeter operator's manual.

Note: If the sensor does not track the pulse reliably, it may be incorrectly positioned—or the sensor site may be too thick, thin, or deeply pigmented, or otherwise deeply colored (for example, as a result of externally applied coloring such as nail polish, dye, or pigmented cream) to permit appropriate light transmission. If any of these situations occurs, reposition the sensor or choose an alternate Nellcor sensor for use on a different site.

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Vous pouvez vous procurer gratuitement des copies supplémentaires de ces instructions auprès de Covidien ou de ses représentants agréés. Par ailleurs, en vertu des droits de copyright de Covidien, les acheteurs de produits obtenus auprès de Covidien ou de ses distributeurs agréés sont autorisés à effectuer des copies supplémentaires de ce mode d'emploi à usage personnel.
COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company.
Covidien Inc, 15 Hampshire Street,
Covidien Ireland Limited,