

Nellcor™

Neonatal SpO₂ Sensor

Non-Adhesive

REF SCNEO	en Instructions For Use
REF SCNEO-I	fr Mode d'emploi
	de Gebruuchsanweisung
	nl Gebruiksaanwijzing
	es Instrucciones de uso
	it Istruzioni per l'uso
	sv Bruksanvisning
	da Bruksanvisning
	pt Instruções de uso
	ru Инструкция по применению
	zh 使用说明

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

This device is not made with natural rubber latex or DEHP.

Neonatal SpO₂ Sensor Non-Adhesive

Directions for Use

Indications/Contraindications	Warnings
The Nellcor™ Neonatal SpO₂ Sensor Non-Adhesive, model SCNEO, is indicated for single patient use when continuous noninvasive arterial oxygen saturation and pulse rate monitoring are required for neonates weighing 1.5 kg to 5 kg.	- Failure to apply the SCNEO properly may cause incorrect measurements. - High oxygen levels may predispose a premature infant to develop retinopathy. Therefore, the upper alarm limit for oxygen saturation must be carefully selected in accord with accepted clinical standards and considering the accuracy range of the oximeter being used. - While the SCNEO is designed to reduce the effects of ambient light, excessive light may cause inaccurate measurements. In such cases, cover the sensor with opaque material. - Circulation distal to the sensor site should be checked routinely. The site must be inspected every 8 hours to ensure skin integrity, and correct optical alignment. If skin integrity changes, move the sensor to another site. - Intravascular dyes or externally applied coloring such as nail polish, dye, or pigmented cream, may lead to inaccurate measurements. - Excessive motion may compromise performance. In such cases, try to keep the patient still, or change the sensor site to one with less motion. - Do not immerse in water or cleaning solutions. Do not resterilize. - If the sensor is wrapped too tightly or supplemental tape is applied, venous pulsations may lead to inaccurate saturation measurements. - Do not use the SCNEO or other oximetry sensors during MRI scanning. Conducted current may cause burns. Also, the SCNEO may affect the MRI image, and the MRI unit may affect the accuracy of oximetry measurements. - Do not alter or modify the SCNEO. Alterations or modifications may affect performance or accuracy. - For additional warnings, cautions or contraindications when using this sensor with Nellcor-compatible instruments, refer to the instrument operator's manual or contact the manufacturer of the instrument. - In the event of damage to the sterile packaging, do NOT resterilize. Follow local governing ordinances and recycling instructions regarding disposal or recycling of sensors. - If the sensor is misapplied with excessive pressure for prolonged periods, a pressure injury can occur.
Instructions for Use	Accuracy Specifications
- Locate transparent windows on the sensor. Windows cover optical components. This is the side that contacts the patient. Note corresponding alignment marks (a) on the outside (b) between the marks A. - Orient the SCNEO so the graphics are on the lateral edge of the site: The preferred site is a foot. Alternatively, use a hand. The window next to the cable goes on the sole of the foot as shown B. Note: When selecting a sensor site, priority should be given to an extremity free of an arterial catheter, blood pressure cuff, or intravascular infusion line. - Wrap the SCNEO firmly, but not too tightly around the foot. Windows must oppose each other. Secure the hook and loop fastener. - Use the cable wrap—attached to the sensor cable—to secure the sensor cable to the patient's ankle. C Secure with the hook and loop fastener. - Plug the SCNEO into the oximeter and verify proper operation as described in the oximeter operator's manual. Note: If the sensor does not track the pulse reliably, it may be incorrectly positioned—or the sensor site may be too thick, thin, or deeply pigmented, or otherwise deeply colored (for example, as a result of externally applied coloring such as nail polish, dye, or pigmented cream) to permit appropriate light transmission. If any of these situations occurs, reposition the sensor or choose an alternate Nellcor sensor for use on a different site.	For the accuracy specification range when used with Nellcor monitors, refer to information provided with the monitor, or (in the U.S.), contact Covidien's Technical Services Department. Outside the U.S., contact your local Covidien representative. For the accuracy specification range of this sensor when used with Nellcor-compatible instruments, refer to the instrument operator's manual or contact the instrument manufacturer.
Reapplication	Additional Copies of Instructions
The SCNEO can be reused <i>on the same patient</i> as long as the sensor attaches without slippage.	Additional copies of these Instructions are available at no charge by calling Covidien or its authorized distributors. Also, permission is hereby granted under Covidien copyrights to purchasers of products obtained from Covidien or its authorized distributors to make additional copies of these instructions for use by such purchasers. Any other use of this sensor is not authorized by Covidien under any patent.

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COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company.
Covidien Ireland Limited, IDA Business & Technology Park, Tullamore.