LISA TRACKER Duo Adalimumab

REF LTA 005

Français

DEFINITION

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BIBLIOGRAPHIE

Ainsworth MA & al. Tumor necrosis factor-alpha binding capacity and anti-infliximab antibodies measured by fluidphase radioimmunoassays as predictors of clinical efficacy of infliximab in Crohn's disease. *Am J Gastroenterol.* 2008 Apr;103(4):944-8.

Afif W & al. Clinical Utility of Measuring Infliximab and Human Anti-Chimeric Antibody Concentrations in Patients With Inflammatory Bowel Disease. *Am J Gastroenterol* 2010; 105:1133–1139.

Arends S & al. The formation of autoantibodies and antibodies to TNF- α blocking agents in relation to clinical response in patients with ankylosing spondylitis. *Clin Exp Rheumatol.* 2010 Sep-Oct;28(5):661-8.

Assa A & al. Proactive Monitoring of Adalimumab Trough Concentration Associated With Increased Clinical Remission in Children With Crohn's Disease Compared With Reactive Monitoring *Gastroenterology* . 2019 Oct;157(4):985-996.e2.

Bartelds GM & al. Development of antidrug antibodies against adalimumab and association with disease activity and treatment failure during long-term follow-up. *JAMA.* 2011 Apr 13;305(14):1460-8.

Bendtsen K. Is There a Need for Immunopharmacologic Guidance of Anti-Tumor Necrosis Factor Therapies. *ARTHRITIS & RHEUMATISM* Vol. 63, No. 4, April 2011, pp 867–870.

Bressler B & al. Clinical Practice Guidelines for the Medical Management of Nonhospitalized Ulcerative Colitis: The Toronto Consensus. *Gastroenterology* 2015;148:1035 – 1058.

Choon Jin Ooi & al. Best practices on immunomodulators and biologic agents for ulcerative colitis and Crohn's disease in Asia. *Intest Res,* May 31, 2019:1-26.

Choy & al. Efficacy of a novel PEGylated humanized anti-TNF (CDP870) in patients with rheumatoid arthritis : a phase II double-blinded, randomized, dose escalating trial. *Rheumatology* 2002;41:1133-1137.

Clinical Guideline: Management of Crohn's Disease in Adults. *Am J Gastroenterol* 2018; 113:481–517.

Desroches M & al. Treatment failure with antagonists of TNF- α : mechanisms and implications for the care of patients. *Eur. Cytokine Netw.,* Vol. 21 n° 4, December 2010, 226-31.

Diagnostics Guidance [DG22] on Therapeutic monitoring of TNF-alpha Inhibitors in Crohn's Disease.
<https://www.nice.org.uk/guidance/dg22/chapter/1-Recommendations>

ECCO-ESGAR Guideline for Diagnostic Assessment in IBD Part 1: Initial diagnosis, monitoring of known IBD, detection of complications. *Journal of Crohn's and Colitis,* 2018, 1–32.

Feuerstein J D. & al. Therapeutic Drug Monitoring in Inflammatory Bowel Disease. *Gastroenterology.* 2017:1-8.

Flanagan E. & al. Infliximab, Adalimumab and Vedolizumab concentrations across pregnancy and vedolizumab concentrations in infants following intrauterine exposure *Aliment Pharmacol Ther.* 2020 Sep 27.

Gibson DJ & al. Review article: determination of the therapeutic range for therapeutic drug monitoring of adalimumab and infliximab in patients with inflammatory bowel disease. *Aliment Pharmacol Ther.* 2020 Mar;51(6):612-628.

Gomollón F & al. 3rd European Evidence-based Consensus on the Diagnosis and Management of Crohn's Disease. *J Crohns Colitis.* 2016:1-23.

Greuter T & al. Therapeutic Drug Monitoring to Guide Clinical Decision Making in Inflammatory Bowel Disease Patients with Loss of Response to Anti-TNF: A Delphi Technique-Based Consensus. *Digestion.* 2019 Aug 28:1-9.

Grossi V & al. Anti-adalimumab and anti-certolizumab antibodies titers after discontinuation of adalimumab: two case reports. *Clin Chem Lab Med.* 2019 Oct 30.

Hässler S & al. Clinicogenomic factors of biotherapy immunogenicity in autoimmune disease: A prospective multicohort study of the ABIRISK consortium *PLoS Med.* 2020 Oct 30;17(10):e1003348

Jaminitski & al. The presence or absence of antibodies to infliximab or adalimumab determines the outcome of switching to etanercept. *Rheum Dis.* 2011 Feb;70(2):284-8.

Kennedy NA & al. Predictors of anti-TNF treatment failure in anti-TNF-naïve patients with active luminal Crohn's disease: a prospective, multicentre, cohort study. *Lancet Gastroenterol Hepatol.* 2019 May;4(5):341-353..

Khan A & al. New Zealand Society of Gastroenterology Guidelines on Therapeutic Drug Monitoring in Inflammatory Bowel Disease. *N Z Med J.* 2019 Mar 8;132(1491):46-62.

Koren & al. Recommendation on risk-based strategies for detection and characterization of antibodies against biotechnology products. *Journal of Immunological Methods,* 333 (2008) 1-9.

Korswagen LA & al. Venous and Arterial Thromboembolic Events in Adalimumab-Treated Patients With Anti-adalimumab Antibodies. *ARTHRITIS & RHEUMATISM* Vol. 63, No. 4, April 2011, pp 877–883.

Labetoulle R & al. Prolonged Persistence of Adalimumab Transferred From Mother to Infant During Pregnancy. *Ann Intern Med.* 2018 Mar 6.

Lamb CA & al. British Society of Gastroenterology consensus guidelines on the management of inflammatory bowel disease in adults. *Gut* 2019;0:1–106.

Lambrescak E & al. Duration of combination therapy and risk of treatment failure in patients with inflammatory bowel disease *Clin Res Hepatol Gastroenterol.* 2020 Sep 3

Little RD & al. Comparison of adalimumab serum drug levels when delivered by pen versus syringe in patients with inflammatory bowel disease. An international, multicentre cohort analysis. *J Crohns Colitis.* 2019 May 15.

Mahil SK & al. Predicting treatment response in psoriasis using serum levels of adalimumab and etanercept: a single-centre, cohort study. *Br J Dermatol.* 2013 Aug;169(2):306-13.

Management of Paediatric Ulcerative Colitis, Part 1: Ambulatory Care—An Evidence-based Guideline. *J Pediatr Gastroenterol Nutr*. 2018 May 30.

Marcus Harbord & al. Third European Evidence-based Consensus on Diagnosis and Management of Ulcerative Colitis. Part 2: Current Management. *Journal of Crohn's and Colitis*, 2017, 1–24.

Martelli L & al. Cost-effectiveness of drug monitoring of anti-TNF therapy in inflammatory bowel disease and rheumatoid arthritis: a systematic review. *J Gastroenterol*. 2017 Jan;52(1):19-25

Marotte H & al. Tapering Without Relapse in Rheumatoid Arthritis Patients With High TNF Blocker Concentrations: Data from STRASS Study *Ann Rheum Dis* . 2020 Jul;79(7):e81..

Mire-Sluis & al. Recommendation for the design and optimization of immunoassays used in the detection of host antibodies against biotechnology products. *Journal of Immunological Methods*, 289 (2004) 1-16.

Mitrev N & al. Consensus statements on therapeutic drug monitoring of anti-tumour necrosis factor therapy in inflammatory bowel diseases. *Aliment Pharmacol Ther*. 2017 Dec;46(11-12):1037-1053.

Papamichael K & al. Appropriate Therapeutic Drug Monitoring of Biologic Agents for Patients With Inflammatory Bowel Diseases. *Clin. Gastroenterology and Hepatology*, 24 March 2019

Papamichael K & al. Therapeutic drug monitoring in patients on biologics: lessons from gastroenterology *Curr Opin Rheumatol*. 2020 Jul;32(4):371-379.

Paul S & al. No Difference of Adalimumab Pharmacokinetics when dosed at 40mg Every Week or 80mg Every Other Week in IBD Patients in Clinical Remission after Adalimumab Dose Intensification *Dig Dis Sci* . 2020 Sep 16

Peyrin-Biroulet L & al. French National Consensus Clinical Guidelines for the Management of Crohn's disease. *Dig Liver Dis*. 2016.

Peyrin-Biroulet L & al. French National Consensus Clinical Guidelines for the Management of Ulcerative Colitis. *Dig Liver Dis*. 2016; 48(7):726-33.

Quistbert J & al. Incidence and risk factors for adalimumab and infliximab anti-drug antibodies in rheumatoid arthritis: A European retrospective multicohort analysis. *Semin Arthritis Rheum*. 2018 Oct 12. pii: S0049-0172(18)30176-8.

Radstake TRDJ & al. Formation of antibodies against infliximab and adalimumab strongly correlates with functional drug levels and clinical responses in rheumatoid arthritis. *Ann Rheum Dis* 2009;68:1739–1745.

Ram Raj Singh, & al. TNF α blockade in human diseases : mechanisms and future. *Clin.Immunol.*, 2008 February ; 126(2):121-136.

Roblin X & al. Association between pharmacokinetics of adalimumab and mucosal healing in patients with inflammatory bowel diseases. *Clin Gastroenterol Hepatol*. 2014 Jan;12(1):80-84.e2.

Roblin X & al. Development of an algorithm incorporating pharmacokinetics of adalimumab in inflammatory bowel diseases. *Am J Gastroenterol*. 2014 Aug;109(8):1250-6.

Roblin X & al. Addition of Azathioprine to the Switch of anti-TNF in Patients With IBD in Clinical Relapse With Undetectable anti-TNF Trough Levels and Antidrug Antibodies: A Prospective Randomised Trial *Gut*. 2020 Jan 24.

Ruemmele F.M. & al. Consensus guidelines of ECCO/ESPGHAN on the medical management of pediatric Crohn's disease. *ournal of Crohn's and Colitis* (2014).

Sparado A & al. Switching from infliximab or etanercept to adalimumab in resistant or intolerant patients with spondyloarthritis: a 4-year study. *Rheumatology (Oxford)*. 2010 Jun;49(6):1107-11.

Therapeutic Drug Monitoring in Inflammatory Bowel Disease. *Gastroenterology*. 2017;153(3):835-857.

Therapeutic Drug Monitoring in Inflammatory Bowel Disease. *Gastroenterology*. 2017;153:858-859.

3rd European Evidence-based Consensus on the Diagnosis and Management of Crohn's Disease. *J Crohns Colitis*. 2017; 11(1):3-25.

3rd European Evidence-based Consensus on the Diagnosis and Management of Ulcerative Colitis. *J Crohns Colitis*. 2017; 11(6):649-670.

Ungar B & al. Addition of an immunomodulator can reverse antibody formation and loss of response in patients treated with adalimumab. *Aliment Pharmacol Ther*. 2017 Jan;45(2):276-282.

Ungar B et al. Prospective Observational Evaluation of Time-Dependency of Adalimumab Immunogenicity and Drug Concentrations: The Poetic Study. *Am J Gastroenterol*. 2018 Jun;113(6):890-898

Ward MG & al. Infliximab and adalimumab drug levels in Crohn's disease: contrasting associations with disease activity and influencing factors. *Aliment Pharmacol Ther*. 2017 Jul;46(2):150-161.

World Health Organization, WHO Expert Committee on Biological Standardization, Wadhwa Meenu, Bird Chris, Atkinson Eleanor. & al. (2019). Report on a collaborative study for proposed 1st international standard for Adalimumab. WHO/BS/2019.2365.

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Dosage Adalimumab	Dosage anti- Adalimumab
----------------------	----------------------------

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