LISA TRACKER[®]

Infliximab

REF LTI 002-48

DEFINITION

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BIBLIOGRAPHIE

Ainsworth MA & al. Tumor necrosis factor-alpha binding capacity and anti-infliximab antibodies measured by fluidphase radioimmunoassays as predictors of clinical efficacy of infliximab in Crohn's disease. *Am J Gastroenterol.* 2008 Apr;103(4):944-8.

Afif W & al. Clinical Utility of Measuring Infliximab and Human Anti-Chimeric Antibody Concentrations in Patients With Inflammatory Bowel Disease. *Am J Gastroenterol* 2010; 105:1133–1139.

Arends S & al. The formation of autoantibodies and antibodies to TNF- α blocking agents in relation to clinical response in patients with ankylosing spondylitis. *Clin Exp Rheumatol.* 2010 Sep-Oct;28(5):661-8.

Attar A & al. Cost savings using a test-based de-escalation strategy for patients with Crohn's disease in remission on optimized infliximab: A discrete event model study. *Dig Liver Dis.* 2019 Jan;51(1):112-119

Avouac J & al. Systematic switch from innovator infliximab to biosimilar infliximab in inflammatory chronic diseases in daily clinical practice: The experience of Cochin University Hospital, Paris, France. *Semin Arthritis Rheum.* 2017 Oct 5.

Baert F & al. Influence of Immunogenicity on the Long-Term Efficacy of Infliximab in Crohn's Disease. *N Engl J Med* 2003;348:601-8.

Ben-Horin S & al. Undetectable anti-TNF drug levels in patients with long-term remission predict successful drug withdrawal. *Aliment Pharmacol Ther.* 2015 Aug;42(3):356-64.

Bendtsen K. Is There a Need for Immunopharmacologic Guidance of Anti-Tumor Necrosis Factor Therapies. *ARTHRITIS & RHEUMATISM* Vol. 63, No. 4, April 2011, pp 867–870.

Benucci M & al. Correlation between HLA haplotypes and the development of antidrug antibodies in a cohort of patients with rheumatic diseases. *Biologics: Targets and Therapy* 2018;12 37–41

Benucci M & al. Safety, efficacy and immunogenicity of switching from innovator to biosimilar infliximab in patients with spondyloarthritis: a 6-month real-life observational study. *Immunol Res.* 2016 Jul 23.

Bor R & al. Clinical role, optimal timing and frequency of serum infliximab and anti-infliximab antibody level measurements in patients with inflammatory bowel disease. *PLoS One.* 2017 Mar 31;12(3):e0172916

Buhl S & al. Therapeutic thresholds and mechanisms for primary non-response to infliximab in inflammatory bowel disease *Scand J Gastroenterol.* 2020 Aug;55(8):884-890

Bressler B & al. Clinical Practice Guidelines for the Medical Management of Nonhospitalized Ulcerative Colitis: The Toronto Consensus. *Gastroenterology* 2015;148:1035 – 1058.

Candon S & al. Clinical and biological consequences of immunization to infliximab in pediatric Crohn's disease. *Clin Immunol.* 2006 Jan;118(1):11-9.

Choon Jin Ooi & al. Best practices on immunomodulators and biologic agents for ulcerative colitis and Crohn's disease in Asia. *Intest Res,* May 31, 2019:1-26.

Choy & al. Efficacy of a novel PEGylated humanized anti-TNF (CDP870) in patients with rheumatoid arthritis : a phase II double-blinded, randomized, dose escalating trial. *Rheumatology* 2002;41:1133-1137.

Clinical Guideline: Management of Crohn's Disease in Adults. *Am J Gastroenterol* 2018; 113:481–517.

Coutzac C & al. Association Between Infliximab Trough Levels and the Occurrence of Paradoxical Manifestations in Patients with Inflammatory Bowel Disease: a Case-Control Study. *J Crohns Colitis.* 2015 Nov;9(11):982-7.

Dave MP & al. Clinical efficacy of infliximab level and anti-infliximab antibody measurement in patients with inflammatory bowel disease: An audit *Indian J Gastroenterol.* 2020 Oct 29.

Desroches M & al. Treatment failure with antagonists of TNF- α : mechanisms and implications for the care of patients. *Eur. Cytokine Netw.,* Vol. 21 n° 4, December 2010, 226-31.

DeVries & al. Inefficacy of infliximab in ankylosing spondylitis is correlated with antibody formation. *Ann Rheum Dis.* 2007 Jan;66(1):133-4.

Diagnostics Guidance [DG22] on Therapeutic monitoring of TNF-alpha Inhibitors in Crohn's Disease.
<https://www.nice.org.uk/guidance/dg22/chapter/1-Recommendations>

ECCO-ESGAR Guideline for Diagnostic Assessment in IBD Part 1: Initial diagnosis, monitoring of known IBD, detection of complications. *Journal of Crohn's and Colitis,* 2018, 1–32.

Edrees & al. Anti-tumor necrosis factor (TNF) therapy in rheumatoid arthritis: correlation of TNF-alpha serum level with clinical response and benefit from changing dose or frequency of infliximab infusions. *Clin Exp Rheumatol.* 2005 Jul-Aug;23(4):469-74.

Fay S & al. The Association Between Drug Levels and Endoscopic Recurrence in Postoperative Patients with Crohn's Disease Treated with Tumor Necrosis Factor Inhibitors. *Inflamm Bowel Dis.* 2017 Aug 22.

Fernandes SR & al. Proactive Infliximab Drug Monitoring Is Superior to Conventional Management in Inflammatory Bowel Disease. *Inflamm Bowel Dis.* 2019 Jun 27

Feuerstein J D. & al. Therapeutic Drug Monitoring in Inflammatory Bowel Disease. *Gastroenterology.* 2017:1-8.

Flanagan E. & al. Infliximab, Adalimumab and Vedolizumab concentrations across pregnancy and vedolizumab concentrations in infants following intrauterine exposure *Aliment Pharmacol Ther.* 2020 Sep 27.

Gecse KB, & al. Efficacy and Safety of the Biosimilar Infliximab CT-P13 Treatment in Inflammatory Bowel Diseases: A Prospective, Multicentre, Nationwide Cohort. *J Crohns Colitis.* 2016 Feb;10(2):133-40.

Gibson DJ & al. Review article: determination of the therapeutic range for therapeutic drug monitoring of adalimumab and infliximab in patients with inflammatory bowel disease. *Aliment Pharmacol Ther.* 2020 Mar;51(6):612-628.

Gomollón F & al. 3rd European Evidence-based Consensus on the Diagnosis and Management of Crohn's Disease. *J Crohns Colitis.* 2016:1-23.

Gonczi L & al. Drug persistence and need for dose intensification to adalimumab therapy; the importance of therapeutic drug monitoring in inflammatory bowel diseases. *BMC Gastroenterol.* 2017 Aug 8;17(1):97

Gonczi L & al. Prediction of Short- and Medium-term Efficacy of Biosimilar Infliximab Therapy. Do Trough Levels and Antidrug Antibody Levels or Clinical And Biochemical Markers Play the More Important Role? *J Crohns Colitis.* 2016 Nov 12

Greuter T & al. Therapeutic Drug Monitoring to Guide Clinical Decision Making in Inflammatory Bowel Disease Patients with Loss of Response

to Anti-TNF: A Delphi Technique-Based Consensus. Digestion. 2019 Aug 28;1-9.

Guidi L & al. Therapeutic drug monitoring is more cost-effective than a clinically-based approach in the management of loss of response to infliximab in inflammatory bowel disease: an observational multi-centre study. *J Crohns Colitis*. 2018 May 31

Hässler S & al. Clinicogenomic factors of biotherapy immunogenicity in autoimmune disease: A prospective multicohort study of the ABIRISK consortium *PLoS Med*. 2020 Oct 30;17(10):e1003348

Jaminitski & al. The presence or absence of antibodies to infliximab or adalimumab determines the outcome of switching to etanercept. *Rheum Dis*. 2011 Feb;70(2):284-8.

Jentzer A & al. Short Communication: Evaluation of infliximab and anti-infliximab LISA-TRACKER immunoassays for the therapeutic drug monitoring of SB2 infliximab biosimilar. *Ther Drug Monit*. 2018 Sep 21

Kennedy NA & al. Predictors of anti-TNF treatment failure in anti-TNF-naïve patients with active luminal Crohn's disease: a prospective, multicentre, cohort study. *Lancet Gastroenterol Hepatol*. 2019 May;4(5):341-353..

Khan A & al. New Zealand Society of Gastroenterology Guidelines on Therapeutic Drug Monitoring in Inflammatory Bowel Disease. *N Z Med J*. 2019 Mar 8;132(1491):46-62.

Koren & al. Recommendation on risk-based strategies for detection and characterization of antibodies against biotechnology products. *Journal of Immunological Methods*, 333 (2008) 1-9.

Krzysiek R & al. Circulating Concentration of Infliximab and Response to Treatment in Ankylosing Spondylitis: Results From a Randomized Control Study. *Arthritis & Rheumatism (Arthritis Care & Research)* Vol. 61, No. 5, May 15, 2009, pp 569–57.

Lamb CA & al. British Society of Gastroenterology consensus guidelines on the management of inflammatory bowel disease in adults. *Gut* 2019;0:1–106.

Lambrescak E & al. Duration of combination therapy and risk of treatment failure in patients with inflammatory bowel disease *Clin Res Hepatol Gastroenterol*. 2020 Sep 3

Lauret A & al. Effects of successive switches to different biosimilars infliximab on immunogenicity in chronic inflammatory diseases in daily clinical practice. *Semin Arthritis Rheum*. 2020 Feb 19.

Lee MW & al. Comparison of infliximab drug measurement across three commercially available ELISA kits. *Pathology*. 2016 Aug 23

Management of Paediatric Ulcerative Colitis, Part 1: Ambulatory Care—An Evidence-based Guideline. *J Pediatr Gastroenterol Nutr*. 2018 May 30.

Management of Paediatric Ulcerative Colitis, Part 2: Acute Severe Colitis—An Evidence-based Consensus Guideline. *J Pediatr Gastroenterol Nutr*. 2018 May 30.

Marcus Harbord & al. Third European Evidence-based Consensus on Diagnosis and Management of Ulcerative Colitis. Part 2: Current Management. *Journal of Crohn's and Colitis*, 2017, 1–24.

Marotte & al. Circulating tumour necrosis factor-alpha bioactivity in rheumatoid arthritis patients treated with infliximab: link to clinical response. *Arthritis Res Ther*. 2005;7(1):R149-55.

Martelli L & al. Cost-effectiveness of drug monitoring of anti-TNF therapy in inflammatory bowel disease and rheumatoid arthritis: a systematic review. *J Gastroenterol*. 2017 Jan;52(1):19-25.

Mire-Sluis & al. Recommendation for the design and optimization of immunoassays used in the detection of host antibodies against

biotechnology products. *Journal of Immunological Methods*, 289 (2004) 1-16.

Mitrev N & al. Consensus statements on therapeutic drug monitoring of anti-tumour necrosis factor therapy in inflammatory bowel diseases. *Aliment Pharmacol Ther*. 2017 Dec;46(11-12):1037-1053.

Nasser Y & al. Comparison of Point-of-Care and Classical Immunoassays for the Monitoring Infliximab and Antibodies Against Infliximab in IBD. *Dig Dis Sci*. 2018 Jun 9

Neveu B & al. An in-Vitro Comparison of Four Different Immunoassays for the Monitoring of Infliximab Biosimilars Drug Levels *Clin Biochem*. 2020 Jan 23.

Papamichael K & al. Appropriate Therapeutic Drug Monitoring of Biologic Agents for Patients With Inflammatory Bowel Diseases. *Clin. Gastroenterology and Hepatology*, 24 March 2019.

Papamichael K & al. Therapeutic drug monitoring in patients on biologics: lessons from gastroenterology *Curr Opin Rheumatol*. 2020 Jul;32(4):371-379.

Paul S & al. Letter: immunogenicity of anti-TNF in elderly IBD patients. *Aliment Pharmacol Ther*. 2019 Aug;50(3):336

Paul S & al. Letter: infliximab de-escalation based on trough levels in patients with inflammatory bowel disease. *Aliment Pharmacol Ther*. 2015 Oct;42(7):939-40.

Paul S & al. Therapeutic drug monitoring of infliximab and mucosal healing in inflammatory bowel disease: a prospective study. *Inflamm Bowel Dis*. 2013 Nov;19(12):2568-76

Petitdidier N & al. Real-world use of therapeutic drug monitoring of CT-P13 in patients with inflammatory bowel disease *Clin Res Hepatol Gastroenterol*. 2020 Jan 7

Peyrin-Biroulet L & al. French National Consensus Clinical Guidelines for the Management of Crohn's disease. *Dig Liver Dis*. 2016.

Peyrin-Biroulet L & al. French National Consensus Clinical Guidelines for the Management of Ulcerative Colitis. *Dig Liver Dis*. 2016; 48(7):726-33.

Poullenot F & al. Severe endoscopic lesions are not associated with more infliximab fecal loss in acute severe ulcerative colitis. *Dig Liver Dis*. 2018 Jul 10.

Radstake TRDJ & al. Formation of antibodies against infliximab and adalimumab strongly correlates with functional drug levels and clinical responses in rheumatoid arthritis. *Ann Rheum Dis* 2009;68:1739–1745.

Ram Raj Singh, & al. TNFα blockade in human diseases : mechanisms and future. *Clin.Immunol.*, 2008 February ; 126(2):121-136.

Ruemmele F.M. & al. Consensus guidelines of ECCO/ESPGHAN on the medical management of pediatric Crohn's disease. *ournal of Crohn's and Colitis* (2014).

Roblin X & al. Combination of C-reactive protein, infliximab trough levels, and stable but not transient antibodies to infliximab are associated with loss of response to infliximab in inflammatory bowel disease. *J Crohns Colitis*. 2015 Jul;9(7):525-31.

Roblin X & al. Cost savings of anti-TNF therapy using a test-based strategy versus an empirical dose escalation in Crohn's disease patients who lose response to infliximab. *Journal of Market Access & Health Policy*, [S.I.], v. 3, oct. 2015.

Roblin X & al. Development and Internal Validation of a Model Using Fecal Calprotectin in Combination with Infliximab Trough Levels to

Predict Clinical Relapse in Crohn's Disease. *Inflamm Bowel Dis*. 2017 Jan;23(1):126-132.

Roblin X & al. Distinct Thresholds of Infliximab Trough Level Are Associated with Different Therapeutic Outcomes in Patients with Inflammatory Bowel Disease: A Prospective Observational Study. *Inflamm Bowel Dis*. 2017 Sep 22.

Roblin X & al. Addition of Azathioprine to the Switch of anti-TNF in Patients With IBD in Clinical Relapse With Undetectable anti-TNF Trough Levels and Antidrug Antibodies: A Prospective Randomised Trial *Gut*. 2020 Jan 24.

Schmitz EM & al. Therapeutic drug monitoring of infliximab: performance evaluation of three commercial ELISA kits. *Clin Chem Lab Med*. 2015 Nov 20.

Seow CH & al. Trough serum infliximab: a predictive factor of clinical outcome for infliximab treatment in acute ulcerative colitis. *Gut*. 2010 Jan;59(1):49-54.

Sparado A & al. Switching from infliximab or etanercept to adalimumab in resistant or intolerant patients with spondyloarthritis: a 4-year study. *Rheumatology (Oxford)*. 2010 Jun;49(6):1107-11.

Steenholdt C, & al. Cut-off levels and diagnostic accuracy of infliximab trough levels and anti-infliximab antibodies in Crohn's disease. *Scandinavian Journal of Gastroenterology*, March 2011, Vol. 46, N°3, Pages 310-318.

Takeuchi & al. Baseline tumour necrosis factor alpha levels predict the necessity for dose escalation of infliximab therapy in patients with rheumatoid arthritis. *Ann Rheum Dis*. 2011 Jul;70(7):1208-15.

Ternant D & al. An enzyme-linked immunosorbent assay for therapeutic drug monitoring of Infliximab. *Ther Drug Monit.*, 2006 April ; 28(2):169-174.

Therapeutic Drug Monitoring in Inflammatory Bowel Disease. *Gastroenterology*. 2017;153(3):835-857.

Therapeutic Drug Monitoring in Inflammatory Bowel Disease. *Gastroenterology*. 2017;153:858-859.

3rd European Evidence-based Consensus on the Diagnosis and Management of Crohn's Disease. *J Crohns Colitis*. 2017; 11(1):3-25.

3rd European Evidence-based Consensus on the Diagnosis and Management of Ulcerative Colitis. *J Crohns Colitis*. 2017; 11(6):649-670.

Ungar B & al. Infliximab therapy intensification upon loss of response: Is there an optimal trough level? *Dig Liver Dis*. 2019 Mar 5.

Valido A & al. Efficacy, immunogenicity and cost analysis of a systematic switch from originator infliximab to biosimilar CT-P13 of all patients with inflammatory arthritis from a single center. *Acta Reumatol Port*. 2019 Nov 22

Van den Bemt BJF & al. Anti-infliximab antibodies are already detectable in most patients with rheumatoid arthritis halfway through an infusion cycle: an open-label pharmacokinetic cohort study. *BMC Musculoskeletal Disorders* 2011.

Vetterlein & al. No antibodies to PEG detected in patients treated with certolizumab pegol. *Ann Rheum Dis* 2008;67(Suppl II):236.

Vincent FB & al. Effect of serum anti-tumour necrosis factor (TNF) drug trough concentrations and antidrug antibodies (ADAb) to further anti-TNF short-term effectiveness after switching in rheumatoid arthritis and axial spondyloarthritis. *Joint Bone Spine*. 2016 Oct;83(5):595-7.

Wolbink GJ & al. Development of Anti-infliximab Antibodies and Relationship to Clinical Response in Patients With Rheumatoid Arthritis. *ARTHRITIS & RHEUMATISM* Vol. 54, No. 3, March 2006, pp711–715.

World Health Organization, WHO Expert Committee on Biological Standardization, Metcalfe, Clive, Bird, Chris, Dougall, Thomas. & al. (2017). Report on a collaborative study for proposed 1st international standard for Infliximab. WHO/BS/2017.2323.

Zufferey P & al. High level of anti-drug antibodies after intra-articular injection of anti-TNF. *Rheumatology (Oxford)*. 2015 Dec;54(12):2291-2

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