### **thermo**scientific

## Technical Data Sheet High Performance General Purpose Refrigerator

Revision- 2

Thermo Fisher Scientific, Asheville, North Carolina

|   | October America   |
|---|---|
|   | Catalog number  |
| Charifications  | TSX5005SV   |
| Specifications  | Model Number  |
|   | TSX5005SV Application, Rating and Electrical Data   |
| Applianting   |   |
| Application   | Storage of General (non-flammable) Laboratory Materials   |
| Storage Volume  | 1447 Liters / 51.1 Cubic Feet   |
| Factory Set Point   | 5℃  |
| Set Point Range / Adjustable Increment  | 3℃ to 7℃ / 0.1℃   |
| Electrical Power / Rated Current  | 208-230 V 50 Hz 1 Phase / 4.6 A   |
| Building Supply Rating  | 15 A dedicated grounded non-GFCI circuit. Protected by circuit breaker rated for inductive loads  |
| Power Plug / Power Cord Length  | 16 Amp (CEE 7/7) / 3.048 Meters (10 feet)   |
| Agency Listings   | CE, UL, cUL   |
| Indoor / Outdoor Usage  | Indoor Use Only   |
| Application Environment   | Non-Corrosive, Non-Flammable, Non-Explosive, Good Air Ventilation   |
| Ambient Operating Temperature   | 15°C - 32°C (60°F - 90°F)   |
|   | Refrigeration Configuration   |
| Refrigeration System  | Vapor Compression System  |
| Compressor / Number / Capacity  | Hermetic Variable Capacity Compressor / 1 / 398 W   |
| Condenser Type / Number   | Enhanced Finned-Tube and Forced-Air Cooled / 1  |
| Expansion Device  | Capillary Tube  |
| Evaporator Type / Defrost Method  | Forced Air Evaporator Coil / Automatic Defrost  |
| Refrigerant Charge / Flammability   | R290, 145 gms, GWP:3 ODP:0 / Flammable  |
| January Company   | Controller/Electrical System Configuration and Features   |
| Controller Level  | Eye Level   |
| Power Switch  | Keyed Off-On - Alarm  |
| Controller Type   | IntrLogic Control System: Digital Display, Graphic Thermometer, Automatic Defrost System  |
| Setpoint Security / Programmable  | Yes / Yes   |
| Compressor Safe Guard   | Current protection  |
| Control Sensor  | Single RTD (1000 ohm shielded Platinum RTD)   |
| Connectivity / Remote Outputs   | Remote Alarm Contacts, RS485 Output and 4-20 mA Output (optional)   |
| Adjustable Warm / Cold Alarms   | Fully Adjustable  |
| Power Failure Alarm   | Standard  |
| Electronic Chart Recorder   | Optional  |
| Electronic Chart Necorder   | Dimensions and Construction   |
| Interior Dimensions (D x W x H)   | 72.4 cm x 133.4 cm x 147.3 cm (28.5" x 52.5" x 58.0")   |
| Exterior Dimensions (D x W x H)   | 96.2 cm x 143.5 cm x 199.4 cm (37.9" x 56.5" x 78.5")   |
| Insulation  | 5.08 cm (2.0") High-Density Water-Blown Polyurethane Foam   |
| Door Perimeter Heater   |   |
|   | Yes (8) Adjustable Wire Shelves / 45 kg (100 lbs) per Shelf   |
| Shelves / Load Limit  | Standard With Locks   |
| All-Direction Casters   |   |
| Shipping Weight   | 282 kg / 621 lb   |
|   | Typical Performance Characteristics   |
| 50 cuft Upright Refrigerator, Pull Down and Warm U ———————————————————————————————————— | 50 cuft Upright Refrigerator at 5C Cycle  MAX — MIN — AVG  Test Unit Series Number or MSO Number: 19585-B-26  |
|   |   |
| 19  | Cabinet Load: Unloaded  Average Cabinet Temp at 5C Cycle (C): 4,8   |
| 17  | Peak Variation from Setpoint (C): +1.47-2.8   |
| ĝ 15  |   |
| 11 13 14 15 15 16 16 17 17 18 18 18 18 18 18 18 18 18 18 18 18 18                       | Stability (C):  1 Stability (C):  2 Duty Cycle at 5C (%):  2 Cycle (on/off) rate at 5C (min):  3 1 Cycle (on/off) rate at 5C (min):  3 2 Cycle (on/off) rate at 5C (min):  3 3 6 Cycle (on/off) rate at 5C (min):  3 3 6 Cycle (on/off) rate at 5C (min):  3 3 6 Cycle (on/off) rate at 5C (min):  3 4 Cycle (on/off) rate at 5C (min):  3 5 7 Cycle (on/off) rate at 5C (min):  3 6 7 Cycle (on/off) rate at 5C (min):  3 7 6 Cycle (on/off) rate at 5C (min):  3 8 7 Cycle (on/off) rate at 5C (min):  3 9 7 Cycle (on/off) rate at 5C (min):  3 1 Cycle (on/off) rate at 5C (min):  4 8 7 Cycle (on/off) rate at 5C (min):  5 9 7 Cycle (on/off) rate at 5C (min): |
| <u>a</u> 11   | Duty Cycle at 5C (%): 32  |
| و ق   | Cycle (on/off) rate at 5C (min): 3 / 6 Energy Consumption (kw-hr/day): 6.2  |
| 7   | Heat Rejection Rate (htulin): 993   |
|   | Pull Down Time (to 5C) (min): 17  |
| 0 50 100 150 200 250  | 0 Warm Up Time (5C to 15C) (min): 258<br>300 0 120 240 360 480 600 720 Sound (dBA): 52  |
| Time, minutes   | 300 0 120 240 360 480 600 720 Sound (dBA): 52  Time, minutes  |
| Performance is nominal and individual.  | •   |

- 1) Performance is nominal and individual units may vary.
- 2) Product performance will differ due to product amount, product size and operating conditions.
- 3) Continuous product enhancements may, without notice, result in amendments or ommisions to this specification. Thermo Scientific cannot accept responsibility for damage, injury, loss or expenses resulting from misapplication of the information herein.

## **thermo**scientific

Typical Cabinet Temperature Map Refrigerator 50cuft , 8 Inner-Shelves + Base, 2 Outer Doors Temperatures are averages during > 20 cycles after reaching a setpoint of 5C


#### **Top View of Shelves**

#### **REAR**

#### FRONT

| М3  | M5  | M4  | M53 | M54<br>M55 |  |  |
|-----|-----|-----|-----|------------|--|--|
| M1  |     | M2  | M51 | M52        |  |  |
|     | M10 |     |     | M60        |  |  |
| M13 |     | M14 | M63 | M64        |  |  |
|     | M15 |     |     | M65        |  |  |
| M11 |     | M12 | M61 | M62        |  |  |
|     | M20 |     |     | M70        |  |  |
| M23 |     | M24 | M73 | M74        |  |  |
|     | M25 |     |     | M75        |  |  |
| M21 |     | M22 | M71 | M72        |  |  |





19585-B-26

**Cabinet Average:** 4.8 C **Probe Average:** 5.4 C

Peak Variation: +1.4 C / -2.8 C

|     | M1  | M2  | M3  | M4  | M5  | M10 |
|-----|-----|-----|-----|-----|-----|-----|
| Avg | 5.4 | 4.7 | 4.6 | 4.6 | 4.7 | 4.7 |
| Max | 6.4 | 6.1 | 6   | 6   | 6   | 6   |
| Min | 4.2 | 3   | 2.8 | 2.6 | 3.1 | 2.9 |

|     | M11 | M12 | M13 | M14 | M15 | M20 |
|-----|-----|-----|-----|-----|-----|-----|
| Avg | 5   | 4.6 | 4.5 | 4.5 | 4.6 | 4.6 |
| Max | 6.2 | 6   | 6.1 | 6.1 | 6   | 6   |
| Min | 3.3 | 2.8 | 2.2 | 2.2 | 2.8 | 2.6 |


| _ |     | M21 | M22 | M23 | M24 | M25 |
|---|-----|-----|-----|-----|-----|-----|
|   | Avg | 4.6 | 4.7 | 4.7 | 4.6 | 4.5 |
|   | Max | 6.1 | 6.2 | 6   | 6.1 | 6.1 |
|   | Min | 2.5 | 2.8 | 3.1 | 2.7 | 2.3 |

|     | M51 | M52 | M53 | M54 | M55 | M60 |
|-----|-----|-----|-----|-----|-----|-----|
| Avg | 5.2 | 5.4 | 4.9 | 4.8 | 5.1 | 4.9 |
| Max | 6.3 | 6.4 | 6.2 | 6.2 | 6.3 | 6.2 |
| Min | 3.8 | 4.1 | 3.1 | 3   | 3.6 | 3.3 |

|     | M61 | M62 | M63 | M64 | M65 | M70 |
|-----|-----|-----|-----|-----|-----|-----|
| Avg | 4.9 | 5.1 | 4.7 | 4.7 | 4.9 | 4.8 |
| Max | 6.2 | 6.3 | 6.2 | 6.2 | 6.2 | 6.1 |
| Min | 3.3 | 3.6 | 2.3 | 2.6 | 3.3 | 3   |

|     | M71 | M72 | M73 | M74 | M75 |
|-----|-----|-----|-----|-----|-----|
| Avg | 4.7 | 4.8 | 4.7 | 4.7 | 4.8 |
| Max | 6.2 | 6.2 | 6.2 | 6.1 | 6.2 |
| Min | 2.6 | 2.9 | 2.4 | 2.8 | 2.7 |

# ONLY CABINET SHOWN FOR CLARITY IN TOP VIEW, TO ILLUSTRATE INTERIOR DIMENSIONS thermoscientific



52.50 [1333.5]

I.D. OF INTERIOR