

## CUTANA™ pAG-Tn5 for CUT&Tag

<b>Catalog No</b>	15-1017	<b>Species</b>	<i>E. coli</i>
<b>Lot No</b>	23082002-C1	<b>Source</b>	<i>E. coli</i>
<b>Pack Size</b>	50 Reactions	<b>Epitope Tag</b>	None
<b>Concentration</b>	20X	<b>MW</b>	191 kDa

### DESCRIPTION

CUTANA pAG-Tn5 is the key reagent for performing efficient mapping of chromatin features in Cleavage Under Targets and Tagmentation (CUT&Tag) [1]. CUT&Tag offers significant improvements in signal to noise at reduced cell inputs and sequencing depth compared to ChIP-seq. As a fusion of Proteins A and G to a highly active *E. coli* transposase mutant (Tn5), CUTANA pAG-Tn5 is compatible with target antibodies from a broad spectrum of host species. It also comes without an epitope tag, making it suitable for tag-mediated CUT&Tag (e.g. FLAG, HA, TY1, V5, etc.). This product is highly purified to remove contaminating *E. coli* DNA, which can be tagmented during storage and complicate analysis at low cell numbers. For superior normalization as well as antibody validation and reaction monitoring, SNAP-CUTANA™ nucleosome spike-ins (e.g. EpiCypher 19-1002) are recommended. The Tn5 comes charged with mosaic adapters and is ready to be used in CUT&Tag.

### RECOMMENDED ACCESSORY REAGENTS

<u>Item</u>	<u>CAT</u>	<u>Item</u>	<u>CAT</u>
Anti-Rabbit Secondary Antibody	13-0047	H3K4me3 Positive Ctrl Antibody	13-0041
Anti-Mouse Secondary Antibody	13-0048	Rabbit IgG Negative Ctrl Antibody	13-0042
Magnetic Separation Rack, 0.2 mL	10-0008	CUTANA™ ConA Beads	21-1401
Magnetic Separation Rack, 1.5 mL	10-0012	SNAP-CUTANA™ K-MetStat	19-1002
CUT&RUN 8-strip Tubes	10-0009	Non-HS 2X PCR Master Mix	15-1018

### TECHNICAL INFORMATION

<b>Storage</b>	Stable for one year at -20°C from date of receipt. The protein is not subject to freeze/thaw under these conditions.
<b>Formulation</b>	50 mM HEPES-KOH pH 7.2, 100 mM NaCl, 0.1 mM EDTA, 1 mM DTT, 0.1% Triton X-100, 50% glycerol
<b>Adapters</b>	Tn5ME-A: 5'-TCGTCGGCAGCGTCAGATGTGTATAAGAGACAG-3' Tn5ME-B: 5'-GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAG-3' Tn5ME-rev: 5'-[phos]CTGTCTCTTATACACATCT-3'

### APPLICATION NOTES

Add 2.5 µL of the supplied enzyme to a 50 µL CUT&Tag reaction (20X dilution). For detailed applications and uses of this product, please see [epicypher.com/protocols](https://epicypher.com/protocols).


### REFERENCES

[1] Kaya-Okur et al. Nat. Commun. (2019). PMID: 31036827


## VALIDATION DATA

### CUT&Tag Methods


CUT&Tag was performed on 100k K562 cells with 0.5 µg of either IgG (EpiCypher 13-0042), H3K4me3 (EpiCypher 13-0041), or H3K27me3 (EpiCypher 13-0055) antibodies using CUTANA™ pAG-Tn5 (1:20 dilution) following the EpiCypher Direct-to-PCR CUT&Tag protocol (epicypher.com/protocols). Libraries were run on an Illumina NextSeq2000 with paired-end sequencing (2x50 bp). Sample sequencing depth was 0.73 million reads (IgG), 7.6 million reads (H3K4me3), and 5.0 million reads (H3K27me3). Data were aligned to the hg19 genome using Bowtie2. Data were filtered to remove duplicates, multi-aligned reads, and blacklist regions.



**FIGURE 1: Protein gel data.** CUTANA pAG-Tn5 for CUT&Tag (1 µg) was resolved via SDS-PAGE and stained with Coomassie blue. The migration and molecular weight of the protein standards are indicated. Uncharged pAG-Tn5 monomer is 78.5 kDa, however once charged with DNA, Tn5 dimerizes to a final complex weight of 191 kDa.



**FIGURE 2: Size distribution of released chromatin.** CUT&Tag was performed as described above. Recovered DNA was directly PCR amplified to produce sequence-ready libraries. Agilent TapeStation traces for libraries derived from negative control IgG (top) and H3K27me3 (bottom) antibodies are shown. Excised DNA is highly enriched for mononucleosomes (peak at ~300 bp reflects ~150 bp insert size).



**FIGURE 3: CUT&Tag data.** CUT&Tag was performed as described above. Representative sequencing tracks obtained using CUTANA pAG-Tn5 show a 229 kb close up view of the LAMC3 gene. CUTANA pAG-Tn5 produced clear peaks with genomic distribution profiles consistent with the known biological functions of H3K4me3 and H3K27me3 as well as minimal background in the IgG negative control.

US Pat. No. 10,689,643, 11,306,307, EU Pat. No. 3,688,157, 2,999,784  
and related patents and pending applications