VarioLux Examination Light

DE	Gebrauchsanweisung Untersuchungsleuchte, Seite 2	BG	Ръководство за работа Лампа за прегледи, стр. 110
en/	Instructions for Use	HR	Upute za rad Svjetlo za pregled, stranica 118
enUS	Examination Light, page 10	DA	Brugervejledning Undersøgelseslampe, side 125
FR	Notice d'utilisation Lampe d'examen, page 17	FI	Käyttöohjeet Tutkimusvalaisin, sivu 132
ES	Instrucciones de uso Lámpara de reconocimiento, página 25	EL	Οδηγίες χρήσης Εξεταστικός φωτισμός, σελίδα 139
IT	Istruzioni per l'uso Lampada da esame, pagina 33	JA	取扱説明書 処置用ライト, 147 ページ
ptBR	Instruções de Uso Foco clínico, página 41	KO	사용 지침서 진료용 조명등, 154 페이지
NL	Gebruiksaanwijzing Onderzoekslamp, pagina 49	RO	Instrucțiuni de utilizare Lampă de examinare, pagina 161
NO	Bruksanvisning Undersøkelseslampe, side 57	SR	Uputstvo za korišćenje Svetlo za pregled, stranica 168
SV	Bruksanvisning Undersökningslampan, sida 64	SK	Návod na použitie Vyšetrovacie svetlo, strana 175
RU	Руководство по эксплуатации Медицинский светильник, стр. 71	SL	Navodilo za uporabo Preiskovalna svetilka, stran 182
PL	Instrukcja obsługi Oświetlenie zabiegowe, strona 79	TR	Kullanma Kılavuzu Muayene lambası, Sayfa 189
CS	Návod k použití Vyšetřovací světlo, strana 87	LT	Apžiūros šviestuvus naudojimo instrukcija, 196 psl.
HU	Használati útmutató Vizsgálólámpa, 95. oldal		
ZH	使用说明书 检查灯, 第 103 页		

WARNING

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VarioLux Examination Light

Trademarks

- DrägerService®
 - VarioLux™
- are trademarks owned by Dräger.

Safety information definitions

WARNING
A **WARNING** statement provides important information about a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION
A **CAUTION** statement provides important information about a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or in damage to the medical device or other property.

NOTE
A **NOTE** provides additional information intended to avoid inconvenience during operation.

Definition of target groups

For this product, users, service personnel, and experts are defined as target groups.

These target groups must have received instruction in the use of the product and must have the necessary training and knowledge to use, install, reprocess, maintain, or repair the product.

The product must be used, installed, reprocessed, maintained, or repaired exclusively by defined target groups.

Users

Users are persons who use the product in accordance with its intended use.

Service personnel

Service personnel are persons who are responsible for the maintenance of the product.

Service personnel must be trained in the maintenance of medical devices and install, reprocess, and maintain the product.

Experts

Experts are persons who perform repair or complex maintenance work on the product.

Experts must have the necessary knowledge and experience with complex maintenance work on the product.

Symbols and Abbreviations

- | | |
|---|--|
|  Atmospheric pressure |  Relative humidity |
|  Do not use if package is damaged |  Consult instructions for use |

- | | |
|--|--|
|  WEEE label, Directive 2002/96/EC |  Storage temperature limitation |
|  Caution |  Protection class II |
|  Device on |  Device off |

For USA: Rx only Caution: Federal law restricts this device to sale by or on the order of a physician.

For your safety and that of your patients

WARNING
To properly use this medical device, the user must obtain a full understanding of the performance characteristics of this medical device prior to use by carefully reading these Instructions for Use.

WARNING
Only the accessories indicated on the order list have been tested and approved for use with the medical device. Therefore, it is strongly recommended that only these accessories are used in conjunction with the medical device. Otherwise, the correct functioning of the medical device may be compromised.

WARNING
Risk of electric shock
Penetrating liquid may cause malfunction of or damage to the device, which may endanger the patient and the user.

WARNING
Impermissible modifications to the medical device lead to malfunction.
This medical device must not be modified without permission from its manufacturer.

WARNING
Use only in rooms that are dry and not subject to explosion hazards.

WARNING
Risk of personal injury
If the light fields of several lights are focussed on a single point, total thermal radiation may exceed 1000 W/m².

CAUTION
Risk of crushing
Movable parts may cause crushing injuries.

General Safety Information

Connection to other devices

Device combinations approved by Dräger (see Instructions for Use of the individual devices) meet the requirements of the following standards:

- IEC 60601-1 (EN 60601-1, 2nd Edition/3rd Edition)
Medical electrical equipment
Part 1: General requirements for safety
- IEC 60601-1-2 (EN 60601-1-2, 2nd Edition/3rd Edition)
Medical electrical equipment
Part 2: General requirements for safety
Collateral standard: Electromagnetic compatibility; Requirements and tests

If Dräger devices are connected to other Dräger devices or third-party devices and the resulting combination is not approved by Dräger, the correct functioning of the devices may be compromised. The owner is responsible for ensuring that the resulting system meets the requirements of the applicable standards.

Strictly observe Assembly Instructions and Instructions for Use for each networked device.

Information on electromagnetic compatibility

General information on electromagnetic compatibility (EMC) according to international EMC standard IEC 60601-1-2:

Medical electrical equipment is subject to special precautionary measures concerning electromagnetic compatibility (EMC) and must be installed and put into operation in accordance with the EMC information provided on page 14.

Portable and mobile RF communications equipment can affect medical electrical equipment.

Training

Training for users is available from the Dräger organization responsible, see www.draeger.com.

Product-specific safety information

WARNING

Risk of electric shock

Do not use the light and the power cable if damaged. The medical device must be checked regularly for damage.

WARNING

Risk of damage to the eyes

Avoid direct illumination of the eyes.

WARNING

Risk of damage to the device due to overheating

- The distance to combustible objects must be at least 0.5 m (19.7 in).
 - Do not store or use the light below a source of heat.
- Do not cover the light head partially or completely when using the light.**

CAUTION

Operating the medical device not in compliance with the intended ambient conditions may lead to damage to the medical device.


NOTE

Do not place any items on the light head or the hinged arms.

Intended use

Medical examination light for local illumination of the human body, to support diagnosis or treatment. The examination light is intended for continuous operation, but not for use in operating rooms.

Overview

- | | |
|--------------------------|---------------------------------|
| A Adjusting knob | B Fastening screw |
| C Rail clamp | D Power cable connection |
| E Power cable | F On/Off switch |
| G Adjusting screw | H Control panel |

Scope of delivery

- VarioLux examination light
- Rail clamp (C)
- Country-specific power cable (E)


Assembly and preparation

- 1 Fasten the rail clamp (C) to the rail.
- 2 Insert the examination light into the rail clamp (C) and fix it in position with the fastening screw (B).
- 3 Check that the mains power supply matches the nominal voltage indicated on the rating plate.
- 4 Connect the power cable (E) to the power cable connection (D) of the examination light.
- 5 Connect the mains plug to a mains power socket.

CAUTION

Risk of personal injury and of damage to the device
Do not pull the arm system beyond the end stops when positioning the examination light.

Getting started

NOTE
Secure the hinged arms when loosening the adjusting knob (A).

- 1 To bring the examination light into the required position, loosen the adjusting knob (A) by rotating it counterclockwise.
- 2 Position the examination light as required. Tighten the adjusting knob (A) by rotating it clockwise.
- 3 To switch the examination light on or off, use the On/Off switch (F): "I" = On, "O" = Off.
- After switching on the light using the On/Off switch (F) an automatic selftest will be performed and some signal lights (N) will flash for a short time.
- After the selftest the indicator light (J) lights up green.

Operation


- 1 Switch on the examination light by pressing the key (M) of the control panel (H) at the light head.
- 2 Adjust the illumination intensity by pressing the keys (K) (adjustment can be performed in steps of 20 % between 6500 lx and 30000 lx).
- The selected illumination intensity is indicated by the signal lights (N).
- 3 Press the keys (L) to select the required color according to the following table:

Color	Color temperature	Recommendation
Blue (cool white)	4700 K	Suitable to detect single vessels, for aftertreatment of wounds, and for minor operations.
White (neutral white)	4100 K	Suitable for general examination.
Red (warm white)	3500 K	Suitable for examination of the skin.

- The selected color is indicated by the signal lights (N).

Position for transport

Recommended position for intrahospital transport:

Fault – Cause – Remedy

Fault	Cause	Remedy
Indicator light (J) does not light up.	On/Off switch (F) is switched off.	Switch on the On/Off switch (F).
	Power cable (E) is not connected to power cable connection (D).	Connect power cable (E) to power cable connection (D).
	Power cable (E) is not connected to mains power supply.	Connect power cable (E) to mains power supply.
	Mains power is not available.	Have mains power checked by service personnel.
Examination light does not light up.	Examination light is faulty.	Call DrägerService.
	Key (M) at light head is switched off.	Switch on key (M) at light head.
Joint is loose.	Adjusting knob (A) cannot be fastened.	Call DrägerService.
Joint at light head is loose.	Adjusting screw (G) is loose.	Call DrägerService.

Cleaning and disinfection

Observe the hygiene regulations of the hospital.

CAUTION
Use only the detergents and disinfectants recommended.

Manual cleaning

Switch off the examination light by pressing the On/Off switch (F) and unplug the power cable (E) from the power cable connection (D).

CAUTION
Risk of burns
Before cleaning, switch off the examination light and let it cool down to room temperature.

NOTE
Keep glass panel clean by cleaning it regularly. Soiling reduces the illumination intensity.

Use glass cleaner to clean the glass panel.
Clean light components with a wipe soaked within household detergent.

Manual disinfection

Disinfectants

The cleaning and disinfection of medical devices has been tested with the following procedures and agents. The following agents showed good material compatibility at the time of the test:

- Surface disinfectants (for device surfaces)
- Dismozon® pur by Bode Chemie GmbH & Co
- Buraton® 10 F by Schülke & Mayr GmbH

Manual disinfection should preferably be carried out with disinfectants based on aldehydes or quaternary ammonia compounds.

Observe the applicable country-specific listings for disinfectants. The list of the German Association for Applied Hygiene (Verbund für Angewandte Hygiene VAH) applies in German-speaking countries.

The composition of the disinfectant is the responsibility of the manufacturer and can change over time.

Strictly observe the manufacturer's information on the disinfectant.

Disinfecting surfaces

- Remove soilings with a wet wipe.

WARNING

Penetrating liquid may cause malfunction of or damage to the device, which may endanger the patient and the user. Only wipe-disinfect device surface and cables and make sure no liquids penetrate into the device.

- 1 Carry out surface disinfection (scrub/wipe disinfection).
- 2 After contact time, remove disinfectant residues.

Checking proper functioning after cleaning and disinfecting

- 1 Connect power cable (E) to power cable connection (D) and switch on the examination light by pressing the On/Off switch (F).
- 2 Check proper functioning according to chapters "Getting started" and "Operation".

Maintenance

Definition of maintenance concepts

Concept	Definition
Maintenance	All measures (inspection, preventive maintenance, repair) intended to maintain and restore the functional condition of a medical device
Inspection	Measures intended to determine and assess the actual state of a medical device
Preventive maintenance	Recurrent specified measures intended to maintain the functional condition of a medical device
Repair	Measures intended to restore the functional condition of a medical device after a device malfunction

Inspection

Perform inspections at regular intervals and observe the following specifications.

Checks	Interval	Personnel responsible
Inspection	Every 2 years	Experts

- 1 Check accompanying documents:
 - Instructions for use are available
- 2 Perform a functional test of the following features according to the instructions for use:

- See "Operation"

- 3 Check that the device is in good condition:

- All labels are complete and legible
- There is no visible damage

- 4 Using the instructions for use, check that all components and accessories needed to use the product are available.

- 5 Check the electrical safety according to IEC 62353.

Repair

Dräger recommends that all repairs are carried out by DrägerService and that only authentic Dräger repair parts are used.

Disposal of the Medical Device

When disposing of the medical device:

- Consult the relevant waste disposal company for appropriate disposal.
- Observe the applicable laws and regulations.

For countries subject to the EU Directive 2002/96/EC

This device is subject to EU Directive 2002/96/EC (WEEE).

It is not registered for use in private households.


In order to comply with its registration according to this directive, this device may not be disposed of at municipal collection points for waste electrical and electronic equipment. Dräger has authorized a company to collect and dispose of this device. To initiate collection or for further information, visit Dräger on the Internet at www.draeger.com. Use the Search function with the keyword "WEEE" to find the relevant information. If access to Dräger's website is not possible, contact the local Dräger Organization.

Technical Data

Ambient conditions

During operation	
Temperature	10 to 40 °C (50 to 104 °F)
Atmospheric pressure	700 to 1060 hPa (10.2 to 15.4 psi)
Relative humidity	10 to 95 %, non-condensing
During storage and transport	
Temperature	–20 to 60 °C (–4 to 140 °F)
Atmospheric pressure	500 to 1060 hPa (7.3 to 15.4 psi)
Relative humidity	10 to 95 %, non-condensing
Operating mode	Continuous operation

Dimensions (mm)

Weight 1.5 kg (3.3 lbs)

Electrical values

Nominal voltage	100 – 240 V ~
Frequency range	50/60 Hz
Effective power (P)	18 W
Apparent power (S)	26 to 34 VA
Classification	Protection class II

Photometric values

Maximum illumination intensity (Ec) at 0.5 m (19.7 in)	30000 lx
Dimming levels	30000 lx, 24000 lx, 18000 lx, 12000 lx, 6500 lx
Light field diameter d10 at 0.5 m (19.7 in)	Ø = 21 cm (8.3 in)
Light field diameter d50 at 0.5 m (19.7 in)	Ø = 9 cm (3.5 in)
Color temperature	Adjustable
Cool white	4700 K ±2 %
Neutral white	4100 K ±2 %
Warm white	3500 K ±2 %

Color reproduction index Ra	≥ 96 (at 3500 K)
Color reproduction index R9	90 (at 3500 K)
Minimum service life of bulb	30000 h
Maximum radiant power at 0.5 m (19.7 in)	105 W/m²
Degree of protection	IP 20
Electrical safety testing and electromagnetic compatibility according to	IEC 60601-1 IEC 60601-2-41 IEC 60601-1-2

Classification according to Directive 93/42/EEC Annex IX Class 1

UMDNS Code 12-347

Universal Medical Device Nomenclature System – Nomenclature for medical devices

GMDN Code 12276

Global Medical Device Nomenclature – Nomenclature for medical devices

EMC Declaration

General information

The EMC compliance of the medical device is also applicable to the external cables, transducers, and accessories specified in the list of accessories. In addition, accessories which do not affect EMC compliance may be used if no other reasons forbid their use (see other sections of the Instructions for Use). The use of non-compliant accessories may result in increased emissions or decreased immunity of the medical device.

Electromagnetic environment

The medical device is intended for use in an electromagnetic environment as specified below. The user must ensure its use in such an environment.

Emissions	Compliance according to	Electromagnetic environment
Radio frequency emissions (CISPR 11)	Group 1	The medical device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
	Class B	The medical device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions (IEC 61000-3-2)	Class A	
Voltage fluctuations/flicker emissions (IEC 61000-3-3)	Complies	

Electromagnetic immunity

The medical device is intended for use in an electromagnetic environment as specified below. The user must ensure its use in such an environment.

Immunity against	Test level (IEC 60601-1-2)	Compliance level (medical device)	Electromagnetic environment
Electrostatic discharge (ESD) (IEC 61000-4-2)	Contact discharge: ± 6 kV	± 6 kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
	Air discharge: ± 8 kV	± 8 kV	
Electrical fast transients/ bursts (IEC 61000-4-4)	Power supply lines: ± 2 kV	± 2 kV	Mains power quality should be that of a typical commercial or hospital environment.
	Longer input / output lines: ± 1 kV	Not applicable	
Surge on AC mains lines (IEC 61000-4-5)	Common mode: ± 2 kV	± 2 kV	Mains power quality should be that of a typical commercial or hospital environment.
	Differential mode: ± 1 kV	± 1 kV	
Power frequency magnetic field (50/60 Hz) (IEC 61000-4-8)	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Voltage dips and short interruptions on AC mains input lines (IEC 61000-4-11)	Dip >95 %, 0.5 periods	Dip >95 %, 0.5 periods	Mains power quality should be that of a typical commercial or hospital environment. If the user of the medical device requires continued operation during mains power supply interruptions, it is recommended that the medical device be powered from an uninterruptible power supply or a battery.
	Dip 60 %, 5 periods	Dip 60 %, 5 periods	
	Dip 30 %, 25 periods	Dip 30 %, 25 periods	
	Dip >95 %, 5 seconds	Dip >95 %, 5 seconds	
Radiated RF (IEC 61000-4-3)	80 MHz to 2.5 GHz: 10 V/m	10 V/m	Recommended minimum distance to portable and mobile RF transmitters with transmission power PEIRP to the medical device including its lines $(1.84 \text{ m} \times \sqrt{\text{PEIRP}})^{1)}$
	80 MHz to 2.5 GHz: 3 V/m	3 V/m	

Immunity against	Test level (IEC 60601-1-2)	Compliance level (medical device)	Electromagnetic environment
Conducted RF (IEC 61000-4-6)	150 kHz to 80 MHz: 3 V	3 V	Recommended minimum distance to portable and mobile RF transmitters with transmission power P_{EIRP} to the medical device including its lines $(1.84 \text{ m} \times \sqrt{P_{EIRP}})^{1)}$

1) For P_{EIRP} the highest possible "equivalent isotropic radiated power" of the adjacent RF transmitter has to be inserted (value in Watts). In the vicinity of equipment marked with the symbol , interference may also occur. Field strengths from fixed, portable or mobile RF transmitters at the location of the medical device should be less than 3 V/m in the frequency range from 150 kHz to 2.5 GHz and less than 1 V/m above 2.5 GHz.

Recommended separation distances to portable and mobile RF telecommunication devices

The separation distances below are in accordance with IEC 60601-1-2.

Max. P_{EIRP} (W)	150 kHz to 2.5 GHz	All other frequencies	Examples
0.03	0.32 m (1.05 ft)	0.96 m (3.15 ft)	e.g., WLAN 5250 / 5775 (Europe)
0.10	0.58 m (1.90 ft)	1.75 m (5.74 ft)	e.g., WLAN 2440 (Europe)
0.17	0.76 m (2.49 ft)	2.28 m (7.48 ft)	e.g., Bluetooth, RFID 2.5 GHz
0.20	0.82 m (2.69 ft)	2.47 m (8.10 ft)	e.g., WLAN 5250 (not in Europe)
0.25	0.92 m (3.02 ft)	2.76 m (9.06 ft)	e.g., UMTS mobiles
0.41	1.18 m (3.87 ft)	3.53 m (11.58 ft)	e.g., DECT devices
0.82	1.67 m (5.48 ft)	5.00 m (16.40 ft)	e.g., RFID 13.56 MHz
1.00	1.84 m (6.04 ft)	5.52 m (18.11 ft)	e.g., WLAN 5600 (not in Europe)
1.64	2.36 m (7.74 ft)	7.07 m (23.20 ft)	e.g., GSM 1800 / GSM 1900
3.28	3.33 m (10.93 ft)	10.00 m (32.81 ft)	e.g., GSM 900 mobiles, RFID 868 MHz

Order List

Description	Part No.
VarioLux EU	MP00601
VarioLux US	MP00602
VarioLux GB	MP00603
VarioLux CHN	MP00604
VarioLux BR	MP00605
VarioLux AUS	MP00606
VarioLux ZA	MP00607
VarioLux JPN	MP00608
Accessories	
Rail clamp for light	MP00615
Adjusting knob, complete	MP00587
Control panel	MP00589
Power cable EU, 1.8 m (70.9 in), 2-pin	MP00591
Power cable US, 1.8 m (70.9 in), 2-pin	MP00592
Power cable GB, 1.8 m (70.9 in), 2-pin	MP00593
Power cable CHN, 1.8 m (70.9 in), 2-pin	MP00594

Description	Part No.
Power cable ZA, 1.8 m (70.9 in), 2-pin	MP00597
Power cable BR, 1.8 m (70.9 in), 2-pin	MP00625
Power cable AUS, 1.8 m (70.9 in), 2-pin	MP00626

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© Dräger Medical GmbH