Hinged Arms

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Instructions for Use of hinged arms

CAUTION

The Instructions for Use of the basic machine must be strictly observed!

WARNING

For a full understanding of the performance characteristics of this medical product, the user should carefully read this manual before use of the medical product.

Definitions

WARNING

A WARNING statement provides important information about a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION

A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or in damage to the equipment or other property.

NOTE

A NOTE provides additional information intended to avoid inconvenience during operation of the medical product.

Symbols


Consult Instructions for Use!

Temperature limit during storage

Relative humidity during storage

Atmospheric pressure during storage

Intended Use

The hinged arms are used for precise positioning and safe fixation of breathing hoses and pressure transducers.

Installation and Operation

A Releasing = turn in a counter-clockwise direction, until end position is reached.
Locking = turn in a clockwise direction.

NOTE

B Traction of the device by means of hinged arm is not allowed.

Maximum load capacity

C Load capacity






Length (m) (inch)	0.2 8	0.4 16	0.6 24	0.8 32	1 40	1.2 48
Weight (kg) max. (lbs)	2 4.4	1 2.2	0.7 1.5	0.5 1.1	0.4 0.9	0.3 0.7

Cleaning and Disinfecting

- Follow the hospital hygiene regulations.
- Remove visible dirt with a disposable cloth.
- After manual cleaning, perform surface disinfection.

Manual disinfection

Manual disinfection can be carried out preferably with disinfectants based on aldehydes or quaternary ammonia compounds. The potency of the disinfectants used must be proven. Observe the applicable country-specific listings. The list of the Association for Applied Hygiene (VAH List) applies in German-speaking countries. The user must strictly observe the manufacturer's information on cleaning agents.

Disinfecting surfaces

WARNING

Penetrating liquid may lead to failure of the product or damage to the product and endanger the patient! Only disinfect parts by wiping and make sure no liquids penetrate into the device.

The following disinfectants are recommended: Buraton 10F (Schülke & Mayr, Norderstedt, Germany).

- Remove disinfectant residues.
- Hinged arms are not suitable for disinfection by immersion or for disinfection in a washer-disinfector.

Disposal of the medical device

When disposing of the medical device:

- Consult the relevant waste disposal company for appropriate disposal.
- Observe the applicable local statutory regulations.

Technical Data

Ambient conditions

During operation:

Temperature	10 to 40 °C (50 to 104 °F)
Air pressure	700 to 1060 hPa (10.2 to 15.4 psi)
Relative humidity	10 to 90 %, no condensation

During storage and transportation:

Temperature	-40 to 70 °C (-40 to 158 °F)
Air pressure	700 to 1100 hPa (10.2 to 16 psi)
Relative humidity	0 to 90 %, no condensation

Order-List

Name	Part No.
Quickstop hinged arm 1	2M85705
Quickstop hinged arm 2	2M85706
Quickstop hinged arm, short	2M86463
Quickstop hinged arm Pod	2M86464
Hinged arm 240 degrees	8409609
Accessories	
Hose Clip, coaxial (2 pieces)	MP00493
Set of Hose Clips (2 pieces)	8406138
Hose Clip (double) (for ball mounting)	8409841
Support for Y-piece (for ball mounting)	8409746
Universal tube holder plate	2M86465
Quickstop Set MultiMed Pod for 2M86464	MP00653
Hinged arm, middle part	MP00583

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© Dräger Medical GmbH
Draeger Medical, Inc.
U.S.A.