Gel-Matratze 2M 20 827

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Intended Use

The gel mattress is intended for use as an anti-bedsore thermal mattress in the Dräger Babytherm 8000 WB and Babytherm 8010.

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Safety instructions

WARNING

Do not use in the Dräger Babytherm 4200 thermal bed or in other thermal beds with a hot plate temperature in excess of 41 °C. Risk of burns.

WARNING

Do not use in the Babytherm 8004. Risk of hypothermia.

WARNING

Observe the heating period of the gel mattress. Risk of hypothermia.

WARNING

Do not switch off the device while an infant is lying on the gel mattress. Risk of hypothermia.

CAUTION

The gel mattress can be damaged! Avoid contact with any sharp objects! Do not fold the gel mattress completely flat. Always store and transport flat and unfolded. Do not store near to strong heat sources (>50 °C).

CAUTION

The gel mattress may become permanently stained if it comes into contact with certain disinfectants (e.g. iodine-containing or peroxide-based). This does not impair the function and the mattress can still be used.

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Storage / transportation:

Temperature	–20 to 60 °C (–4 to 140 °F)
Atmospheric pressure	700 to 1060 hPa (10.15 to 15.37 psi)
Rel. humidity	0 to 90 %, without condensation

Dimensions

Length x width x height	745 x 485 x 13 mm 29.3 x 19.1 x 0.5 inches
-------------------------	---

Construction

The gel mattress has a plasticiser-free silicone core with high elasticity and restoring force. The outer cover consists of polyurethane plastic and is skin-friendly.

Materials used

Membrane	POLYETHYLENE
Core	Silicone gel (mercury-free; plasticiser-free)

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Gel Mattress 2M 20 827 (only for USA)

WARNING

For a full understanding of the performance characteristics of this accessory, the user should carefully read this supplement as well as the Instructions for Use of the device to which it attaches.

Intended Use

The Gel Mattress is intended for use as an anti-bedsore heated mattress only in the Dräger Babytherm 8000 WB and Babytherm 8010.

Preparing for Use

- Lay the mattress flat in the cot and cover it with a cotton sheet.

Safety instructions

WARNING

Do not use in the Dräger Babytherm 4200 thermal bed or in other thermal beds with a hot plate temperature in excess of 41 °C (106 °F). The gel mattress cannot tolerate these high temperatures and will over heat thereby risking burns to the patient.

WARNING

Do not use in the Babytherm 8004, because the Gel Mattress will not be heated. Risk of hypothermia.

WARNING

Observe the warm-up time of the Gel Mattress (1 hour). Risk of hypothermia.

WARNING

Do not switch off the device while an infant is lying on the Gel Mattress as the mattress will cool and there is risk of hypothermia.

CAUTION

Keep sharp objects away from the mattress and do not bend the gel mattress since the mattress may be damaged and its correct functioning may be compromised. For the same reason, always store and transport flat. Do not store near to strong heat sources (>50 °C (>122 °F)).

CAUTION

Although the gel may become permanently stained if it comes into contact with liquid, alcoholic iodine-containing or iodine-like disinfectants, this will not impair the function and the mattress can still be used.

Maintenance

The gel mattress is highly resistant to a multitude of detergents and to UV light. However, the gel mattress is subject to different degrees of wear, depending on the way it is used. Frequent and prolonged contact with urine, secretions or chemicals may affect the elasticity and serviceability of the gel mattress. The gel mattress must be replaced in the case of visible material changes or cracks. It is recommended to replace the gel mattress after 2 years, at the latest.

Reprocessing/cleaning

- Remove the mattress from the bed. Store the gel mattress flat.

WARNING

To avoid risks to hospital staff and other patients, disinfect and clean the device after use. When disinfecting contaminated parts, follow the hospital hygiene regulations (protective clothing, eye protection, etc.)

Disinfect and clean the gel mattress:

- after each change of patient
- at least once a week
- as required

Wipe disinfect with commercial disinfectants (e.g. Incidur, Ecolab, Germany). Observe the correct exposure time for the specific disinfectant in order to prevent excessive disinfection. Note the disinfectant manufacturer's instructions. In case of doubt, check material compatibility (polyurethane film/silicone core) or obtain confirmation from the disinfectant manufacturer.

CAUTION

Only clean and wipe disinfect the gel mattress.

CAUTION

In order to avoid irreversible material damage, do not wash the mattress in a washing machine and do not autoclave.

Following cleaning/disinfection, completely remove disinfectant residues by washing, drying and airing. Consideration of patient compatibility with any disinfectant is the responsibility of the medical personnel and the hygiene specialist.

Information on how to extend the service life of the Gel Mattress.**NOTE**

When using surface disinfectants other than those recommended, observe the disinfectant manufacturers' recommendations. In addition to the main active agents, disinfectants often contain constituents that can cause material damage. In case of doubt, check the compatibility of the used material or obtain confirmation from the disinfectant manufacturer. The user is liable for any damage caused by the use of unsuitable products.

For reasons of material compatibility, we recommend the use of aldehyde-based surface disinfectants (e.g. Incidur® from Ecolab). Normally highly suitable are disinfectants based on:

- aldehydes

Normally unsuitable or only conditionally suitable are disinfectants based on:

- halogen-releasing compounds
- strong organic acids/bases
- O₂-releasing compounds (peroxide)
- alcohol

CAUTION

The use of unsuitable disinfectants (e.g. peroxide-based or iodine-containing) may cause irreversible material change.

NOTE

The components to be reprocessed must only be subjected to brief contact with disinfectant. Do not immerse. Select the lowest possible active agent concentration. After allowing the disinfectant time to take effect (see manufacturer's specifications), wipe the surface with a clean damp disposable cloth and then wipe dry. Note the disinfectant manufacturer's instructions.

Technical Data**Ambient conditions**

During operation:

Temperature	15 to 35 °C (59 °F to 95 °F)
Atmospheric pressure	900 to 1060 hPa (13.05 to 15.37 psi)
Rel. humidity	0 to 75 %, without condensation

Storage / transportation:

Temperature	–20 to 60 °C (–4 to 140 °F)
Atmospheric pressure	700 to 1060 hPa (10.15 to 15.37 psi)
rel. humidity	0 to 90 %, without condensation

Dimensions

Length x width x height	745 x 485 x 13 mm 29.3 x 19.1 x 0.5 inches
-------------------------	---

Construction

The gel mattress has a plasticizer-free silicone core with high elasticity and restoring force. The outer cover consists of polyurethane plastic and is skin-friendly.

Materials used

Membrane	POLYETHYLENE
Core	Silicone gel (mercury-free; plasticizer-free)

Classification

(according to EC Directive 93/42/EEC Appendix IX)

Class I

UMDNS code

(Universal Medical Device Nomenclature System)

12-473

Directive 93/42/EEC
concerning Medical Devices

 Manufacturer:

Dräger Medical AG & Co. KG

Moislinger Allee 53 – 55
D-23542 Lübeck
Germany
+49 451 8 82-0

FAX +49 451 8 82-20 80
 <http://www.draeger.com>

For USA only distributed by

Draeger Medical, Inc.

3135 Quarry Road
Telford, PA 18969
U.S.A.

(215) 721-5400
(800) 4DRAGER
(800 437-2437)

FAX (215) 723-5935
 <http://www.draeger.com>

90 39 195 - de/en/fr/es/nl/it/sv/da/fi/ru/pt/
pl/no/hu/cs/sk/tr/zh/ja/enUS

© Dräger Medical AG & Co. KG

4th edition – July 2008

Dräger reserves the right to make
modifications to the equipment without
prior notice.