Flowmeter

de	Durchflussmesser Gebrauchsanweisung, Seite 2	sv	Flödesmätare Bruksanvisning, sida 66	ro	Debitmetru Instrucțiuni de utilizare, pagina 131
en/	Flowmeter	ru	Расходомер Руководство по эксплуатации, на стр. 74	sr	Merač protoka Uputstvo za korišćenje, na stranici 139
enUS	Instructions for use, page 10				
fr	Débitmètre Notice d'utilisation, page 18	pl	Przepływomierz Instrukcja obsługi, strona 83	bg	Ротаметър Ръководство за работа, на страница 147
es	Flujómetro Instrucciones de uso, página 26	cs	Průtokoměr Návod k použití, strana 91	el	Ροόμετρο Οδηγίες χρήσης, στη σελίδα 155
it	Flussometro Istruzioni per l'uso, pagina 34	sk	Prietokomer Návod na použitie, strana 99	tr	Akış Ölçer Kullanma kılavuzu στη σελίδα 163
ptBR	Fluxômetro Instruções de uso, página 42	sl	Merilnik pretoka Navodilo za uporabo, stran 107	zh	流量计 使用说明书、第 171 页
nl	Flowmeter Gebruiksaanwijzing, pagina 50	hu	Áramlásmérő Használati útmutató, 115. oldal	ko	유량계 사용 지침서, 178 페이지
no	Flowmåler Bruksanvisning, side 58	hr	Mjerač protoka Upute za rad, na stranici 123	ja	流量計 取扱説明書、185 ページ

WARNING

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Instructions for use Flowmeter

For your safety and that of your patients	11
Intended use	12
Overview	12
Preparation	12
Functional test	13
Operation	13
Operation in magnetic fields	13
Measured value	13
Shut down	13
Fault-Cause-Remedy	14
Maintenance	14
Cleaning and disinfection	14
Manual cleaning and disinfection	15
Disposal of the medical device	15
Technical data	15
Order list	16
Replacement parts	16
Recommended accessories for operation	17

Safety information definitions

WARNING

A WARNING statement provides important information about a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION

A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or in damage to the medical device or other property.

NOTE

A NOTE provides additional information intended to avoid inconvenience during operation.

Definition of target groups

For this product, users, service personnel, and experts are defined as target groups.

These target groups must have received instruction in the use of the product and must have the necessary training and knowledge to use, install, reprocess, or maintain the medical device.

The product must be used, installed, reprocessed, or maintained exclusively by defined target groups.

Users

Users are persons who use the product in medical rooms in accordance with its intended use.

Service personnel

Service personnel are persons who are responsible for the maintenance of the product.

Service personnel must be trained in the maintenance of medical devices and install, reprocess, and maintain the product.

Experts

Experts are persons who perform repair or complex maintenance work on the product.

Experts must have the necessary knowledge and experience with complex maintenance work on the product.

Symbols

Keep free of oil and grease

Atmospheric pressure

Relative humidity

Storage temperature limitation

Manufacturer

Date of manufacture

Quantity

Part number

Serial number

When using, do not exceed the 20 mT line

For USA: Rx only

Caution: Federal law restricts this device to sale by or on the order of a physician.

For your safety and that of your patients

Strictly follow these instructions for use

WARNING

Any use of the medical device requires full understanding and strict observation of all sections of these instructions for use. The medical device must only be used for the purpose specified under "Intended use" on page 12. Strictly observe all WARNING and CAUTION statements throughout these instructions for use and all statements on medical device labels. Failure to observe these safety information statements constitutes a use of the medical device that is inconsistent with its intended use.

WARNING

Do not modify the medical device. Modification may damage or impair the proper functioning of the device which may lead to patient injury.

Maintenance

CAUTION

Risk of medical device failure and of patient injury. The medical device must be inspected regularly by service personnel. Repair and complex maintenance carried out on the medical device must be performed by experts.

If the above is not complied with, medical device failure and patient injury may occur. Observe chapter "Maintenance".

Dräger recommends that a service contract is obtained with DrägerService and that all repairs are performed by DrägerService. For maintenance Dräger recommends the use of authentic Dräger repair parts.

Accessories

CAUTION

Only the accessories indicated on the "Order list" on page 16 have been tested and approved for use with the medical device.

Therefore, it is strongly recommended that only these accessories are used in conjunction with the medical device. Otherwise, the correct functioning of the medical device may be compromised.

Safety notes for handling oxygen

WARNING

Risk of fire

Do not allow the O₂ metering valve to come into contact with oil, grease, or flammable liquids. Avoid open flames, sparks or other possible ignition sources.

Patient safety

The design of the medical device, the accompanying documentation, and the labeling on the medical device are based on the assumption that the purchase and the use of the medical device are restricted to persons familiar with the most important inherent characteristics of the medical device.

Instructions and WARNING and CAUTION statements are therefore largely limited to the specifics of the Dräger medical device.

The instructions for use do not contain any information on the following points:

- Risks that are obvious to users
- Consequences of obvious improper use of the medical device
- Potentially negative effects on patients with different underlying diseases

Medical device modification or misuse can be dangerous.

Intended use


The flowmeter is intended to control the flow of oxygen or Air for inhalation or insufflation in conjunction with the accessories needed for operation.

Only for use on patients with spontaneous breathing.

WARNING

Set the O₂ flow necessary for the patient. Monitoring of oxygen saturation is required, especially for neonates.

Overview

- A** Probe of flowmeter (DIN)
- B** Flow tube with flow display
- C** Hand wheel for flow adjustment
- D** Connector for humidifier, socket or connection hoses
- E** Terminal unit
- F** Probe of flowmeter, British Standard (BS)
- G** Probe of flowmeter, Norme Française (NF)
- H** Probe of flowmeter, Japan Industrial Standard (JIS)
- I** Rail clamp
- J** Probe of supply hose
- K** NIST connector for supply hose
- L** Operating pressure

Preparation

CAUTION

Risk of patient injury
Always check that the delivered gas type is appropriate for the application.

- Prepare accessories in accordance with the respective instructions for use.
- 1** Screw Dräger humidifier, nebulizer, or adapter to the connector (D) of the flowmeter using the union nut. Alternatively, screw on commercial sterile water system. Connect application device.
- 2** If necessary, connect the supply hose to the connector (K) on the rail clamp using the NIST connection.
- 3** Place rail clamp (I) on rail and fix in position by pressing it down.

- 4 Connect probe of flowmeter, country-specific (A), (F), (G) or (H), or probe of supply hose (J) to terminal unit (E).

Functional test

- Open hand wheel by turning it counterclockwise. The ball in the flow tube must rise evenly.
- Read flow at upper rim of ball.

Operation

WARNING

Risk of malfunction

Obstructions, damage, and foreign matter can lead to malfunction.

Check all system components for obstructions, damage, and foreign matter before installation.

WARNING


Risk of breakage

Check the medical device before use, particularly for cracks in the pressure dome. Unscrew the depressurized pressure dome and check the thread for irregularities and cracks at regular intervals (at least once a year). Do not use the medical device if it is damaged.

WARNING

Use the medical device only within the operating pressure range specified on the flow tube and within the specified environmental conditions. Observe section "Technical data". Otherwise, the set flow and the actual flow do not correspond.

Do not expose flowmeter to direct sunlight.

- Set the flow necessary for the patient on the device (A).
- Application and delivery only according to physician's instructions.
- SpO₂ monitoring is recommended.

Operation in magnetic fields

Due to the very small amount of ferromagnetic material they contain, the flowmeters can be used in magnetic fields, i.e., in the vicinity of magnetic resonance imaging systems.

When using, do not exceed the 20 mT line.

Measured value

The measured value refers to an ambient temperature of 20 °C (68 °F) and an ambient pressure of 1013 hPa.

To convert this value to the flow at ambient conditions $\dot{V}_2$, proceed as follows (gas equation):

$$\dot{V}_2 = \frac{\dot{V}_1 \times T_2 \times p_1}{p_2 \times T_1}$$

$\dot{V}_1$ = Flow displayed in L/min (at 1013 hPa and 20 °C/68 °F)

$\dot{V}_2$ = Flow rate at

p_2 = Ambient pressure (hPa)

t_2 = Ambient temperature (°C/°F)

p_1 = Reference ambient pressure 1013 hPa

t_1 = Reference temperature 20 °C (68 °F)

$T_1 = t_1 + 273 \text{ °C (523.4 °F)}$

$T_2 = t_2 + 273 \text{ °C (523.4 °F)}$

Example:

Displayed value 10 L/min, temperature 24 °C, ambient pressure 960 hPa

$$\dot{V}_2 = \frac{10 \text{ L/min} \times (24 \text{ °C} + 273 \text{ °C}) \times 1013 \text{ hPa}}{960 \text{ hPa} \times (20 \text{ °C} + 273 \text{ °C})}$$

$\dot{V}_2 = 10.7 \text{ L/min}$ (at 960 hPa and 24 °C)

Shut down

- 1 Close the hand wheel by turning it clockwise.

NOTE

Use only slight pressure when closing the hand wheel, otherwise the seal may be damaged.

- 2 Remove accessories.
- 3 Remove rail flowmeter from the rail by pressing the toggle.

Fault-Cause-Remedy

Fault	Cause	Remedy
Ball sticks after flow has been terminated	Static charge	Disconnect flowmeter and set down briefly
	Soiling inside the flow tube	Replace flow tube, call DrägerService
Cracks and/or clouding of the pressure dome	Incorrect disinfectant or exposure to sunlight	Replace the pressure dome. Call DrägerService
Leakage between flowmeter and connected accessories	Missing or worn sealing ring at the gas outlet	Call DrägerService

Maintenance

This chapter describes the maintenance measures required to maintain the proper functioning of the medical device. Maintenance measures must be performed by the personnel responsible.

WARNING

Risk of infection

Users and service personnel can become infected with pathogenic germs. Disinfect and clean device or device parts before any maintenance measures and also before returning the medical device for repair.

Definition of maintenance concepts

Concept	Definition
Maintenance	All measures (inspection, preventive maintenance, repair) intended to maintain and restore the functional condition of a medical device
Inspection	Measures intended to determine and assess the actual state of a medical device
Preventive maintenance	Recurrent specified measures intended to maintain the functional condition of a medical device
Repair	Measures intended to restore the functional condition of a medical device after a device malfunction

Inspection

Perform inspections at regular intervals and observe the following specifications.

Checks	Interval	Personnel responsible
Inspection	Every 5 years	Service personnel

Performing inspection

- 1 Check completeness and readability of labels
- 2 Visual inspection of
 - connectors
 - gas type code conformity
 - the sealing ring at the gas outlet

- the CS connection hose, if available
- the pressure dome and the respective sealing ring for cracks and contamination

3 Functional check

Repair

Dräger recommends that all repairs are carried out by DrägerService and that only authentic Dräger repair parts are used.

Cleaning and disinfection

WARNING

Risk of infection

Reusable products must be reprocessed. Otherwise there is an increased risk of infection and the functioning of the products may be impaired.

- Comply with the hospital's hygiene regulations.
- Use validated procedures for reprocessing.
- Reprocess reusable products before each use.
- Observe the manufacturer's instructions on cleaning agents and disinfectants.

WARNING

Risk due to faulty products

Reusable products may show signs of wear, such as cracks, deformations, discolorations, or peeling. Check products for signs of wear and replace if necessary.

WARNING

Risk of stress cracking
Do not use disinfectants or cleaning agents
containing alcohol.

Observe the hygiene regulations of the hospital.

Testing of procedures and agents

The cleaning and disinfection of medical devices has been tested with the following procedures and agents. The following agents showed good material compatibility at the time of the test:

- Dismozon pur® by Bode Chemie GmbH & Co
- Buraton® 10 F by Schülke & Mayr

Manual cleaning and disinfection

Manual disinfection should preferably be carried out with disinfectants based on aldehydes.

Observe the applicable country-specific listings for disinfectants. The list of the German Association for Applied Hygiene (Verbund für Angewandte Hygiene VAH) applies in German-speaking countries.

The composition of the disinfectant is the responsibility of the manufacturer and can change over time.

Strictly observe the manufacturer's information on the disinfectant.

Carry out manual cleaning and disinfection

- Remove soiling with a wipe soaked in disinfectant.

WARNING

Penetrating liquid may cause malfunction or damage to the device, which may endanger the patient. Only wipe-disinfect device surface and make sure no liquids penetrate into the device.

- 1 Carry out surface disinfection.
- 2 After contact time, remove disinfectant residues.

Disposal of the medical device

When disposing of the medical device:

- Consult the relevant waste disposal company for appropriate disposal.
- Observe the applicable laws and regulations.

Technical data**Ambient conditions**

During operation

Temperature	10 °C to 40 °C (50 °F to 104 °F)
Atmospheric pressure	700 to 1200 hPa (10.2 to 17.4 psi)
Relative humidity	0 to 95 %, no condensation

During storage and transport

Temperature	–20 °C to 70 °C (–4 °F to 158 °F)
Atmospheric pressure	800 to 1200 hPa (11.6 to 17.4 psi)
Relative humidity	0 to 90 %, no condensation

Performance characteristics

Operating pressure	4.5 bar (65.3 psi)
Setting ranges	0.5 to 4 L/min 2 to 10 L/min 2 to 16 L/min 4 to 32 L/min

Accuracy of flow rate at: ±15 % (corresponding to smallest setting value)
 4.5 bar (65.3 psi)*
 (1013 hPa/20 °C/68 °F)** ±10 % (all other setting values)

Device characteristics

Screw thread for accessories

Oxygen	M 34 x 1.5 or 9/16" adapter
Air	M 32 x 1.5

Dimensions (H x W x D) 160 x 40 x 130 mm
(6.3 x 1.6 x 5.1 in)

Material

Housing	Aluminum
Pressure dome	PA PACM 12
Flow tube	Glass

Weight

CS version	approx. 500 g (1.1 lbs)
Rail-mounted device	approx. 950 g (2.09 lbs)
Double flowmeter	approx. 1200 g (2.65 lbs)

Classification

as per EC Directive 93/42/EEC Annex IX Class IIa

UMDNS code

11-748
 Universal Medical Device
 Nomenclature System –
 Nomenclature for medical
 devices

* 1 bar = 1 kPa x 100


** If flowmeter is used outside the specified ambient conditions, see page 13.

Order list

Description	Part no.
O ₂ flowmeter 4 L/min, central supply, DIN	MP04551
O ₂ flowmeter 4 L/min, rail, NIST	MP04552
O ₂ flowmeter 4 L/min, rail, 90°, NIST	MP04553
O ₂ flowmeter 10 L/min, central supply, Carbamed	MP04554
O ₂ double flowmeter 10 L/min, central supply, Carbamed	MP04555
O ₂ flowmeter 10 L/min, rail, NIST	MP04556
O ₂ double flowmeter 10 L/min, rail, NIST	MP04557
O ₂ flowmeter 16 L/min, central supply, DIN	MP04558
AIR flowmeter 16 L/min, central supply, DIN	MP04559
O ₂ double flowmeter 16 L/min, central supply, DIN	MP04560
O ₂ flowmeter 16 L/min, central supply, DIN, long	MP04561
AIR flowmeter 16 L/min, central supply, DIN, long	MP04562
O ₂ flowmeter 16 L/min, central supply, Carbamed	MP04563
AIR flowmeter 16 L/min, central supply, Carbamed	MP04564
O ₂ flowmeter 16 L/min, central supply, British Standard (BS)	MP04565
O ₂ flowmeter 16 L/min, central supply, Norme Française (NF)	MP04566
O ₂ flowmeter 16 L/min, rail, NIST	MP04567
AIR flowmeter 16 L/min, rail, NIST	MP04568

Description	Part no.
O ₂ flowmeter 32 L/min, rail, NIST	MP04580
AIR flowmeter 32 L/min, rail, NIST	MP04581
O ₂ double flowmeter 16 L/min, rail, NIST	MP04569
O ₂ flowmeter 16 L/min, rail, 90°, NIST	MP04570
O ₂ flowmeter 32 L/min, central supply, DIN	MP04571
AIR flowmeter 32 L/min, central supply, DIN	MP04572
O ₂ flowmeter 16 L/min, hose, DIN ISO	MP04573
O ₂ flowmeter 16 L/min, hose, NIST-NIST ISO, anesthesia	MP04576
O ₂ flowmeter 16 L/min, hose neutral, NIST-NIST, anesthesia	MP04577
O ₂ flowmeter 16 L/min, 90°, hose, NIST-NIST ISO, anesthesia	MP04578
O ₂ flowmeter 16 L/min, 90°, hose neutral, NIST-NIST, anesthesia	MP04579
O ₂ flowmeter 16 L/min, hose, DIN neutral	MP04575

Replacement parts

00537050

WARNING

Do not disassemble the flowmeter while it is under pressure.

	Description	Part no.
A	Flow tube replacement set, 4.5 bar (65.3 psi)	
	O ₂ 4 L/min	MP04584
	O ₂ 10 L/min	MP04585
	O ₂ 16 L/min	MP04586
	AIR 16 L/min	MP04587
	O ₂ 32 L/min	MP04588
	AIR 32 L/min	MP04589

	Description	Part no.
B	Sealing ring	M04237
C	Pressure dome replacement set with O-ring	
	AIR	2M86557
	O ₂	2M86558
D	Hand wheel replacement set	MP04583
	Spindle replacement set, 16 L/32 L	2M86560
	Spindle replacement set, 4 L	MP00434

Recommended accessories for operation

Description	Part no.
Humidifier	
Humidifier bubbler	2M85834
Humidifier nebulizer	2M85835
Medication nebulizer	
Nebulizer element	2M85571
Nebulizer mask	M11015
Mouthpiece	M12394
O ₂ /Air connection hose (1.5 m/59.05 in)	M17716
O ₂ /Air connection hose (3 m/118.11 in)	M17717
O ₂ /Air spiral hose (3 m/118.11 in)	MP03322

Description	Part no.
Spiral hose Air/O ₂ replacement set, M15x1 (3 m/118.11 in)	MP03323
Adapter for Inhalette O ₂ (M34 x 1.5 to M15 x 1)	MP01184
Adapter for Inhalette Air (M32 x 1.5 to M15 x 1)	MP01185
Adapter	
Connector Ø 5 (M32 x 1.5; Air)	2M85575
Connector Ø 6 (M34 x 1.5; O ₂)	M06258
Adapter M34 x 1.5 to 9/16"	2M86680
Adapter 9/16" to 6 mm nozzle	MP00416
Hose configurator	8601697
The flowmeter is also compatible with the O2Star™ oxygen therapy accessories of Dräger. For further information contact Dräger.	

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© Drägerwerk AG & Co. KGaA